

PCR Purification:

Requirements:

- Omega E.Z.N.A.® Cycle Pure Kit
- PCR product
- 100% ethanol
- Microcentrifuge capable of at least 13,000 x g
- Nuclease-free 1.5mL microcentrifuge tubes
- Sterile deionized water
- Electric dry oven of 65°C
- Water bath of 65°C

Before Starting:

- Heat sterile deionized water to 65°C using water bath
- Add 100mL 100% ethanol to the bottle of DNA Wash Buffer if there's no mark on the bottle and store at room temperature

Protocol:

1. Determine the volume of PCR product, and transfer the product into a clean 1.5mL microcentrifuge tube.
2. Add 4-5 volumes CP Buffer. For PCR products smaller than 200bp, add 6 volumes CP Buffer.
3. Vortex to mix thoroughly.
4. Insert a HiBind® DNA Mini Column into a 2mL Collection Tube.
5. Add the sample from Step 3 to the HiBind® DNA Mini Column.
6. Centrifuge at maximum speed ($\geq 13,000 \times g$) for 1 minute at room temperature.
7. Discard the filtrate and reuse collection tube.
8. Add 700 μ L DNA Wash Buffer.
9. Centrifuge at maximum speed for 1 minute.
10. Discard the filtrate and reuse collection minute.
11. Centrifuge the empty HiBind® DNA Mini Column at maximum speed for 2 minutes to dry the column.
12. Transfer the HiBind® DNA Mini Column into a clean 1.5mL microcentrifuge tube. Open the lid and put it in the electric dry oven for 10 minutes to volatilize alcohol.
13. Add 40 μ L sterile deionized water directly to the center of column matrix.
14. Let sit at room temperature for 2 minutes.
15. Centrifuge at maximum speed for 1 minute.
16. Store DNA at -20°C.