

Protocol of RNA extraction

Kit content:

RNeasy Mini Spin Columns

Collection Tubes (2 ml)

Buffer RLT Buffer RW1 Buffer RPE

RNase-Free Water

Process:

1. Harvest the cells appropriate number of cells (maximum 1×10^9) by centrifuging at $5000 \times g$ for 3-5 min at 4°C . Discard supernatant ensuring all media is completely removed.
2. Resuspend cells in 100 μl lysozyme-containing TE buffer by vortexing.

	Lysozyme concentration in TE buffer	Incubation time at room temperature
Gram-negative bacteria	400 $\mu\text{g}/\text{ml}$	3-5 min
Gram-positive bacteria	3mg/ ml	5-10 min

3. Add 350 μl Buffer RLT and vortex vigorously to lyse the spheroplasts. If insoluble material is visible, centrifuge for 2 min at full speed, and use only the supernatant in subsequent steps.
4. Add 250 μl ethanol (96–100%) to the homogenized lysate, and mix well by pipetting without centrifugation.
5. Apply sample (usually 700 μl), including any precipitate that may have formed, to an RNeasy spin column placed in a 2 ml collection tube (supplied). Close the lid gently, and centrifuge for 15 s at $\geq 8000 \times g$ ($\geq 10,000$ rpm).
6. Add 700 μl Buffer RW1 to the RNeasy spin column, and centrifuge for 15 s at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to wash

7. Transfer RNeasy column into a new 2-ml collection tube. Pipet 500 μ l Buffer RPE onto the RNeasy column, and centrifuge for 15 s at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to wash
8. Centrifuge for 2 min at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to wash the spin column membrane.
9. Place the RNeasy spin column in a new 1.5 ml collection tube (supplied). Add 30–50 μ l RNase-free water directly to the spin column membrane. Close the lid gently, and centrifuge for 1 min at $\geq 8000 \times g$ ($\geq 10,000$ rpm) to elute the RNA.