

P1 lysates

- To a 150ml flask add:
 - o 25µl 1M CaCl₂
 - o 25µl overnight culture
 - o 5ml LB
 - o 2µl P1 stock (approximately 10⁸ pfu/ml)
 - o 50µl 1M MgSO₄
- Let grow shaking at 37°C until you see complete lysis (usually 4 hours)
- Add 10 drops of chloroform and vortex
- Let stand 10 minutes at room temperature
- Transfer lysate to Eppendorf tubes
- Centrifuge 2min
- Transfer supernatant to a new tube
- Add 2 drops of chloroform and vortex
- Let stand 10 min at RT
- Centrifuge 2min
- Transfer supernatant to a new tube, leaving the chloroform behind
- Centrifuge 2min
- Transfer supernatant to new tube.

P1 transduction

Note: if using P1 lysates from the fridge add 3 drops of chloroform, vortex and let stand 10min at RT then centrifuge 2mins and transfer to a new tube.

Remember to do controls where you don't add any lysates!

- Centrifuge overnight culture
- Resuspend cells in 1/2 volume (e.g. 1ml overnight, resuspend in 500µl) MC buffer
- Add 10-100µl P1 lysate (20µl mostly works) to 100µl cells
- Incubate 20mins at 37°C
- Add 150µl 1M NaCitate
- Add 1ml LB
- Incubate 1 hour at 37°C
- Centrifuge
- Resuspend pellet in 0,1ml LB (okay to leave some supernatant and resuspend in that) and spread all on a selective plate

Also spot 5µl P1 lysate, 5µl MC buffer, 5µl LB and 5µl NaCitate onto a LB plate containing no antibiotics to check for contamination.

MC buffer:

- 0,1M MgSO_4
- 0,005M CaCl_2
- H_2O

100 μl 1M MgSO_4 for 1ml MC

50 μl 0,1M CaCl_2 for 1ml MC

850 μl for 1ml MC