

## Genomic DNA Miniprep

- \* grow 5 ml cells over night at 30 °C
- \* spin, wash once with 1 ml H<sub>2</sub>O
- \* resuspend in 500 µl lysis buffer
- \* add acid washed glass beads to 1.25 ml
- \* vortex 2 min
- \* recover liquid phase with blue tip or  
punch whole into the bottom of the tube and centrifuge into another tube
- \* add 275 µl 7 M Ammonium acetate pH 7.0
- \* incubate 5 min at 65 °C, then 5 min on ice
- \* add 500 µl Chloroform, vortex, spin 2 min in microfuge
- \* take supernatant and precipitate with 1 ml Isopropanol
- \* incubate 5 min at RT, spin 5 min
- \* wash pellet with 70 % EtOH, dry and dissolve in 50 µl H<sub>2</sub>O
  
- \* for Southern: digest 5 µl DNA
- \* for PCR: use 0.5 to 1 µl DNA

### Buffers:

#### lysis buffer:

100 mM Tris pH 8.0; 50 mM EDTA; 1 % SDS.

[for 50 ml: 5 ml 1 M Tris; 5 ml 0.5 M EDTA; 5 ml 10 % SDS]