

Splitting HeLa / CHO cells

1. Detachment of adherent cells with trypsin
 1. Aspirate the medium from the cell culture dish
 2. Add 2 ml PBS and gently mix
 3. Aspirate the PBS
 4. Add 2 ml trypsin solution to the dish and incubate for 5 min at 37°C
 5. Detach cells by pipetting up and down
 6. Transfer cells to a 15 ml Falcon tube, filled with 5 ml medium
 7. Centrifuge for 3 min at 300g
 8. Discard (aspirate) the supernatant and resuspend the pellet in 5-10 ml medium