

## **Leishmania and DC2.4 cells culture**

In our project, we used wild-type clone 12-1 of *leishmania amazonensis* (RAT/BA/74/LV78) which doubly transfected with pX-*alad* and p6.5-*pbgd* (DT; 8–10) as model. *Leishmania* were routinely grown as promastigotes at 25 °C in medium 199 (Sigma) supplemented with 10% heat-inactivated fetal bovine serum (HIFBS) and 25mM HEPES to pH 7.4. Through two different antibiotics, G418 (100 µg/ml) and tunicamycin (20 µg/ml) to select the correct content of PBGD and ALAD, which are third and second enzymes in the heme biosynthetic pathway. The cells need to be passaged in a week or when reached the peak cell density of  $\sim 5 \times 10^7$  cells/ml at the stationary phase. Use T25 flask at 4ml/flask with solid screw-cap tightened to expand the cells. Stock cultures were set up by inoculation of a turbid culture at one drop ( $\sim 50$  µl) per ml into the 15 ml culture tubes. Before exposure of these DT cultures to ALA for inducing cytoplasmic accumulation of uroporphyrin 1 (URO), they were grown in drug-free medium to stationary phase to avoid the potential cytotoxicity of the drugs effect carried over to the host cells.

DC2.4 line, mouse dendritic cells were grown in RPMI-1640 supplemented with 10% FBS at 37 °C, 5% CO<sub>2</sub> incubator.