

E. coli heat-shock competent cells

Introduction

Protocol for preparing E. coli heat-shock competent cells from <http://www.genomearchitecture.com/protocols/ecoli-heat-shock-competent-cells-preparation.html>.

When preparing cells, it is important to keep everything chilled for use. For example, you should pre-chill your pipette tip boxes on ice before using.

Materials

- › E. coli strain
- › LB medium
- › LB plates with proper antibiotic
- › 0.1 M CaCl₂ solution (ice cold)
- › 0.1 M CaCl₂ solution containing 15 % glycerol (ice cold)

Procedure

Competent cells preparation

1. The day before: put CaCl₂ solutions at 5 °C, inoculate one single colony of your E. coli strain in 5 mL LB medium. Shake at +37 °C overnight.
 - When preparing cells; filter solutions before use.
2. Put 1.5 mL eppendorf tubes at -80 °C in advance.
3. Subculture at 37 °C with shaking till OD₆₀₀ reaches 0.25-0.3 (about 2 hrs subculture time)
4. Chill the culture on ice for 15 minutes.
5. Separate 100 mL chilled bacterial culture into two 50 mL falcon tubes and centrifuge at 4 °C at 4000 rpm for 10 minutes.
6. Discard the supernatant and resuspend the pellet with 40 mL ice-cold CaCl₂ solution.
7. Keep cells on ice again for 30 minutes.
8. Centrifuge cells at 4000 rpm at 4 °C for 10 minutes.
9. Discard the supernatant and resuspend pellet with 5 mL ice cold 0.1 M CaCl₂ solution containing 15 % glycerol.
10. Pipet 50 µL (or 100µL, or 200µL) of cell suspension into -80 °C frozen eppendorfs and directly transfer them to -80 °C freezer.