

# Restriction Digest

---

## Introduction

Basic protocol for restriction digest.

## Materials

- › DNA
- › Restriction enzymes
- › Buffer
- › dH<sub>2</sub>O

## Procedure

### Digest

1. Combine the materials:
  - >500 ng DNA (500 ng for diagnostic digest, >1000ng for restriction cloning)
  - 0.5-1 uL each restriction enzyme
  - Buffer (appropriate buffer indicated by enzyme manufacturer, to a final concentration of 1x)
  - dH<sub>2</sub>O to total volume of 20 ul.
2. Mix gently by pipetting.
3. Incubate tube for an appropriate temperature at an appropriate time (usually 37 C, time varies; generally 1 h for NEB enzymes, 10-30 min for FastDigest enzymes. Always follow the manufacturer's instructions.).
4. Visualize the results of your digest, conduct gel electrophoresis.