

Ligation

Rationale:	
Special Observations:	
Results:	
Interpretation :	

Experiment Date:	Source: NEB, Dujduan Waraho
Experiment Time:	
Primary Experimenter (contact):	Assembled: 6/27/2012
Other Experimenters:	

Reagent	Details	Quantity
ddH2O (nuclease-free)		Bring total vol. to 20 µL
10X T4 DNA Ligase Buffer*		2 µL
Vector DNA		**Vector : Insert = 1 : 3
Insert DNA		**Vector : Insert = 1 : 3
T4 Ligase		1 µL

*Thawed and resuspended at room temperature

** 1:3 molar ratio of vector : insert (but it can vary from 1:2 to 1:6);

Calculations:

Start with 10 ul vector and calculate volume of insert, then scale sum to ≤17uL.

$$Y_{ng\ insert} = (3) * (X_{ng\ vector}) * \frac{bp\ length\ of\ insert}{bp\ length\ vector}$$

$$Y_{ng\ insert} = (3) * (X_{ng\ vector}) * \frac{bp\ length\ of\ insert}{bp\ length\ vector}$$

$$Y_{ng\ insert} = (3) * (X_{ng\ vector}) * \frac{bp\ length\ of\ insert}{bp\ length\ vector}$$

$$\mu\text{L insert to use} = \frac{Y \text{ ng insert}}{\text{insert conc. in } \frac{\text{ng}}{\mu\text{L}}}$$

$$\mu\text{L vector to use} = \frac{(X \text{ ng vector})}{\text{vector conc. in } \frac{\text{ng}}{\mu\text{L}}}$$

Procedure:

Critical Steps:

Label microcentrifuge tubes, put on ice

Add components in listed order to microcentrifuge tubes on ice

For numerous ligations, prepare mastermix with water and buffer

Incubate at 4°C overnight and/or room temperature for at least 30 minutes

Chill on ice and transform 5 µl of the reaction into 50 µl competent cells.

- o See electroporation protocol