## Mowiol 4-88

Art.-Nr. 0713

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## Mowiol 4-88

Art. No. 0713

Mowiol 4-88 consists of colourless to light yellow granules. It is a partly hydrolysed polyvinyl alcohol with high adhesive power, which is principally used in industry working with adhesives.

In histology Mowiol 4-88 is used as hydrophilic mounting medium. Usually it is applied in combination with the anti-bleaching reagent DABCO (1.4-Diazabicyclo-(2.2.2)octan), which conserves the colours and the fluorescence of durable preparations.

### Instruction for hydrophilic embedding with Mowiol 4-88

#### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5  
76185 Karlsruhe  
Postfach 100121  
76231 Karlsruhe  
Telefon: +49 (0) 721/ 5606-0  
Telefax: +49 (0) 721/ 5606-149  
E-Mail: info@carlroth.de  
Internet: www.carlroth.de

e.d. 11/2014

#### Ingredients:

- 6.0 g Glycerol (Art. No. 3783)
- 2.4 g Mowiol 4-88
- 6.0 ml distilled water, sterile (Art. No. T143)
- 12.0 ml 0.2M Tris-HCl (pH 8.5) (Art. No. 9090)
- 25 mg per ml solution DABCO (Art. No. 0718) **(optional)**

#### Preparation:

- Mix 6.0 g Glycerol with 2.4 g Mowiol and dissolve with frequent agitation for one hour at room temperature (RT).
- Add 6.0 ml water and stir for one more hour at RT.
- Add 12.0 ml 0.2 M Tris-HCl (pH 8.5) and incubate for two hours at 50 °C under periodical stirring (every 20 minutes for 2 minutes).
- Note: In many cases Mowiol does not dissolve completely. We recommend centrifugation for 15 minutes at 5000 g. Continue with the supernatant.
- Optional:** For each ml solution, add 25 mg DABCO and stir until complete dissolution.
- Pour the medium into 1.5 ml- or 2 ml-tubes and store at -20 °C.

#### Application:

Thaw a Mowiol- or Mowiol/DABCO-aliquot and keep it on ice until use.

**Note: Mowiol/DABCO-aliquots should undergo 2 cycles of freeze-and-thaw in maximum.**

Apply the tissue specimen onto the slide and remove the excessive buffer. Don't allow the sample to fall dry.

Pipet 50-100 µl Mowiol or Mowiol/DABCO without bubbles onto a cover slip – preferably forming a „T“. Apply the cover slip bubble-free to the slide by lowering with caution (s. Fig.1). Let the samples harden over night.

**Completely hardened specimens can be stored at room temperature.**

Fig. 1

Lower the cover slip wet side down onto the slide starting at the cross of the „T“, in order to disperse the embedding medium optimally.

<b>Mowiol</b>	0713.1	50 g
	0713.2	100 g