

## RNeasy® Plus Mini Kit

The RNeasy Plus Mini Kit (cat. nos. 74134 and 74136) can be stored at room temperature (15–25°C) for at least 9 months.

For more information, additional and more detailed protocols, and safety information, please refer to the *RNeasy Plus Mini Handbook*, which can be found at [www.qiagen.com/handbooks](http://www.qiagen.com/handbooks).

For technical assistance, please call toll-free 00800-22-44-6000, or find regional phone numbers at [www.qiagen.com/contact](http://www.qiagen.com/contact).

### Notes before starting

- If purifying RNA from cell lines rich in RNases, or tissue, add either 10  $\mu$ l  $\beta$ -mercaptoethanol ( $\beta$ -ME), or 20  $\mu$ l 2 M dithiothreitol (DTT), to 1 ml Buffer RLT Plus before use. Buffer RLT Plus containing DTT or  $\beta$ -ME can be stored at room temperature for up to 1 month.
- Add 4 volumes of ethanol (96–100%) to Buffer RPE for a working solution.
- Foaming can be reduced by adding Reagent DX (cat. no. 19088) at a final concentration of 0.5% (v/v) before disruption and homogenization.\*  
\* This option not included in handbook; handbook to be updated.

1. **Cells:** Harvest a maximum of  $1 \times 10^7$  cells, either as a cell pellet, or lysed directly in the vessel. Add the appropriate volume of Buffer RLT Plus (see Table 1). Vortex for 30 s, or homogenize.  
**Tissues:** Disrupt the tissue ( $\leq 30$  mg) and homogenize the lysate in the appropriate volume of Buffer RLT Plus (see Table 1). Centrifuge the lysate for 3 min at maximum speed. Carefully remove the supernatant by pipetting and use it in step 2.
2. Transfer the homogenized lysate to a gDNA Eliminator spin column placed in a 2 ml collection tube (supplied).
3. Centrifuge for 30 s at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm). Discard the column, and save the flow-through. Add 1 volume (usually 350  $\mu$ l or 600  $\mu$ l) of 70%

January 2011



ethanol to the flow-through, and mix well by pipetting. Do not centrifuge. Proceed immediately to step 4.

4. Transfer up to 700  $\mu$ l of the sample, including any precipitate, to an RNeasy spin column placed in a 2 ml collection tube (supplied). Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
5. Add 700  $\mu$ l Buffer RW1 to the RNeasy Mini spin column (in a 2 ml collection tube). Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
6. Add 500  $\mu$ l Buffer RPE to the RNeasy spin column. Close the lid, and centrifuge for 15 s at  $\geq 8000 \times g$ . Discard the flow-through.
7. Add 500  $\mu$ l Buffer RPE to the RNeasy spin column. Close the lid gently, and centrifuge for 2 min at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm).  
**Optional:** Place the RNeasy spin column in a new 2 ml collection tube (supplied). Centrifuge at full speed for 1 min to further dry the membrane.
8. Place the RNeasy spin column in a new 1.5 ml collection tube (supplied). Add 30–50  $\mu$ l RNase-free water directly to the spin column membrane. Close the lid, and centrifuge for 1 min at  $\geq 8000 \times g$  to elute the RNA.  
**Optional:** Repeat elution with another volume of water or with RNA eluate.

**Table 1. Volumes of Buffer RLT Plus for sample disruption and homogenization**

| Sample         | Amount               | Dish     | Buffer RLT Plus* Disruption and homogenization |  |
|----------------|----------------------|----------|--|--|
| Pelleted cells | $< 5 \times 10^6$    | $< 6$ cm | 350 $\mu$ l                                    | Add Buffer RLT Plus, vortex ( $\leq 1 \times 10^5$ cells);                       |
|                | $\leq 1 \times 10^7$ | 6–10 cm  | 600 $\mu$ l                                    | or use QIAshredder, TissueRuptor®, or needle and syringe                         |
| Animal tissues | $< 20$ mg            | –        | 350 $\mu$ l                                    | TissueLyser LT; TissueLyser II;  |
|                | 20–30 mg             | –        | 600 $\mu$ l                                    | TissueRuptor, or mortar and pestle followed by QIAshredder or needle and syringe |

\* Use 600  $\mu$ l Buffer RLT Plus for tissues stabilized in RNAlater® Reagent, or for difficult-to-lyse tissues.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual.

“RNAlater®” is a trademark of AMBION, Inc., Austin, Texas and is covered by various U.S. and foreign patents.

Trademarks: QIAGEN®, RNeasy®, TissueRuptor® (QIAGEN Group);

1067551 01/2011 © 2011 QIAGEN, all rights reserved.

