

Preparation of nano-cellulose fiber aerogel

Experiment material: KP blanched pine leaves pulp, TEMPO, NaClO solution (contains 50% effective Cl), HCl solution (12mol/L), NaOH solution (1mol/L)

Equipment: magnetic stirrer, lyophilizer, liquid nitrogen

- 1, Pipette 100 ml 1% fiber pulp into a 250 ml beaker, then pipette 0.1g NaClO solution and certain amount of TEMPO solution. Stir the mixture on the magnetic stirrer.
- 2, Again, pipette certain amount of NaClO solution into the mixture. Use HCl solution to stabilize the PH of the system around 10.
- 3, During the reaction, continuously add NaOH solution to the system to maintain the PH as 10.
- 4, At the end of the reaction, filter the sample through a membrane which has a 0.22um bore diameter on the sanding core. Rinse the sample several times, and keep the sample in the refrigerator.
- 5, Disperse the oxidize fiber over distilled water. Adjust the concentration of the suspension to 0.2%. then treat the suspension with high pressure homogenization.
- 6, Freeze the nano-cellulose fiber suspicion with liquid nitrogen. Uses lyophilizer to dry the sample 48 hours at -55 centigrade.