

iGEM2014 – Microbiology – BMB – SDU	
<b>Title:</b> TSB transformation <b>SOP number:</b> SOP0009_v01 <b>Version number:</b> 01	<b>Date issued:</b> 2013.06.17 <b>Review date:</b> 2013.06.17 <b>Written by:</b> Patrick Rosendahl Andreassen

### 1. Purpose

To transform *E. coli* cells with plasmid using TSB buffer

### 2. Area of application

*All E. coli cells*

### 3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
Pipettes (p1000,200,10)	Micro Storage	•	
Heating block	Laboratory 1. Floor	•	
Ice	Across V18-403b-2	•	

### 4. Materials and reagents – their shelf life and risk labelling

Name	Components	Supplier / Cat. #	Room (hallway storage)	Safety considerations
Purple pipette tips		Contact lab-manager	Micro storage	
Green pipette tips		Contact lab-manager	Micro storage	
Blue pipette tips		Contact lab-manager	Micro storage	
Fort. LB		The new Anne-mette	Autoclave room	

<b>Polyethylene glycol (PEG) 3,350</b>		Sigma Aldrich	Micro chemical room	
<b>Dimethyl sulfoxid (DMSO)</b>		Sigma Aldrich	Micro chemical room	
<b>Magnesium chloride (MgCl<sub>2</sub>) 1M</b>		The new Anne-mette?	Autoclave room	
<b>Sterile filter (Pref. Blue)</b>		Contact lab-manager	Micro storage	
<b>Plasmid</b>				
<b>15mL falcon tube</b>		Contact lab-manager	Micro storage	
<b>10mL syringe</b>		Contact lab-manager	Micro storage	
<b>Long needle for syringe</b>		Contact lab-manager	Micro storage	

## 5. QC – Quality Control

Colony PCR on transformed cells using primers for the plasmid.

## 6. List of other SOPs relevant to this SOP

## 7. Environmental conditions required

## 8. Procedure

1. Preparation of E. coli culture:
  - 1.1. Add at least 5mL fort. LB (depending on amount of transformation to perform) to a bulb
  - 1.2. With a blue pipette tip add E. coli culture from agar plate to the LB media
  - 1.3. Grow culture to a OD600 of 0.3 to 0.5
2. Preparation of TSB buffer
  - 2.1. Add the following components to a 15mL falcon tube:
    - 2.1.1. PEG 3,350 1g
    - 2.1.2. DMSO 500µL
    - 2.1.3. MgCl<sub>2</sub> (1M) 200µL
    - 2.1.4. Fort. LB →10mL
  - 2.2. When everything is completely dissolved, transfer it to a new (sterile) falcon tube through  
  
a sterile filter using a syringe
3. TSB Transformation
  - 3.1. Spin 0.5-1.0mL culture for 5 min. at 4000 rpm.

- 3.2. Remove supernatant
- 3.3. Dissolve pellet in 200µL TSB buffer
- 3.4. Add plasmid (varying amount)
- 3.5. Keep at ice for 30 min.
- 3.6. Transfer directly to a heating block at 42°C for 2 min.
- 3.7. Add 1 mL fort. LB
- 3.8. Phenotypical expression at 37°C (0-2 hours)
- 3.9. Spin for 5 min. at 4000 rpm.
- 3.10. Remove most supernatant and dissolve pellet in the remaining supernatant (50-150µL)
- 3.11. Plate on agar plate with appropriate antibiotic

## 9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks
ON Culture		Liquid bacterial waste	
Once use plastic		GMO yellow waste	

## 10. Time consumption

- Total-time 4-6 hours.
- Hands-on-time 45 min.

## 11. Scheme of development

Date / Initials	Version No.	Description of changes
13.06.18 / PRA	01	The SOP has been written
13.06.18 / AK	01	The SOP has been approved

## 12. Appendices