

Virus Transduction to 293FT cells

1. Seed cells to be 80% confluent at a 75mL culture flask.
2. Dilute 1ml lentiviral vector (packaged) in 10ml high glucose DMEM medium containing 10% FBS
3. Withdraw culture medium from the 75mL culture flask.
4. Add vector-DMEM complex to cells
5. Incubate for 12 hours in 37C°.
6. Withdraw vector-DMEM complex from culture flask.
7. Add 10ml DMEM medium containing 10% FBS to cells and incubate.
8. Observe the cells under Inverted fluorescence microscope.
9. Collect the medium containing virus vector after 48 hours and 72 hours, respectively.