

## Digestion

1.Enzyme A 1 $\mu$ l

Buffer (10X) 0.2 $\mu$ l

BSA(100X) 2 $\mu$ l

ddH<sub>2</sub>O rest of the volume

incubate at recommended temperature (37°C) for at least 2 hour, and temperature (68°C) for 10 mins

2.Purify the digestion product;

3.Enzyme B 1 $\mu$ l

Buffer (10X) 0.2 $\mu$ l

BSA(100X) 2 $\mu$ l

ddH<sub>2</sub>O rest of the volume

incubate at recommended temperature (37°C) for at least 2 hour, and temperature (68°C) for 10 mins

4.Purify the digestion product;

Notes:

The enzymes used here are NEB enzymes ( EcoRI / XbaI / SpeI / PstI), and buffer4 is suitable for most of double digestion;

For 50  $\mu$ l reaction system, the suggested amount of each restriction enzyme is 2 $\mu$ l;

According to personal experience, 50 $\mu$ l reaction system has lower efficiency than 20  $\mu$ l reaction system, so 20 $\mu$ l reaction system is recommended.