

Protocol for transformation with ligation reaction system

1. Get the competent cell from -70°C and wait for its fusion;
2. Add $20\mu\text{l}$ TFBII and $20\mu\text{l}$ into a 1.5ml centrifuge tube;
3. Add the entire ligation product ($10\mu\text{l}$) into the tube;
4. Mix and incubate on ice for 30 min;
5. Heat pulse for 1min at 42°C ;
6. Put back the tube on ice and incubate for 2-3min;
7. Add $900\mu\text{l}$ LB non-antibiotic agar plates and incubate at 37°C for 60min;
8. Plate the culture on LB plate containing corresponding antibiotics