

3A ASSEMBLY (DIGESTION AND LIGATION)

Aim

To assemble two BioBricks into a new vector.

Procedure

Digestion

1. Prepare following Master Mixes

Table 1: Master Mixes		
Enzyme Master Mix for Plasmid Backbone (25 μ l total, for 5 runs)	Enzyme Master Mix for Part A (BioBrick on the 5' end)	Enzyme Master Mix for Part B (BioBrick on the 3' end)
5 μ l NEB Buffer 2	5 μ l NEB Buffer 2	5 μ l NEB Buffer 2
0.5 μ l BSA	0.5 μ l BSA	0.5 μ l BSA
0.5 μ l EcoRI-HF	0.5 μ l EcoRI-HF	0.5 μ l XbaI
0.5 μ l PstI	0.5 μ l SpeI	0.5 μ l PstI
18 μ l dH ₂ O	18.5 μ l dH ₂ O	18.5 μ l dH ₂ O

2. Digest Plasmid Backbone

- Add 4 μ l linearized plasmid backbone (25 ng/ μ l for 100 ng total)
- Add 4 μ l of Enzyme Master Mix

3. Digest Part A

- Add 4 μ l Part A (25 ng/ μ l for 100 ng total)
- Add 4 μ l of Enzyme Master Mix

4. Digest Part B

- Add 4 μ l Part B (25 ng/ μ l for 100 ng total)
- Add 4 μ l of Enzyme Master Mix

5. Digest all three reactions at 37 °C for 60 min, heat kill 80 °C at 20 min

Ligation

1. Add 2 μ l of digested Plasmid Backbone (25 ng)
2. Add equimolar amount of Part A (EcoRI-HF SpeI digested) fragment (< 3 μ l)*
3. Add equimolar amount of Part B (XbaI PstI digested fragment) (< 3 μ l)*
4. Add 1 μ l T4 DNA ligase buffer.
5. Add 0.5 μ l T4 DNA ligase
6. Add water to 10 μ l

7. Ligate at RT for 10 min, heat kill 80 °C/20 min

8. Transform with 1-2 μ l of product

* Easiest to calculate using http://www.insilico.uni-duesseldorf.de/Lig_Input.html

Use a 1:3 vector:insert ratio

Sources

This protocol is modified from the 3A Assembly protocol at igem.org:

http://parts.igem.org/Help:Assembly/3A_Assembly

N.B. We used exceedingly more DNA than iGEM did - DNA volumes of 500 - 1000 μ l.