

PCR AMPLIFICATION USING Q5[®] HIGH-FIDELITY 2X MASTER MIX

Aim

Selectively amplify regions of interest from a linear or circular DNA molecule.

Procedure

1. Use of high quality, purified DNA templates greatly enhances the success of PCR. Recommended amounts of DNA template for a 50 μ l reaction are as follows:

Table 1: Amount of DNA template

DNA	AMOUNT
DNA genomic	1ng - 1 μ l
Plasmid or viral	1pg - 1ng

2. Gently mix the reaction mix described below. Collect all liquid to the bottom of the tube by a quick spin if necessary.

Table 2: Reaction Setup

Component	25 μ l Reaction	50 μ l Reaction	Final concentration
Q5 High-Fidelity master mix	12.5 μ l	25 μ l	1x
10 μ M Forward Primer	1.25 μ l	2.5 μ l	0.5 μ M
10 μ M Reverse Primer	1.25 μ l	2.5 μ l	0.5 μ M
Template DNA	Variable	Variable	<1.000 ng
Nuclease-Free Water	to 25 μ l	To 50 μ l	

Table 3: Thermocycling Conditions for a Routine PCR

STEP	TEMP	TIME
Initial Denaturation	98°C	30 sec
25-30 Cycles	98°C	5-10 sec
	50-72°C	10-30 sec
	72°C	20-30 sec
	72°C	2 min
Final Extension	72°C	2 min
Hold	4-10°C	

3. Transfer PCR tubes to a PCR machine and begin thermocycling.

Sources

To access the original protocol, [click here](#).