

Nuclease

Week 3

Summarized below are the experiments conducted this week in chronological order. Click on the experiment name to view it. To go back to this summary, click **Summary** in the footer.

Summary

1	Digestion of T7 & Nuc	2
2	PCR of Nuc from Glycerol Stocks	4

Digestion of T7 & Nuc

Responsible

Aman Mebrau, Oscar Frisell, Sigrun Stulz, Bethel Tesfai Embaie

Protocols used

- 3A Assembly protocol
- Gel purification

Modifications and comments to protocols

Not enough Nuc samples were available and the PstI single digestion was not performed.

Table 1: Modifications to digestion reactions

N4	N - (neg)
2 μ l DNA	2 μ l DNA
0.6 μ l 10x cutsmart	0.6 μ l 10x cutsmart
3.4 μ l ddH ₂ O	3.4 μ l ddH ₂ O
2 μ l enzyme MM	

Experimental setup

The digestion followed the protocol. See table below for Sample structure.

Table 2: Sample structure for T7 and Nuc 2

T7	T1	T2	T3	T4	T5	T -
Conc 35.8 ng/ μ l	sd EcoRI	sd SpeI	dd	dd	dd	negative control
Nuc 2	N1	N2	N3	N4	N -	
Conc 30.05 ng/ μ l	dd (XbaI & PstI)	dd	dd	sd XbaI	negative control	

Sample Calculation

$$\text{Volume of DNA sample} = \frac{\text{Final amount of DNA}}{\text{Concentration of DNA in sample}} \quad (1)$$

Results and Conclusions

Gel 1: ran for too long; samples were lost

Gel 2: gel was repeated with remaining T7 digestion products

Table 3: Gel labelling

Well	Label	Comment
1	DNA ladder	
2	T7 sd EcoRI	successful
3	T7 sd SpeI	successful
4	T7 dd	similar length; could be used for ligation
5	T7 dd	similar length; could be used for ligation
6	T7 dd	similar length; could be used for ligation
7	T7 negative control	looks right
8	ladder	




Figure 1: Successful gel of T7 digestion products.

Discussion and Troubleshooting

The Nuc gel ran for too long so we lost it. Hence we need to amplify more Nuc and run a new gel.

2 PCR of Nuc from Glycerol Stocks

Responsible

Bethel Tesfai Embaie, Oscar Frisell

Protocols used

- PCR protocol

Experimental Data

The aim was to amplify Nuc DNA for 3A assembly.

The total reaction volume was 50 μ l for each of the 4 per tubes.

The table below shows how much DNA template that was used:

Table 4: Amount of DNA template per reaction

Label	DNA template [μ l]
1	2
2	1
3	1
4	1

Results and Conclusions

Table 5: Gel description

Well	Component	Expected bands
1	DNA ladder	-
2	Nuc PCR product 2ul	Between 800-900 bp
3	Nuc PCR product 1ul	Between 800-900 bp
4	Nuc PCR product 1ul	Between 800-900 bp
5	Nuc PCR product 1ul	Between 800-900 bp




Figure 2: Gel of Nuc PCR products.

The bands in Figure 2 seem to have a size just below 1000 bp, however it is hard to tell as the bands are large and a bit smeared. Considering the expected size of between 800-900 bp, the bands might correspond to the PCR product of Nuc and thus will be PCR purified for further analysis.

Discussion and Troubleshooting

PCR purification will be followed by NanoDrop and digestion for 3A assembly.