

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU_Eindhoven

Digestion

Table of contents

Digestion

1

Digestion

3

1.1

Materials

3

1.2

Setup & Protocol

Error! Bookm

1 Digestion

Estimated bench time: 30 minutes

Estimated total time: 5 hours

Purpose: Giving sticky ends to the vector and the insert(s) so that they can be ligated.
It is essential to work with gloves at all times to protect the DNA from DNase activity.

1.1 Materials

- 10X NEB Cut-Smart buffer
- Autoclaved H₂O
- Autoclaved PCR tubes
- Bucket with ice
- Insert(s) which is/are to be digested (pcr amplification product)
- Pipettes and tips
- Thermal cycler
- Vector which is to be digested
- Restriction enzyme BamHI-HF
- Restriction enzyme HindIII-HF

1.2 Setup & Protocol

Properties of the restriction enzymes can be found at www.neb.com

- Please check the stock concentration of the restriction enzyme you use
- Use PCR tubes. (Evt. distribute the digestion mixtures over >1 tube)
- Add all components, end with restriction enzymes (transport in cold block)
- mix well

Reaction mixture insert		
Component	Quantity/mass/final concentration	Volume (µl)
DNA (PCR product ng/µL)	PCR product: ~3 ug	
10x NEB CutSmart buffer	1x	5
Restriction enzyme BamHI-HF (20 U/µL stock)	10 U	1
Restriction enzyme HindIII-HF (20 U/µL stock)	10 U	1
H ₂ O		
Total		50

Reaction mixture vector		
Component	Quantity/mass/final concentration	Volume (µl)
DNA pET28a (vector ng/uL)	Vector: 5 ug	
10x NEB CutSmart buffer	1x	5
Restriction enzyme BamHI-HF (20 U/uL stock)	20 U	1
Restriction enzyme HindIII-HF (20 U/uL stock)	20 U	1
H ₂ O		
total		50

- Start PCR machine using program shown below

Step	Temp (°C)	Time (min)
Incubation	37	240 (4 hours)
Cooling	4	hold