

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU_Eindhoven

Gibson Assembly

Table of contents

Gibson Assembly

1	NEBuilder HiFi DNA Assembly Reaction	Error! Bookm
1.1	Materials	Error! Bookm
1.2	Setup & Protocol	Error! Bookm
2	Transformation into NEB 5-alpha competent E. Coli	4
2.1	Materials	4
2.2	Setup & Protocol	4
3	References & Acknowledgements	4

1 NEBuilder HiFi DNA Assembly Reaction

Estimated bench time: 30 minutes

Estimated total time: 50 minutes

Purpose: The assembly of multiple inserts within the vector of choice.

1.1 Materials

- Autoclaved dH₂O
- Bucket with ice
- gBlocks
- Linearized vector
- NEBuilder HiFi DNA Cloning Kit
- PCR machine

1.2 Setup & Protocol

- Look up the following parameters to complete the table. The molecular weight of the double stranded vector and gBlocks can be determined using a DNA molecular weight calculator online.

Component	Molecular weight (Da)	Concentration (ng/μl)
<i>Linearized vector</i>		
<i>gBlock 1</i>		10
<i>Additional gBlocks</i>		10

- The amount of pmol for the Gibson assembly has to be between the following ranges:
 1 or 2 gBlocks: 0.02 – 0.5 pmol (with a vector : gBlock ratio of 1:2 or 1:3)
 4 to 6 gBlocks: 0.02 – 1.0 pmol (with a vector : gBlock ratio of 1:2 or 1:3)

$$\text{pmols} = (\text{weight in ng}) \times 1,000 / (\text{base pairs} \times 650 \text{ daltons})$$

- Set up the following reaction on ice:

Component	Value in mixture	Volume (μl)
<i>Linearized vector</i>	50-100 ng	
<i>gBlock 1 (10 ng/μl)</i>	3-fold excess (3x pmol vector)	
<i>Additional gBlocks (10 ng/μl)</i>	3-fold excess	
<i>NEBuilder HiFi DNA Assembly Master Mix</i>	10 μl	10
<i>Autoclaved dH₂O</i>	Till 20 μl	
<i>Total volume</i>		20

The volume of linearized vector and gBlock can be calculated with:

$$\text{Volume } (\mu\text{l}) = \frac{\text{amount (pmol)} * \text{molecular weight (Da)}}{\text{concentration } (\frac{\text{g}}{\text{l}})}$$

- For 2-3 gBlocks: Incubate the samples at 50 °C for 15 minutes using the PCR machine.
- For 4-6 gBlocks: Incubate the samples at 50 °C for 15 minutes using the PCR machine.
- After incubation, store the samples on ice for subsequent transformation.

2 Transformation into NEB 5-alpha competent E. Coli

Estimated bench time: 1 hour

Estimated total time: 2 hours

Purpose: Transformation of the assembled vector into competent cells. These cells can be used for subsequent plasmid amplification.

During transformation with NEB 5-alpha you are working with bacterial cells. Therefore, you need to work near the Bunsen burner flame, prohibiting the use of gloves during this step.

2.1 Materials

- Bucket with ice
- Bunsen Burner
- Heat/shaking-block
- Incubator
- LB-agar plates supplemented with the correct antibiotic
- NEB 5-alpha competent cells
- Pipettes and tips
- Plasmids to be transformed
- SOC solution
- Water bath

2.2 Setup & Protocol

- Thaw the NEB 5-alpha competent E. Coli on ice such that all ice crystals disappear.
- Add 2 µl of the chilled samples to the competent cells. Mix gently by stirring with your pipette. Do not vortex nor pipette up and down.
- Place mixture on ice for 30 minutes. Do not mix.
- Heat shock the cells for exactly 30 seconds at 42 °C.
- Transfer the tubes directly to ice for 2 minutes.
- Add 950 µl of room-temperature SOC media to the tube.
- Incubate for 60 minutes at 37 °C and 300 rpm.
- Spread 100 µl of the cells onto the agar plates while working near the Bunsen burner flame. Incubate overnight at 37 °C.

3 References & Acknowledgements

This protocol was adapted from the NEBuilder HiFi Assembly Cloning Kit's manual with help from 2015 team of iGEM TU Eindhoven.