

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU_Eindhoven

Mutagenesis

Table of contents

Mutagenesis

1 Site-Directed Mutagenesis

1.1 Materials

1.2 Setup & Protocol

2 DpnI digestion

2.1 Materials

2.2 Setup & Protocol

3 Transformation into E. Coli XL10-GOLD

4 Plating cells on agar

Error! Bookm

Error! Bookm

Error! Bookm

Error! Bookm

Error! Bookm

Error! Bookm

4

4

1 Site-Directed Mutagenesis

Estimated bench time: 15 minutes

Estimated total time: 2 hours, but much longer if there is a big insert in pET28a

Purpose: Introducing a mutation into the plasmid

1.1 Materials

- 10X buffer
- Autoclaved H₂O
- Bucket with ice
- dNTP mix
- Pair of primers which yield the mutation for the insert
- PCR tubes
- Pipettes and tips
- Plasmid DNA
- QuikChange Lightning Multi enzyme blend
- QuikSolution
- PCR machine

1.2 Setup & Protocol

- Keep dNTP mix on ice before addition to reaction mixture.
- Add the components to reaction mixture in the order shown below:

Component	Final Concentration	Volume (μL)
<i>H₂O</i>	Add up to 25μL	
<i>10X buffer</i>	1X	2.5
<i>QuikSolution</i>		0.5
<i>Plasmid DNA</i>	100ng (stock ..ng/μL)	
<i>Primer 1</i>	100ng (stock 10 μM)	1
<i>dNTP mix</i>		1
<i>QuikChange Lightning Multi enzyme blend</i>		1
<i>Total</i>		25

- Take the DNA polymerase out of the -20 °C freezer only after all components have been added to the reaction mixture and the PCR machine has been programmed. Do not touch the bottom of the polymerase tube.
- Mix reaction mixture well before start of PCR.
- Run the following PRC program:

Step	Temperature (°C)	Time (sec.)	Cycles
<i>Initial Denaturation</i>	95	120	1
<i>Denaturation</i>	95	20	30
<i>Annealing</i>	55	30	
<i>Extension</i>	65	30/kb	
<i>Final extension</i>	65	300	1
<i>Hold</i>	4		

2 DpnI digestion

Estimated bench time: 5 minutes

Estimated total time: 65 minutes

Purpose: Digestion of the template plasmid from the PCR product mixture.

2.1 Materials

- DpnI restriction enzyme
- Bucket with ice
- PCR product
- PCR machine

2.2 Setup & Protocol

- Thaw the DpnI restriction enzyme on ice.
- Add 1 µL of DpnI restriction enzyme directly to each amplification reaction.
- Incubate for 1 hour at 37 °C, in PCR machine.

3 Transformation into E. Coli XL10-GOLD

Estimated bench time: 40 minutes

Estimated total time: 2.5 hours

Purpose: Placing the plasmid into bacteria for amplification.

For more information, see our general Transformation into E. Coli XL10-GOLD protocol.

4 Plating cells on agar

Estimated bench time: 15 minutes

Estimated total time: 18 hours

Purpose: Amplification of the ligation product.

For more information, see our general Plating protocol.

- Be sure to dilute the bacteria 10x with SOC!
- SOC : bacteria = 10:1 - you have 50 µL bacteria, so add 500 µL SOC
- Plate 100 µL