

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU_Eindhoven

PCR Amplification

Table of contents

PCR Amplification	1	PCR Amplification	Error! Bookm
	1.1	Materials	Error! Bookm
	1.2	Setup & Protocol	Error! Bookm

1 PCR Amplification

Estimated bench time: 30 minutes

Estimated total time: 2 hours

Purpose: Amplification of DNA with the possibility of expanding the DNA sequence at the beginning and/or end with the primers.

It is essential to work with gloves at all times to protect the DNA from DNase activity.

1.1 Materials

- Autoclaved H₂O
- Autoclaved PCR tubes
- Bucket with ice
- DNA to be amplified
- Forward primer
- Pipettes and tips
- Reverse primer
- 5x phusion buffer
- Phusion DNA polymerase
- Thermal cycler

1.2 Setup & Protocol

- Construct 3 PCR mixtures in the following way:

Component	Quantity/mass/final concentration	Volume (μL)
DNA	30 ng (stock n ng/μL)	30/n
Primer FW	0.5 uM (10 uM stock)	2.5
Primer RW	0.5 uM (10 uM stock)	2.5
5x Phusion buffer	1x	10
dNTPs	200 uM (10 mM stock)	2.5
Phusion DNA polymerase	0.02 U/μL (2 U/μL stock)	0.5
H₂O		30-30/n
Total		50

- Use PCR tubes
- Add all components, end with polymerase (transport in cold block)
- Mix well by pipetting up and down.
- Run the following PCR program:

Step	Temperature (°C)	Time (sec)	Cycles
Denaturation	98	30	1
Denaturation	98	10	30
Annealing*	...	15	
Extension*	72	.. (20 sec/kb)	
Final extension	72	600	1
Cooling	4	hold	1

*differs for each PCR dependent on T_m primers

The annealing temperature can be calculated for the set of primers using New England Biolabs T_m calculator. An annealing temperature of 3°C lower than the lowest melting temperature was used to increase yields.