

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU_Eindhoven

Vector Linearization

Table of contents

Vector Linearization

1	Vector Linearization through PCR	
1.1	Materials	
1.2	Setup & Protocol	
2	DpnI digestion	
2.1	Materials	
2.2	Setup & Protocol	
3	PCR Purification	4
4	NanoDrop	4
5	Gel Electrophoresis	5
6	References & Acknowledgements	5

Error! Bookm
Error! Bookm
Error! Bookm

Error! Bookm
Error! Bookm
Error! Bookm

1 Vector Linearization through PCR

Estimated bench time: 45 minutes

Estimated total time: 5-7 hours (depends on the vector)

Purpose: Preparing a linear vector which can be used in the Gibson Assembly reaction.

1.1 Materials

- Autoclaved H₂O
- Bucket with ice
- Pair of primers which yield the necessary overlaps for the insert
- PCR tubes
- Pipettes and tips
- Q5 High-Fidelity 2X Master Mix
- PCR machine
- Vector
- Vortex

1.2 Setup & Protocol

- Thaw the Q5-HF 2X Master Mix on ice. If the Master Mix contains a pellet, briefly vortex the tube until the pellet disappears.
- Set up a PCR with the following reaction components for the vector to be amplified. Add the Q5-HF 2X Master Mix lastly! Quickly transfer the PCR tube to the PCR machine adding the polymerase:

Component	Stock Concentration	Final Concentration	Volume (ul)
<i>Plasmid</i>	1 ng/ul	1 ng	2
<i>Forward primer</i>	10 uM	0.5 uM	2.5
<i>Reverse primer</i>	10 uM	0.5 uM	2.5
<i>Q5 High-Fidelity 2X Master Mix</i>	2X	1X	25
<i>H₂O</i>			18
<i>Total</i>			50

- Run the following PCR program:

Step	Temperature (°C)	Time (sec.)	Cycles
<i>Initial Denaturation</i>	98	120	1
<i>Denaturation</i>	98	15	35
<i>Annealing*</i>	72 ¹	20	
<i>Extension*</i>	72	30/kb → 162	
<i>Final extension*</i>	72	120	1
<i>Hold</i>	4		

*Differs per PCR dependent on T_m of primers!

¹ The annealing temperature can be calculated for the set of primers using New England Biolabs T_m calculator.

2 DpnI digestion

Estimated bench time: 5 minutes plus 1 minute per sample

Estimated total time: 1.5 hours

Purpose: Digestion of the template vector from the PCR product mixture. This will remove the number of background colonies which do not carry the desired insert after Gibson Assembly.

2.1 Materials

- 10X CutSmart buffer from New England Biolabs
- Bucket with ice
- DpnI restriction enzyme
- PCR product
- PCR machine

2.2 Setup & Protocol

- Thaw the 10X CutSmart buffer at room temperature.
- Thaw the DpnI restriction enzyme on ice. Setup the following reaction:

Component	Stock concentration	Final concentration	Volume
<i>PCR product</i>			50
<i>10X CutSmart buffer</i>	10X	1X	5.7
<i>DpnI</i>	20U/ul	20U	1
<i>Total</i>	2X	1X	56.7

- Digest the vector for 1 hour at 37 °C. Heat inactivate DpnI for 20 minutes at 80 °

3 PCR Purification

Estimated bench time: 45 minutes

Estimated total time: 45 minutes

Purpose: If the PCR product is less than 90% pure, large volumes of unpurified PCR product could significantly inhibit the Gibson Assembly. PCR purification may be performed to increase the efficiency.

For more information, see our general PCR Purification protocol.

4 NanoDrop

Estimated bench time: 5 minutes start-up and 2 minutes per sample

Estimated total time: 5 minutes start-up and 2 minutes per sample

Purpose: Measuring the concentration of the PCR product which is necessary to set up the Gibson Assembly reaction.

For more information, see our general NanoDrop protocol.

5 Gel Electrophoresis

Estimated bench time: 40 minutes

Estimated total time: 1.5 hours

Purpose: Agarose gel electrophoresis may be used to verify the purity of your PCR product. If the product is pure, a single band will show up during the gel electrophoresis.

For more information, see our general Gel Electrophoresis protocol.

6 References & Acknowledgements

This protocol was based on the information provided by the 2015 iGEM team of the TU/Eindhoven, New England Biolabs NEBuilder HiFi DNA Assembly Cloning Kit manual and on Integrated DNA Technologies' g Blocks Gene Fragments Cloning Protocols.