

# PCR

## *Materials and equipment*

- Thermocycler
- Microtubes (0.2 ml)
- Micropipette tips
- Micropipettes

## *Reagents*

- Template DNA (cassette plasmid extraction)
- Standard reaction buffer 5X
- dNTPs
- Forward primer (VF2)
- Reverse primer (VR)
- DNA polymerase
- Nuclease-free water

## *Methodology*

1. The procedure is as follows:
2. Thaw the template DNA, Standard reaction buffer 5X, dNTPs and primers (FP, RP) on ice.
3. Add in a 0.2 ml tube (for a 25 $\mu$ L reaction):
4. Up to 25  $\mu$ L nuclease-free water
5. 5 $\mu$ L Standard Reaction buffer
6. 0.5  $\mu$ L of 10 mM dNTPs
7. 0.5  $\mu$ L of 10  $\mu$ M Forward Primer
8. 0.5  $\mu$ L of 10  $\mu$ M Reverse Primer
9. 1 $\mu$ L of Template DNA
10. 0.125  $\mu$ L DNA Polymerase
11. Mix gently each tube (by pipetting up and down several times).
12. Place at thermoblock in a program, based on the polymerase and primers specifications:
13. Denaturation step: 94°C, 5 minutes
14. Second Denaturation: 94°C, 30 seconds
15. Annealing: 56°C, for 30 seconds
16. Extension: 68°C, 1 minute
17. Final extension: 68°C, 5 minutes
18. Hold: 4°C, infinite time
19. Repeat for 30-35 cycles from b-d.
20. Store PCR products at -20°C.