

Plasmid digestion

Materials and Equipment

- Alcohol lamp
- Microtube rack
- Micropipette
- Micropipette tips
- Microtubes (2 ml)
- Thermoblock

Reagents

- Nuclease free water
- NEB Buffer 10X
- Plasmid
- Restriction Enzymes

Methodology

Each reaction will be added in a 0.2ml tube (for a 50 μ L reaction).

Control 1:

- A. Add up to 50 μ L nuclease free water (32 μ L)
- B. 5ul Buffer 10X
- C. 500 ug plasmid (12 μ L)
- D. 1 μ L Restriction Enzyme 1

Control 2:

- A. Add up to 50 μ L nuclease free water (32 μ L)
- B. 5ul Buffer 10X
- C. 500 ug plasmid (12 μ L)
- D. 1 μ L Restriction Enzyme 2

Negative Control:

- A. Plasmid without enzymes

Digestion:

- A. Add up to 50 μ L nuclease free water (32 μ L),
- B. 5ul Buffer 10X
- C. 500 ug plasmid (12 μ L)
- D. 1 μ L Restriction Enzyme 1
- E. 1 μ L Restriction Enzyme 2

1. Mix gently each tube.
2. Place at thermoblock, 37°C, 1 hour.
3. Inactivate enzymes at 80°C for 20 minutes.
4. Store digestions at -20°C.