

Synthesis resuspension

Materials and Equipment

- Micropipette
- Micropipette tips
- Centrifuge
- Vortex
- Termoblock

Reagents

- Nuclease free water
- Oligonucleotides (Reverse and Forward)

Methodology

1. From lyophilized synthesis tube:
2. Centrifuge tube at 5,000 rpm for 3 minutes.
3. Resuspend each biopart with nuclease-free water
4. Mix each tube with the Vortex for 5 seconds. (3 times)
5. Place microtubes in a boiling water bath at 50°C for 10 minutes
6. Store bioparts at -20°C.