

Glucose cloning protocol

First step:

Primers arrived => Dilute them as described from the company

Second step:

1. 2-stage PCR reaction for:

variant	F-primer	R-primer	Tm1	Tm2
1	Glu_BA_For	Glu_AA_Rev	48.6	65.7
2	Glu_BB_For	Glu_AB_Rev	55.6	68.3
3	Glu_BB_For	Glu_AC_Rev	55.6	68.3
4	Glu_BC_For	Glu_AD_Rev	55.6	66.4
5	Glu_BD_For	Glu_AB_Rev	55.6	68.3
6	Glu_BE_For	Glu_AE_Rev	48.6	67
7	Glu_BF_For	Glu_AF_Rev	55.6	66.9
8	Glu_BG_For	Glu_AF_Rev	55.6	66.2
control	Glu_BD_For	Glu_AB_Rev	55.6	68.3

Set up reaction:

Component	Volume	control
5x PCR Buffer	10 µl	10 µl
dNTPs (10mM)	1 µl	1 µl
Forward primer [10 µM]	2.5 µl	2.5µl
Reverse primer [10 µM]	2.5 µl	2.5µl
Template DNA - Tar2 CM (2- 10 ng)	1 µl	-
Phusion DNA Polymerase	0.5 µl	0.5
UPW	32.5 µl	33.5 µl
Total	50	50

PCR program:

cycles	Step	Temperature	Time
1	Initial denaturation	98	30 sec
10	Denaturation	98	10 sec
	Annealing	grad of 48.6-55.6	30 sec
	Extension	72	2 min
20	Denaturation	98	10 sec
	Annealing	grad of 65.7-68.3	30 sec
	Extension	72	2 min

1	Final extension	72	10 min
1	Hold	4	hold

30 cycles

- Gel electrophoresis (p.10-12) - 1% agar (500 bp to 10 kb)

5 µl of ladder

4 µl of sample + 1 µl dye

ladder	1	2	3	4	5	6	7	8	control
--------	---	---	---	---	---	---	---	---	---------

Expected length-

THE WHOLE TAR- 3930bp

- Purification of PCR product (Equal volume = 45 µl)

variant	conc.(ng/µl)	260/280	260/230
1			
2			
3			
4			
5			
6			
7			
8			

- DpnI for DNA (in order to digest methylated plasmid)

Set up reaction:

component	Volume(ul)
DNA	30
10X CutSmart Buffer	4
DpnI	1
UPW	5
Total	40ul

Incubate for 1 hour in 37°C

- Purification of DpnI product

variant	conc.(ng/μl)	260/280	260/230
1			
2			
3			
4			
5			
6			
7			
8			

6. Phosphorylation

Set up reaction:

component	Volume(μl)
T4 DNA ligase buffer	2.5 μl
T4 kinase (PNK)	1 μl
DNA	X~100 ng= 3μl
PEG 4000	2.5 μl (5%)
MB H ₂ O	19-X μl
Total	25 μl

Incubate for 30 minutes in 37°C

Heat inactivation prior to ligation for 20 minutes in 65°C

Cool on ice for 10 min

7. Self Ligation:

Add 1 μl of T4 DNA ligase to phosphorylation reaction tube

Incubate for 2 hours in 30°C

8. Transformation to UU

9. Starters + mini prep

10. Sequencing

11. Chemotaxis checks: Swarming assay, microscope