

PctA/ NarX protocol

1. **Extraction of Terminator part** from DNA 2016 extraction kit.
Part B0015- 2016 kit, plate 3, 3F.
2. **Transformation to *E. coli* Top 10.**
Glycerol stock + miniprep.
3. **Restriction of Terminator part with EcoRI & XbaI.**

Materials:

- Restriction Enzymes: EcoRI & XbaI
- NEB buffer: cutSmart
- Molecular biology water

Procedure:

Set up the following reaction as described below:

Component	Volume
Buffer (X10)	4µl
DNA 2.5µg	33µl
Enzyme 1 -HF	1µl
enzyme 2 -HF	1µl
MB water	1 µl
Total volume	40µl

[B0015]= ng/µl

1µg (B0015)= 13.5 µl (~µl)

2.5µg (B0015)= 33 µl (~µl)

→ MB water=1 µl

1. Add appropriate amount of MB water to sterile 1.5ml tube.
2. Add restriction enzyme buffer to the tube. Vortex buffer before pipetting to ensure that it is well-mixed.
3. Add appropriate amount of DNA to be cut to the tube. Vortex DNA before pipetting to ensure that it is well-mixed.
4. Add 1µl of each enzyme
5. Place in a water bath at 37°C for 1hr.
6. Purify digested inserts using gel extraction.

5. Run gel agarose (1%)- 5ul of the sample.

Goal- Purification inserts using gel before ligation.

Gel to 5ul of reaction →

Purification with Kit → Check Concentration of product (nano drop).

store at -20C or continue to next part.

6. Gibson assembly on 3 parts (Terminator, 2 parts of gBlock)

Gibson assembly of the products according to the lab manual

Calculation:

$$DNA\ Con'_{\left[\frac{pmol}{\mu l}\right]} = \left(\frac{DNA\ Con'_{\left[\frac{ng}{\mu l}\right]} \times 1000}{Part\ size\ \times 650} \right)$$

Recommended Amount of Fragments Used for Assembly	
	2–3 Fragment Assembly
Total Amount of Fragments	0.02–0.2 <u>pmols</u> *X μl
Gibson Assembly Master Mix (2X)	10 μl
Deionized H ₂ O	10-X μl
Total Volume	20 μl***

control- vector only

Use 0.05pmol of plasmid and 0.15pmol of insert

Component	Volume
Plasmid (TT)	X
Insert (gBlock 1)	Y
Insert (gBlock 2)	Z

Gibson assembly mix	15µl
MBW	5-X-Y-Z
Total	20µl

Incubate samples in a thermocycler at 50°C for 60 minutes. Following incubation, store samples on ice or at –20°C for subsequent transformation.

7. Transformation

of the Gibson assembly products and control to TOP 10 (5ul)

8. plating on amp plates--> O/N

9. Verification

-PCR with primers to Mid of Tar only.

- Glycerol stock & miniprep

-Restriction of product with restriction enzymes

Materials:

- Restriction Enzymes Sall- HF, PstI- HF, clal)
- NEB buffer (Cutsmart Buffer)
- Molecular biology water

Procedure:

Set up the following reaction as described below:

For PctA-Tar

Component	Volume	Volume	Volume
Buffer (X10)	2µl	2µl	2µl
DNA up to 2µg	Xµl	Xµl	Xµl
PstI- HF	1µl	-	1µl
Sal I- HF	1µl	1µl	-
MB water	16-X µl	17-X µl	17-X µl
Total volume	20µl	20µl	20µl


[cloned PctA-Tar]= Y ng/ μ l
 1 μ g (cloned PctA-Tar)= μ l (~ μ l)
 2 μ g (cloned PctA-Tar)= μ l (~ μ l)
 → MB water= μ l

For NarX-Tar

Component	Volume	Volume	Volume
Buffer (X10)	2 μ l	2 μ l	2 μ l
DNA up to 2 μ g	X μ l	X μ l	X μ l
PstI-HF	1 μ l	-	1 μ l
Clal	1 μ l	1 μ l	-
MB water	16-X μ l	17-X μ l	17-X μ l
Total volume	20μl	20μl	20μl

[cloned NarX-Tar]= Y ng/ μ l
 1 μ g (cloned NarX-Tar)= μ l (~ μ l)
 2 μ g (cloned NarX-Tar)= μ l (~ μ l)
 → MB water= μ l

- Gel
Expected sizes:



MW: 1 Kb DNA Ladder

- ▼ 1: Cloned PctA+tar full plasmid
SalI + PstI
 - 1. 2546 bp
 - 2. 1644 bp
- ▼ 2: Cloned PctA+tar full plasmid
SalI
 - 1. 4190 bp
- ▼ 3: Cloned PctA+tar full plasmid
PstI
 - 1. 4190 bp
- ▼ 4: Cloned NarX+tar full plasmid
PstI + ClaI
 - 1. 2108 bp
 - 2. 1710 bp
- 5: Cloned NarX+tar full plasmid
ClaI
 - 1. 3818 bp
- 6: Cloned NarX+tar full plasmid
PstI
 - 1. 3818 bp

- Sending to sequence