

## WESTERN BLOT PROTOCOL

Goal:

Confirm that the protein of interest is expressed in the cell.

### SOLUTIONS

#### Lysis Buffer (for whole bacteria lysis)

65mM Tris pH 6.8 (7.4)->pH depends on the stacking gel

2% SDS

10% Glycerol

#### Loading Buffer

65mM Tris-HCl pH 6.8 (7.4)-> pH depends on the stacking gel

2% SDS

10% Glycerol

5% mercaptoethanol (add under the HOOD)

#### 1X Tris-Glycine Running Buffer

25mM Tris

190mM glycine

0.1% SDS

Adjust to pH 8.3 (if necessary or 7.4 depending on the gel)

#### 1X Transfer Buffer

25mM Tris

190mM glycine

20% Methanol

Adjust to pH 8.3 (if necessary); if protein is >80kDa then add 0.1% of SDS

#### PorcEAU Staining Buffer

0.2% (w/v) PorcEAU S

5% glacial acetic acid

#### 1X TBST WashingBuffer

20mM Tris

150mM NaCl

0.1% Tween-20

#### Blocking Buffer Milk

(5%) in TBST OR

BSA (3%) in TBST (if you will use a phospho-Antibody detection)

#### Antibody Solution

Antibody is diluted in Blocking Buffer (e.g. 1:5000)

### Stripping Buffer

20mL 10% SDS

12.5 mL 0.5M Tris-HCl

67.5mL MilliQ water

0.8mL  $\beta$ -mercaptoethanol

### PROCEDURE

#### 1. Whole Cell Lysis

- Grow bacteria to and OD600 about 0.5. Take 5 mL of the bacterial solution and centrifuge at 16,000 rpm for 5 min
- Remove the supernatant carefully
- Resuspend the pellet in 50  $\mu$ L of cold *Lysis Buffer*
- Boil the sample for 10 min at 95°C and freeze them at -20°C for 10min. Repeat this 3 times.

*Note:* You can aliquot the sample and store it at -20°C or -80°C for long-term storage

#### 2. Protein Separation by Gel electrophoresis

- Add to 10  $\mu$ L of the sample to 10  $\mu$ L of Loading Buffer. Add to this 1  $\mu$ L of bromophenol blue.
- Boil the samples at 95°C for 5 min.
- Quick-spin the sample down.
- Set-up the SDS-PAGE chamber by putting the pre-cast gel in the construct and then fills the tank with *Running Buffer*. Remove the comb in the tank.
- Load equal amounts (20  $\mu$ L) into the wells of a mini pre-cast gel (Biorad)
- Run the gel at 200V (max. mA) and let the blue marker run down shortly before it reaches the bottom (in the last third of the gel)- Approx 35-45min. Stop the run and get the gel out of it cast. (CAREFUL!!! Gel is very fragile! Don't break it!)
- Cut one corner to remember in which well you have loaded which sample.

#### 3. Transferring proteins from the gel to the membrane

- Let the gel stand in the Transfer buffer for 15 min.
- Prepare the transfer sandwich that consists of 2x sponge, then filter paper, then gel, then nitrocellulose membrane (also cut), then filter paper, then 2x sponges. All steps are done IN A TRANSFER BUFFER BATH.*Note:* Make clear to press out all bubbles and close the transfer cassette
- Put the transfer cassette in the transfer tank. Verify that the blot is on the

cathode and the membrane directed to the anode.

- Place also a block in the tank.
- Transfer the proteins at 100V for 30min (this step can immensely influence the bands that you will see later). Alternatively you transfer at 10mA during overnight in the cold room.
- 
- 
- 
- 

#### 4. Antibody incubations

- Briefly rinse the blot in water and stain it with Ponceau solution to check the transfer quality.
- Rinse off the Ponceau S stain with three washes with TBST.
- Block in 3% BSA in TBST (required for phopho--Antibodies) or 5% milk in TBST at room temperature for 1 hour.
- Wash three times with TBST.
- Cut the membrane into two pieces: one which will be incubated with the anti- body targeting the target protein; one against a control protein.
- Incubate overnight in the primary antibody solution against the target protein at 4°C.
- *Note:* As control you can use GAPDH-targeting antibody
- *Note:* The antibody should be diluted in the blocking buffer according to the manufacturer's recommended ratio. At room temperature: primery antibody needs about 1-3 hours depending on the antibody quality and performance.
- Recollect the antibody!
- Rinse the blot 3-5 times with TBST for 5 min.
- Incubate the secondary antibody (e.g. HRP-conjugated) for 1 hour at room temperature
- Recollect the antibody!
- Rinse the blot 3-5 times for 5 min with TBST

#### 5. Imaging & Data Analysis

- Apply the chemoluminescent substrate to the blot according to the manufacturers recommendation
- Capture the signal!
- Use the program imagej to read the band intensity of the target protein (in comparison to the loading control GAPDH)

#### 6. Stripping and reprobing

(if the protein sizes are quite different not necessary!!!)

Warm the Stripping Buffer to 50°C

- Add the Stripping Buffer to the membrane in a container designated for stripping. Incubate at 50°C for up to 45min with some agitation.
- Rinse the blot under running water for 1 hour
- Transfer the membrane to a clean container, wash 5 times for 5 min with TBST
- Block in 3% BSA or 5% Milk in TBST at room temperature for 1 hour.
- Incubate with the primary antibody overnight. (Recollect the antibody!)
- Rinse the blot for 3-5 times for 5 min with TBST.
- Incubate with secondary antibody (e.g. HRP-conjugated) for 1 hour at room temperature.
- Rinse the blot 3-5 times for 5 min with TBST.

## 7. Imaging & Data Analysis

(if the protein sizes are quite different not necessary!!!)

- Apply the chemiluminescent substrate to the blot following the manufacturer's suggestions.
- Capture the chemiluminescent signals
- Use imagej to read the band intensity of the loading control proteins.
- Use the loading proteins levels to normalize the target protein levels.