

6/3/2016 -- Gen. Notes

Friday, June 3, 2016

1:26 PM

General Plasmid Process

Build plasmid --> cloning strain amplification --> purify plasmid --> uptake plasmid into expression strain

Plasmid Purification Process

1. Spin/pellet the cells and remove supernatant
2. Lyse and extract cells
 - a. 250 uL of Buffer 1 (resuspends DNA) pipette into 1.5 mL tube
 - b. 250 uL of Buffer 2 (lyses cells) then invert tube
 - c. 350 uL of Buffer 3 (neutralizes cells) then invert tube
3. Spin and purify cells
 - a. Spin for 10 minutes (DNE = supernatant)
 - b. Load supernatant into column
 - i. Bind DNA
 - ii. Buffer A (inactivates enzymes from DNA degradation)
 - iii. Wash Buffer
 - iv. Elute with water or elution buffers

Cloning strains --> Top10 and DH5(alpha) (high transformation efficiencies)

Expression strains --> BL21 (need to overexpress proteins)

Expression and Genomic data --> MG1655

6/8/2016 -- Digests

Tuesday, June 14, 2016

3:39 PM

14A, 14C -- EcoRI, SpeI -- Cutsmart

-- EcoRI, XbaI, SAP -- Cutsmart

Note: Google NEB Double Digest Finder

General Protocol

250 ng DNA

2 uL Cutsmart

0.5 EcoRI

0.5 SpeI

--> water up to 20 uL

14A

7.12 uL DNA

8 uL CS
2 uL EcoRI
2 uL SpeI

60.88 uL water
14C

8.76 DNA
8 uL CS
2 uL EcoRI
2 uL SpeI

59.24 uL water
1P

7.00 uL DNA
8 uL CS
2 uL EcoRI
2 uL SpeI

61 uL water

Incubate 1 hour (37 degrees C)
Heat kill 20 min (80 degrees C)

8% gel, 100V, 40 min

6/9/2015 -- Enzyme Tests

Same protocol as 6/8

Needed to test enzyme due to abnormal results

Ran smaller digests

Tuesday, June 14, 2016

3:45 PM

<i>Plasmid</i>	14A	14C	1P
<i>DNA</i>	1.79	2.19	1.75
<i>Cutsmart</i>	2	2	2
<i>SpeI/XbaI</i>	0.5 SpeI	0.5 SpeI	0.5 XbaI
<i>EcoRI</i>	0.5	0.5	0.5
<i>Water</i>	15.22	14.81	15.25

Incubate 37 degrees C for 1 hour
Heart shock 80 degrees C for 15 min

6/14/2016 -- DNA Isolation


Tuesday, June 14, 2016

3:53 PM

Ran protocol from [UltraClean Microbial DNA Isolation Kit Instruction Manual](#)


Used detailed protocol

UltraClean[®] Microbial DNA Isolation



Screen clipping taken: 6/14/2016 3:56 PM

BBa_K1604010



Screen clipping taken: 6/29/2016 12:54 PM

7/11/2016 -- PCR (Primer+GFP)

Monday, July 11, 2016

1:56 PM

General PCR Protocol

<https://www.neb.com/protocols/1/01/01/taq-dna-polymerase-with-standard-taq-buffer-m0273>

- Setting the PCR thermocycler settings
 - Login as iGEM (password = igem)
 - Click on the iGEM folder, select a protocol (or make a new one), and adjust the settings
 - A temperature gradient can be created to cover a range of annealing temperatures
- Preparing PCR solution
 - Get a container and fill it with ice. PCR components should be mixed in test tubes on ice.
 - Mix the following concentrations according to the NEB PCR Kit
 - Ex:

Component	25 μ L Rxn	50 μ L Rxn
10 μ M Forward Primer	1.25 μ L	2.5 μ L
10 μ M Reverse Primer	1.25 μ L	2.5 μ L

Template DNA		
Q5 High Fidelity 2X Master Mix	12.5 µL	25 µL
Nuclease Free Water	to 25 µL	to 50 µL

Promoter PCR list for Viraat and I:

iGEM8	For Prom 2531	Add 308 uL for 100 uM
iGEM9	Rev Prom 2531	Add 279 uL for 100 uM
iGEM12	For Prom 1224	Add 305 uL for 100 uM
iGEM13	Rev Prom 1224	Add 288 uL for 100 uM
iGEM22	For Prom 3517	Add 276 uL for 100 uM
iGEM23	Rev Prom 3517	Add 288 uL for 100 uM
Template DNA		119.2 ng/uL Concentration
GFP Block		0.1 ng/uL
GFP_For		Add 279 uL for 100 uM
GFP_Rev		Add 292 uL for 100 uM

Centrifuge the 100 uM primers and resuspend/dilute with DI water in 1:10 ratio to get 10 uM.

Today's PCR Protocol

Component	25 uL rxn
GFP_For	1.25 uL
GFP_Rev	1.25 uL
2D Diluted GFP_Block (0.1 ng/uL)	
Master Mix	
DI Water	

Component	25 uL rxn
iGEM8 (F2531 Primer)	1.25 uL
iGEM9 (R2531 Primer)	1.25 uL
Template DNA (119.2 ng/uL)	1 uL
Master Mix	12.5 uL
DI Water	9 uL

Component	25 uL rxn
iGEM12 (F1224 Primer)	1.25 uL
iGEM13 (R1224 Primer)	1.25 uL
Template DNA	1 uL
Master Mix	12.5 uL
DI Water	9 uL
Component	25 uL rxn
iGEM22 (F3517 Primer)	1.25 uL
iGEM23 (R3517 Primer)	1.25 uL
Template DNA	1 uL
Master Mix	12.5 uL
DI Water	9 uL

<-- Need to add 1uL or less of GFP Block Template, but adding 1 uL at current concentration will yield 0.1 ng of GFP Block which seems too low??

Thermocycler:

Temp Gradient: 57-63 degrees C
Primers Tm: [59,61] degrees C