

Page No. \_\_\_\_\_

7/16 - IGEM

50- Pick colonies from proteins A(D), and B(D) and streaked out on new LB Kan 50 plates. Incubated overnight at 37°C

50- PCR reaction to amplify H out of Labrenzia CP44 ~~had~~ have not been able to amplify out of mixed (MC) community

H is heterostaggered, because it contains an internal cut site

PCR reaction:

H1	H2	Neg 1	Neg 2
25 $\mu$ L 2x master mix	25 $\mu$ L 2x master mix	25 $\mu$ L 2x master mix	25 $\mu$ L 2x master mix
25 $\mu$ M H Fwd (#71)	1 $\mu$ L 25 $\mu$ M H Fwd (#71)	1 $\mu$ L 25 $\mu$ M H Fwd (#71)	1 $\mu$ L 25 $\mu$ M H Fwd (#71)
25 $\mu$ M H Rev1 (#72)	1 $\mu$ L 25 $\mu$ M H Rev1 (#73)	1 $\mu$ L 25 $\mu$ M H Rev1 (#72)	1 $\mu$ L 25 $\mu$ M H Rev1 (#72)
23 $\mu$ L NF H <sub>2</sub> O	23 $\mu$ L NF H <sub>2</sub> O	23 $\mu$ L NF H <sub>2</sub> O	23 $\mu$ L NF H <sub>2</sub> O
Labrenzia CP colony	1 Labrenzia CP colony	1 Labrenzia CP colony	1 Labrenzia CP colony

amplifier Protocol:

Step	Temp	Time
Initial Denaturation	98°C	3 minutes
35 cycles	98°C	10 seconds
	70°C - 72°C	30 seconds
	72°C	45 seconds
Final Extension	72°C	2 minutes
Hold	10°C	50

Primers contained N's, so ~~their~~ their annealing temperatures cannot be determined (they range from 70°C to 72°C). H1 and H2 3 replicates of H1 and H2 were made, so each temp could be tested

Moved to -20°C when done

To Page No. \_\_\_\_\_

Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_


Date \_\_\_\_\_

TLE

om Page No. \_\_\_\_\_

8/18/16 - IGEM

2:00 - Ran 1% Agarose gel of H's PCR products at 110V for 40min.



analyzing F's transformations

Antibiotic Resistance	Volume Plate	~# of colonies	Notes
Kan	10ul	12	4 small, 8 large
Kan	100ul	~100	1/4 small
Kan	400ul	~1000	1/4 small
chloro	10ul	field	cannot select good isolated colonies
chloro	100ul		
chloro	400ul		

Will continue on with the 10ul ~~chloro~~ Kan plate.

Pick 8 colonies (4 small and 4 large) and streak onto one <sup>square</sup> grid of a 3x3 grid on a LB Kan50 plate. Put rest of colony left on plate into a PCR tube for colony PCR. Placed plate in 37C incubator.

Per PCR tube

12.5ul 2x Green Mix  
0.5ul Fwd1 (Primer # 68)  
0.5ul FRev (Primer # 70)  
11.5ul NF H<sub>2</sub>O  
1 colony

Thermocycler

1) 95°C for 3:00  
2) 95°C for 0:30  
3) 67°C for 0:30  
4) 72°C for 3:00  
5) 2-4 x 34 hms  
6) 72°C for 5:00  
7) 10°C oo

Drawing of 3x3 Grid



Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_

Date \_\_\_\_\_


Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

Page No. \_\_\_\_\_

18/16 - I GEM

TITLE \_\_\_\_\_



← Gel from PCR amplification of  
PB1 H from *Laobrenzia* CP4  
No PCR products visible  
↳ 4th attempt at PCR w/o  
success → will analyze primers  
again

00 - <sup>03</sup> PCR check of Gene A (D1 and B D2) and B (D3)

Picked 1/2 of a colony w/ P10 pipette tip from plates streaked on  
OB/17/16 and used for colony PCR (outlined below)

Picked second half of colony w/ wooden stick and used grew up  
in 5mL LB Kan for H-G in 37°C shaker

Colony PCR:  
ran a negative control of each without a colony

Gene A	Gene B
12.5 µL 2x Green Mix	12.5 µL 2x Green Mix
0.5 µL A Fwd (Primer G2)	0.5 µL B Fwd (Primer G5)
0.5 µL A Rev (Primer G4)	0.5 µL B Rev (Primer G6)
11.5 µL NF H <sub>2</sub> O	11.5 µL NF H <sub>2</sub> O
1 colony	1 colony

Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_

Date \_\_\_\_\_

To Page No. \_\_\_\_\_

From Page No. \_\_\_\_\_

08/18/16 - I GEM

Thermocycler for Genes A and B PCR

STEP	TEMP	TIME
Initial Denaturation	95°C	3:00
35 cycles	95°C	0:30
	66°C	0:30
	72°C	1:00
Final Extension	72°C	2:00
Hold	10°C	∞

Ran 1% Agarose Gel at 110V for 20min of A and B's PCR and F's PCR

-A		AD1	2-105 100bp Ladder	AD2	BD3	Ladder	-B		
FA1	FA2	FA3	Ladder	FB1	FB2	Ladder	FB3	FC1	FC2
10 $\mu$ L			7 $\mu$ L	10 $\mu$ L		7 $\mu$ L	10 $\mu$ L		



More likely that something went wrong with the PCR runs than all 3 colonies did not have the insertions

15:30 - Split each liquid culture (AD1, AD2) into 3 15 mL tubes  
 3 mL LB  
 3 mL Kan (50 mg/mL)  
 30 $\mu$ L culture  
 grow up overnight in 37°C shaker

Verified by: \_\_\_\_\_

Date \_\_\_\_\_

Project No. \_\_\_\_\_  
Book No. \_\_\_\_\_ TITLE \_\_\_\_\_

Page No. \_\_\_\_\_  
9/16 - I GEM

Freezer stocks (800  $\mu$ L cells, 250  $\mu$ L 50% glycerol) and dot  
prepped remaining with overnight growths of AD1, AD2, BD3  
eluted w/ 60  $\mu$ L Buffer EB

Stock #	Insert	Plasmid	Label
89	A	Kan	AD1
90	A	Kan	AD2
91	B	Kan	BD3

Not dropped

Sample ID	User name	Date and Time	Nucleic Acid Conc.	Unit	A260	A280	260/280	260/230
A D1	SS Gabrielle Tender	8/19/2016 11:43:07 AM	328.1	ng/ $\mu$ L	6.561	3.460	1.90	2.29
A D2	SS Gabrielle Tender	8/19/2016 11:44:01 AM	285.5	ng/ $\mu$ L	5.709	3.007	1.90	2.29
B D3	SS Gabrielle Tender	8/19/2016 11:44:48 AM	516.2	ng/ $\mu$ L	10.324	5.462	1.89	2.24
1MCS-2::Undecaheme D1	SS Gabrielle Tender	8/19/2016 11:45:31 AM	94.3	ng/ $\mu$ L	1.886	0.971	1.94	2.17
1MCS-2::Undecaheme D2	SS Gabrielle Tender	8/19/2016 11:46:09 AM	82.6	ng/ $\mu$ L	1.653	0.854	1.93	2.18
EB	SS Gabrielle Tender	8/19/2016 11:46:44 AM	-0.3	ng/ $\mu$ L	-0.006	-0.005	1.14	0.35

Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_

Date \_\_\_\_\_

To Page No. \_\_\_\_\_

TITLE		Project No.	Book No.
From Page No.			
08/22/16 - IGEN!			
Streaked FB3 and FCI for singles on a LB Kan50 plate Incubated at 37°C overnight			
08/23/16 - IGEN!			
Picked 1/2 of a colony from F plates (one from FB3 and one from FCI) w/ a p10 pipette tip. Colony PCR (outlined below) to verify insert. Picked second 1/2 of colony w/ wooden stick and grew up in 3mL of LB Kan50 in 37°C shaker for 4-6 hrs			
PCR		3 PCR tubes	"master mix"
12.5 µL 2x Green Mix		• neg (no colony)	37.5 µL 2x Green Mix
0.5 µL Fwd Primer (#60)		• FB3 (FB3 colony)	1.5 µL Fwd
0.5 µL FRev Primer (#70)		• FCI (FCI colony)	1.5 µL FRev
14.5 µL MFH20			34.5 µL MFH20
1 colony			25 µL per tube
Thermocycler			
1) 95°C for 3:00			
2) 95°C for 0:30			
3) 67°C for 0:30			
4) 72°C for 3:00			
5) 2~4 x 34 times			
6) 72°C for 5:00			
7) 10°C for ∞			
Digestion for tomorrow's ligation (small scale to verify correct digestion) Plan is to dig B w/ EcoRI and SpeI, Fw/XbaI and PstI, and A w/ EcoRI and PstI Ligate together B and F into plasmid from A (or find linearized backbone)			
AD1		AD2	BD3
0.5 µL EcoRI	0.5 µL EcoRI	0.5 µL PstI/EcoRI	37°C for 1 hr
0.5 µL PstI	0.5 µL PstI	0.5 µL PstI/SpeI	
1 µL AD1	1 µL AD2	0.5 µL BD3	
1 µL 10x usimant	1 µL 10x usimant	1 µL 10x usimant	
Recorded by: Tal MFH20	7 µL MFH20	Verified by: 7.5 µL MFH20	Date





Project No. \_\_\_\_\_

TITLE \_\_\_\_\_


Book No. \_\_\_\_\_

From Page No. \_\_\_\_\_

08/33/16-IGEM!

Run gel of PCR and RE digests at 110V for 40min

10kb	neg	FB3	FC1	10kb	AD1	AD2	BOS	10kb
7 $\mu$ L	PCR rxn			7 $\mu$ L	RE digest			7 $\mu$ L
	5 $\mu$ L				10 $\mu$ L + 2 $\mu$ L Gx dye			



Split FB3 and FC1 liquid cultures into 3 liquid cultures each  
3mL LB  
3mL Kan 50mg/mL  
30mL previous liquid culture  
incubate in 37°C shaker overnight

To Page No. \_\_\_\_\_

Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_

Date \_\_\_\_\_

Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

Page No. \_\_\_\_\_

24/10 - iGEM

de freezer stocks (800  $\mu$ L cells, 250  $\mu$ L 50% glycerol)

Freezer stocks #

DNA

92

iGEM F- B3

93

iGEM F- C1

nipped DNA (eluted in 50  $\mu$ L Buffer EB) and nanodropped

Sample ID	User name	Date and Time	Nucleic Acid Conc.	Unit	A260	A280	260/280	260/230
EM FB3	SS Gabrielle Tender	8/24/2016 12:23:08 PM	166.4	ng/ $\mu$ L	3.327	1.752	1.90	2.28
EM FC1	SS Gabrielle Tender	8/24/2016 12:24:03 PM	247.9	ng/ $\mu$ L	4.957	2.613	1.90	2.19

## Restriction Enzyme Digest

DB3	FB3	FC1	2 $\mu$ L BD3	Chloro Vector
BD3	9 $\mu$ L FB3	4 $\mu$ L FC1	Kan vector	5 $\mu$ L pSB1C3
Cutsmart Buffer	5 $\mu$ L cutsmart	5 $\mu$ L cutsmart	5 $\mu$ L cutsmart	5 $\mu$ L cutsmart
EcoRI	1 $\mu$ L EcoRI	1 $\mu$ L EcoRI	1 $\mu$ L EcoRI	1 $\mu$ L EcoRI
SpeI	1 $\mu$ L SpeI	1 $\mu$ L SpeI	1 $\mu$ L SpeI	1 $\mu$ L SpeI
NFHa0	28 $\mu$ L NFHa0	40 $\mu$ L NFHa0	42 $\mu$ L NFHa0	29 $\mu$ L NFHa0

37  $^{\circ}$ C for 1 hr

ed to separate backbone from insert

on 1% Agarose gel of all but chloro backbone (doesn't contain insert)  
+ 110V for 40 min

BD3	2 $\mu$ L BD3	FB3	FC1	2 $\mu$ L BD3	Kan vector
2 $\mu$ L	2 $\mu$ L	50 $\mu$ L	20 $\mu$ L	50 $\mu$ L	10 $\mu$ L
10 $\mu$ L Gdye		10 $\mu$ L Gdye			10 $\mu$ L Gdye

To Page No. \_\_\_\_\_

Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_

Date \_\_\_\_\_



Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

From Page No. \_\_\_\_\_

08/24/16 - iGEM

PCR cleaned up Clono vector (eluted w/ 50 $\mu$ L Buffer EB)

Gel purified BD3 and Kpn vector (eluted w/ 50 $\mu$ L Buffer EB)

No insert in either F vector, both

picked 5 colonies from each F<sup>+</sup> B3 and F<sup>+</sup> C1 plates

Grew up in 3 mL LB Kana 50

08/25/16 - iGEM

Made freezer stocks of all 10 F<sup>+</sup> growths (FB3A, FB3B, ..., FC1D, FC1E) (800 $\mu$ cells, 25 $\mu$ L)

Labeled each one 99, but then put a label w/ tape on side

Miniprep DNA (eluted w/ 25 $\mu$ L Buffer EB)

↳ accidentally labeled FC1C as FC1D, so can't differentiate - relabeled the two

"FC1D" as FC1C/D and FC1D/L

Nanodropped

Sample ID	User name	Date and Time	Nucleic Acid Conc. Unit	A260	A280	260/280
psB1K3 digested and purified	SS Gabrielle Tender	8/25/2016 10:43:37 AM	10.5	ng/ $\mu$ L	0.210	0.092 2.29
Buffer EB	SS Gabrielle Tender	8/25/2016 10:44:25 AM	0.4	ng/ $\mu$ L	0.008	-0.010 -0.79
psB1K3 digested and purified	SS Gabrielle Tender	8/25/2016 10:45:16 AM	27.67	ng/ $\mu$ L	0.551	0.285 1.93
iGEM BD3 digested and purified	SS Gabrielle Tender	8/25/2016 10:46:13 AM	5.7	ng/ $\mu$ L	0.113	0.035 3.21
psB1K3 digested and purified	SS Gabrielle Tender	8/25/2016 10:47:05 AM	3.6	ng/ $\mu$ L	0.072	0.009 8.04
iGEM FB3 A	SS Gabrielle Tender	8/25/2016 10:48:03 AM	152.5	ng/ $\mu$ L	3.050	1.573 1.94
iGEM FB3 B	SS Gabrielle Tender	8/25/2016 10:48:51 AM	149.7	ng/ $\mu$ L	2.994	1.550 1.93
iGEM FB3 C	SS Gabrielle Tender	8/25/2016 10:49:42 AM	151.3	ng/ $\mu$ L	3.027	1.566 1.93
iGEM FB3 D	SS Gabrielle Tender	8/25/2016 10:50:28 AM	125.9	ng/ $\mu$ L	2.519	1.286 1.96
iGEM FB3 E	SS Gabrielle Tender	8/25/2016 10:51:29 AM	129.7	ng/ $\mu$ L	2.594	1.340 1.94
1 iGEM FC1 A	SS Gabrielle Tender	8/25/2016 10:52:22 AM	107.3	ng/ $\mu$ L	2.145	1.096 1.96
2 iGEM FC1 B	SS Gabrielle Tender	8/25/2016 10:53:22 AM	115.5	ng/ $\mu$ L	2.310	1.186 1.95
3 iGEM FC1 C/D	SS Gabrielle Tender	8/25/2016 10:54:20 AM	122.5	ng/ $\mu$ L	2.450	1.267 1.93
4 iGEM FC1 D/C	SS Gabrielle Tender	8/25/2016 10:57:57 AM	103.1	ng/ $\mu$ L	2.061	1.037 1.99
5 iGEM FC1 E	SS Gabrielle Tender	8/25/2016 10:58:54 AM	119.6	ng/ $\mu$ L	2.393	1.222 1.96
6 Buffer EB	SS Gabrielle Tender	8/25/2016 10:59:47 AM	0.5	ng/ $\mu$ L	0.011	-0.005 -2.29

Recorded by: \_\_\_\_\_

Date: \_\_\_\_\_

Verified by: \_\_\_\_\_

Date: \_\_\_\_\_

To Page No. \_\_\_\_\_

Project no. \_\_\_\_\_ TITLE \_\_\_\_\_  
 Book No. \_\_\_\_\_  
 Form Page No. \_\_\_\_\_

B125/Kc - iGEM1

~~incubated~~  
 digested each F vector w/ XbaI and PstI (BstNotw  
 only did 10uL digestion (didn't want to waste enzyme or Buffer)

Per tube

2uL DNA		"Master Mix"
1.5uL XbaI		5uL XbaI
0.5uL PstI		5uL PstI
1uL CutSmart Buffer		10uL CutSmart
7.5uL NF H <sub>2</sub> O		75uL NF H <sub>2</sub> O
		9uL per tube

37°C for 1hr

Found digested and undigested plasmid backbones.  
 NanoDrop and image to determine length, purity and conc.

Sample ID	User name	Date and Time	Nucleic Acid Conc.	Unit	A260	A280	260/280
GEM clean cut plasmid K 7/12/16	SS Gabrielle Tender	8/25/2016 12:04:34 PM	15.1	ng/ul	0.303	0.176	1.72
digested plasmid C	SS Gabrielle Tender	8/25/2016 12:05:45 PM	26.2	ng/ul	0.523	0.305	1.72
GEM K plasmid 7/11/16 KC	SS Gabrielle Tender	8/25/2016 12:06:48 PM	288.1	ng/ul	5.762	3.033	1.90

Ran 1% Agarose Gel at 110V for 20min

FB3 A	FB3 B	FB3 C	FB3 D	FB3 E	2log 10kb	GEM A-Vector 200kb	clean cut plasmid 372kb	plasmid clean 100kb	2log 10kb
FC1 A	FC1 B	FC1 C	FC1 D	FC1 E	2log 10kb	digested plasmid C	very dirty digest without boiler	plasmid from 100kb	2log 10kb

RE digest: 10uL + 2uL Qsdye

7uL

2uL + 4uL NF H<sub>2</sub>O + 1uL Gage

Recorded by: \_\_\_\_\_ Date \_\_\_\_\_ Verified by: \_\_\_\_\_ Date \_\_\_\_\_ To Page No. \_\_\_\_\_


||||

Project No. \_\_\_\_\_  
Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

From Page No. \_\_\_\_\_

OB 12/14 - IGEM



Sol'n: amplify F out of MCL extract, FCID, FCIA. amplify B out of BD3.

F - MCL

7.  $\mu$ L 350 extract  
1.  $\mu$ L 250M Fwd  
1.  $\mu$ L 250M Rev  
30  $\mu$ L BD3  
1.  $\mu$ L NT 120

F - FCID

7.  $\mu$ L 350 extract  
1.  $\mu$ L 250M Fwd  
1.  $\mu$ L 250M Rev  
30  $\mu$ L BD3  
1.  $\mu$ L NT 120

F - FCIA

7.  $\mu$ L 350 extract  
1.  $\mu$ L 250M Fwd  
1.  $\mu$ L 250M Rev  
30  $\mu$ L BD3  
1.  $\mu$ L NT 120

B - BD3

7.  $\mu$ L 350 extract  
1.  $\mu$ L 250M Fwd  
1.  $\mu$ L 250M Rev  
30  $\mu$ L BD3  
1.  $\mu$ L NT 120

Made 1:10 dilutions of FCID, FCIA, BD3 (use these in PCR)  
Made 1:4 dilutions of 350 extract #5 and 6 Rev primer #66  
Master mix

To Page No. \_\_\_\_\_

Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_

Date \_\_\_\_\_

From Page No.

08/22/14 - T-GEM

PCR from pg 71

F master mix

1  
 0.75  $\mu$ L 100  $\mu$ M FFH1 (primer #68)  
 0.75  $\mu$ L 100  $\mu$ M FFH2 (primer #70)  
 0.75  $\mu$ L Q5 2x master mix  
 4.00  $\mu$ L NF H<sub>2</sub>O

B master mix

1  $\mu$ L 100  $\mu$ M BFwd (primer #65)  
 1  $\mu$ L 100  $\mu$ M BRev (primer #66)  
 6  $\mu$ L NF H<sub>2</sub>O

F-MCL

1  $\mu$ L MCL SE extract  
 49  $\mu$ L F master mix

F-FC1D

1  $\mu$ L 1:10 FC1 D/C  
 49  $\mu$ L F master mix

F-FC1A

1  $\mu$ L 1:10 FC1A  
 49  $\mu$ L F master mix

B BD3

1  $\mu$ L 1:10 BD3  
 2  $\mu$ L B master mix  
 25  $\mu$ L Q5 2x master  
 22  $\mu$ L NF H<sub>2</sub>O

F neg

1  $\mu$ L NF H<sub>2</sub>O  
 49  $\mu$ L F master mix

Thermocycler:

STEP	TEMP	TIME
Initial Denaturation	98°C	30 sec
35 cycles	98°C	10 sec
	72°C	30 sec
	72°C	1 min
Final Extension	72°C	2 min
Hold	10°C	00

Recorded by:

Date

Verified by:

Date



Project No. \_\_\_\_\_

Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

From Page No. \_\_\_\_\_

DB/26/Ka-igEM

I do not trust any of the PCR'd or digested Kan vector (psBK3.m1) left for me. I will instead digest it out of psBK3.m1::BD3

~~digestion~~  
8/26/11

Run 1% agarose gel of PCR products




→ Digestion of psBK3.m1::BD3 to recover psBK3.m1

Did 4 50 $\mu$ L reactions

Per tube (did 4x master mix and split it into 4 tubes)

1 $\mu$ L EcoRI-HF1 $\mu$ L PstI2 $\mu$ L igEM BD35 $\mu$ L 10X cutsmart Buffer40 $\mu$ L NF H<sub>2</sub>O

4x master mix

4 $\mu$ L EcoRI4 $\mu$ L PstI8 $\mu$ L igEM BD320 $\mu$ L 10X cutsmart104 $\mu$ L NF H<sub>2</sub>O37 $^{\circ}$ C incubation for one hour

I am very confused, because it amplified out of F-PCD but not MCL ~~gen~~ genomic DNA

PCR clean up F-PCD

Labeled ~~kan~~ "igEM F"

Gel purified B-BD3 by cutting out 1 kb band and purifying (look on next page for more info)

Recorded by: \_\_\_\_\_

Date: \_\_\_\_\_

Verified by: \_\_\_\_\_

Date: \_\_\_\_\_

To Page No. \_\_\_\_\_

Project No. \_\_\_\_\_  
Book No. \_\_\_\_\_

TITLE \_\_\_\_\_

Page No. \_\_\_\_\_

26/10-IGEM

purified pSB1K3+ and gene B 600  
through 1% agarose gel for 40 min  
digested pSB1K3+ : 1.005  
200 1000  
15 µL  
50 µL 100 µL  
edge  
← stopped loading b/c gel leaky

was all squishy...

it amplified B, will do gel purification of digested product

2x50 µL rxns:

rxn tube:

1 µL EcoRI

1 µL SpeI

1 µL IGEM BD3

1 µL 2x cutsmart

1 µL NPH20

200 µL Master Mix

2 µL EcoRI

1 µL SpeI

2 µL IGEM BD3

5 µL 2x cutsmart

42 µL NPH20

C for 1 hr

cleaned up: F from MCL and F from FCIA too (added w/ 50 µL Buffer BB cut)

Sample ID	User name	Date and Time	Nucleic Acid Conc.	Unit	A260	A280	260/280	260/230
F PCRd from FCID/C	SS Gabrielle Tender	8/26/2016 2:35:21 PM	8.6	ng/µl	0.173	0.083	2.08	2.68
F PCRd from MCL	SS Gabrielle Tender	8/26/2016 2:36:21 PM	6.4	ng/µl	0.128	0.053	2.42	2.26
F PCRd from FCIA	SS Gabrielle Tender	8/26/2016 2:37:06 PM	6.5	ng/µl	0.130	0.055	2.38	2.04
err EB	SS Gabrielle Tender	8/26/2016 2:37:53 PM	-1.0	ng/µl	-0.021	-0.015	1.37	0.35

Recorded by: \_\_\_\_\_

Date \_\_\_\_\_

Verified by: \_\_\_\_\_

Date \_\_\_\_\_

To Page No. \_\_\_\_\_