

Week 4

June 27, 2016

What we did today:

PCR

Liquid cultures

Drafted Email for Human Practices

Transformed RFP, InterLab Device 3, pRETRO-ON plasmids

Streaked GFP from Glycerol Stock

Made glycerol stock of pRETRO-ON

PCR:

Mixed two tubes of PCR reaction:

25 μ L 2xPCR mastermix

1.25 μ L dCas9 BamH1 primer

1.25 μ L dCas9 Pac1 primer

22.5 μ L DI water

1 picked colony of pdCas9

PCR protocol:

2 minutes at 95°C

30 seconds at 95°C

45 seconds at 50°C


4 minutes 30 seconds at 72°C

----- repeat 30 times-----

10 minutes at 72°C

Held at 10°C until taken out of PCR

Ran on Gel:



Cut out two pieces of gel circled above
Weighed each gel piece
2 μL of NT for every 1 μg of reaction
Allowed gel pieces to melt in 42°C water bath
Added to gel extraction columns
Spun for 1 minute
Dumped out flow through
Added 700 μL of NT3
Spun for 1 minute
Dumped out flow through
Spun for 2 minutes
Dumped out flow through
Put column into microfuge tube
Added 50 μL of NE
RT for 1 minute
Spun for 1 minute

Restriction:
50 μL of PCR extracted
6 μL cutsmart
2 μL Pac1
2 μL BamH1
Incubated in 37°C water bath over night

Transformations into DH10Beta cells

- GFP
- RFP
- pRetro-On
 - Made Glycerol Stocks of Each

PCR

PCR of dCas9

dCas9 grew better on 2x amp plates

Double digest of Pac1 and BamH1 on plasmid

10ug plasmid

5uL Cutsmart

2uL each enzyme

31uL water (50uL total)

Liquid cultures 3 tubes dCas9 in amp and LB

Gel wide wells (.6g agarose/60mL 1x TAE)

50uL sample + 10uL load dye + SYBR green

Lane 1: ladder

Lane 2: 30uL sample

Lane 3: 30uL sample

Cut out band with razor blade, place in tube and do tube kit

- Plasmid prep: pRetro
 - Total = 50 μ L
 - 5 μ L CutSmart buffer
 - 10 μ g pRetro
 - 2 μ L each enzyme
 - 1 μ L SAP after 2 hours
- Liquid culture:
 - dCas9-humanized, 2X Amp LB (10 μ L/tube), 3 tubes
- Gel
 - pRetro: BamHI, PacI digest
- DNA purification from gel
 - 350 mg gel/100 mg * 200 μ L NTI
 - Water bath to melt
 - Spin, add 600 μ L NT3
 - 50 μ L elution buffer

June 28, 2016

What We Did Today:

Transformed pRETRO-ON

Transformed RFP from iGEM 2016 competency pack

Restreaked ptet-RFP, pLac-RFP, and GFP

PCR dCas9

Emailed Mr. Sontheimer and Mr. Zhang about integrated human practices

(Sontheimer will not be in the office until July 18, email again around the 20?)

Heated each PCR reaction in 80°C water bath for 20 minutes

PCR Clean-up:

100 µL of NT and 50 µL of reaction into PCR clean-up columns

Spun for 1 minute

Dumped out flow through

Added 700 µL of NT3

Spun for 1 minute

Dumped out flow through

Spun for 2 minutes

Dumped out flow through

Put column into microfuge tube

Added 50 µL of NE

RT for 1 minute

Spun for 1 minute

Ligation:

Created 4 separate ligation reactions

1. 1 µL vector (pRetro-ON), 1 µL Insert (dCas9), 15 µL sterile water, 2 µL 10x ligation buffer, and 1 µL T4 ligase
2. 1 µL vector, 2.5 µL Insert, 13.5 µL sterile water, 2 µL 10x ligation buffer, and 1 µL T4 ligase
3. 1 µL vector, 7 µL Insert, 9 µL sterile water, 2 µL 10x ligation buffer, and 1 µL T4 ligase
4. 1 µL vector, 0 µL Insert, 16 µL sterile water, 2 µL 10x ligation buffer, and 1 µL T4 ligase (negative control)

Incubated at RT for 30 minutes

Transformed 5 µL of each ligation with 50 µL of DH10β cells

Plated onto 1xAMP plates

None grew

PCR


Using dCas9 mini-preps

3 mini-preps:

1.) 1.002 (diluted 1:8)

2.) 0.897 (diluted 1:8)

3.) 0.078



PCR Gel

- Lane 1: Log2 Ladder
- Lane 2: Empty
- Lane 3: Mini-prep 1
- Lane 4: Mini-prep 2
- Lane 5: Mini-prep 3

Analysis: All of the bands for the three min-preps are at the same spot in the gel around the 6kb mark.

Cut out three pieces of gel where the bands were located above

Weighed each gel piece

- 1.) 0.57g -> 570mg -> 1140 μ L
- 2.) 0.46g -> 460mg -> 920 μ L
- 3.) 0.81g -> 810mg -> 1620 μ L

2 μ L of NT for every 1 μ g of reaction

Allowed gel pieces to melt in 42°C water bath

Added to gel extraction columns

Spun for 1 minute

Dumped out flow through

Added 700 μ L of NT3

Spun for 1 minute

Dumped out flow through

Spun for 2 minutes

Dumped out flow through
Put column into microfuge tube
Added 50 μ L of NE
RT for 1 minute
Spun for 1 minute

Restriction:
50 μ L of PCR extracted
6 μ L cutsmart
2 μ L Pac1
2 μ L BamH1
Incubated in 37°C water bath over night

Made glycerol stock of dCas9, in the #5 slot
Miniprep dCas9 humanized

Digest confirmation of pgRNA
Digested with Xho1, BamH1, EcoR1, Not1, and double digest of Not1 and EcoR1, and negative control

Ran gel
Lanes:
1: Ladder
2: Negative Control
3: Xho1
4: BamH1
5: EcoR1
6: Not1
7: Not1 EcoR1 double digest

Ligation of dCas9 into pRetro-ON
-4 types: 1 μ L insert, 2.5 μ L insert, and 7 μ L insert
Incubated at RT for 30 min
Transformed into DH10Beta, (5 μ L ligation + 50 μ L cells), plated on 1x AMP

- Miniprep: pgRNA-humanized (3 samples)
 - 1: 100.7 ng
 - 2: 107.7 ng
 - 3: 833.4 ng
- Digest & Gel
 - pgRNA-humanized (used 833.4 ng stock, 1.2 μ L/sample)
 - Lanes (L to R): Ladder, Negative, XhoI, BamHI, ECoRI, NotI, ECoRI+NotI
 - 12 μ L/well, 4 μ L dye/tube

June 29, 2016

What We Did Today:

PCR Cleanup

Transformed ligations from digest (using both DH5(alpha) cells and DH10(Beta) cells)

Made liquid cultures of RFP and GFP

Streaked and made liquid culture of new DH5(alpha) cells

Results from yesterday's transformations:

RFP from cell competency pack worked

PRETRO-ON was overgrown (possible problem with competent cells)

Restreaks all worked

PCR Clean-up

100 μ L of NT and 50 μ L of reaction into PCR clean-up columns

Spun for 1 minute

Dumped out flow through

Added 700 μ L of NT3

Spun for 1 minute

Dumped out flow through

Spun for 2 minutes

Dumped out flow through

Put column into microfuge tube

Added 50 μ L of NE

RT for 1 minute

Spun for 1 minute

Ligation:

Using protocol

Made different dilution (neg control, 1 μ L, 2.5 μ L, and 7 μ L)

Transformed using both Dh5(alpha) and DH10(beta) cells

Made liquid cultures of P-Tag RFP

Transformation of ligation reaction into EGM2K lambda cells.

7.5 μ L of reaction into 50 μ L cells, shaker for 2 hours plated on 1x AMP plates

New insert

-80 C 20min and let cool to rt,

PCR clean up. Ligation

Transformed into DH10Beta and DH5Alpha cells

- Transformation

Cells: DH5 α , EMG2K λ

DNA: ligation products (pRetro-dCas9), 1, 2.5, 7 μ L samples

- 7.5 μL DNA/tube
- DH5 α in 1X Amp, EMG2K λ in 2X Amp
- PCR set up
 - dCas9, 1.002 $\mu\text{g}/\mu\text{L}$ dilute 1:8
 - 22.5 μL dH₂O
 - 1 μL DNA
 - 1.25 (1.3) μL 3' (PacI) & 5' (BamHI) primer each
 - 25 μL 2X PCR Mastermix
 - PCR cycle * 30

June 30, 2016

What we did today:

Checked cultures – RFP did not glow red, GFP did

Unable to make glycerol stocks because didn't glow red **(USE CHLOR NEXT TIME NOT AMP!!)**

Results from Transformations:

Transformations using DH5(alpha) cells were very overgrown

Transformations using DH10(beta) barely grew

DH5a Competent Cell Protocol:

- 1.) 2ml of the liquid culture into flask with 200mL of LB (no antibiotics)
- 2.) Shaker to grow for 1.5hrs
- 3.) Incubate until OD_{600nm} of 0.3. Continue incubation until OD reading is correct
- 4.) Pre-chill 4 50 mL conical tubes and large centrifuge to 4 degrees.
- 5.) Transfer 20mL of cells to each tube
- 6.) Incubate on ice for 20 mins.
- 7.) Combine two tubes into one (results in two tubes with 40mL each)
- 8.) Centrifuge at 4 degrees at 3000g for 15 Mins
- 9.) Pour out supernatant into waste.
- 10.) Resuspend pellet in 5mL of ice cold CCMB80 buffer
- 11.) Incubate for 20 mins on ice
- 12.) Put 200µL into chilled microcentrifuge tubes
- 13.) Store at -80
- 14.) Test Competency by transforming with GFP.

PCR run in gel

5uL PCR rxn + 1uL dye

Lane 1: Ladder

Lane 2: -

Lane 3: PCR rxn 1

Lane 4: -

Lane 5: PCR rxn 2

Lane 6: -

Lane 7: PCR rxn 3

Lane 8: -

Lane 9: PCR rxn 4

All had bands except #4, 1 the brightest

45uL PCR cleanup

90uL NT1

spun in yellow column

add 700uL NT3 spin 1 min

dump spin 2 min

add 50uL elution buffer

digestion 50uL of above
6uL smartcut buffer
2uL Pac1
2uL BamH1
digested 1.5 hours
Heat kill 20 mins 80 degrees
Ran another gel
#1 & 3 had bands

Ligation:
EMG2k lambda cells 2xamp plates

Insert: 0, 2, 7, 15uL
water: 15, 13, 8, 0uL
2uL vector
2uL ligase buffer
1uL ligase

sits room temp 30 mins
Transform 5uL into 50uL cells
Transformation protocol
Wait full 2 hours on the shaker

APOBEC1 PCR

- Gel
 - PCR product (dCas9)
 - 5 μ L/well + 1 μ L dye
 - Lanes (L to R): Ladder, 1, 2, 3, 4
- PCR clean up
 - PCR product (dCas9)
 - ~45 μ L PCR + 200 μ L NT1/100 μ L PCR = 90 μ L NT1
 - 700 μ L NT3
 - 50 μ L elution buffer
- Digest
 - PacI & BamHI
 - 60 μ L final volume
- Gel
 - Digested PCR clean up product (dCas9)
 - 5 μ L/well + 1 μ L dye
 - Lanes (L to R): Ladder, 1, 2, 3, 4
- Heat kill enzyme in the rest of the digest (80 °C, 20 min)
- Ligation
 - Vector: pRetro
 - Insert: PCR product (dCas9)
 - Water: 15, 13, 8, 0 μ L
 - Insert: 0, 2, 7, 15 μ L

- 2 μ L buffer
- 2 μ L vector
- 1 μ L ligase
- 30 min in RT
- Transformation
 - Cells: EMG2K λ
 - 50 μ L cell
 - 5 μ L ligation product
 - 2X Amp plates