

# Week 9

August 1, 2016

Gibson Plate Counts




Plate cells today for transfection Tuesday with Apobec 2x and 3x

Wiki day

For the test digests of last week


The only two enzymes that gave a small amount of background were BamHI and HindIII

HEK293T cell maintenance

Created 2 12 well plates for fluorescence test and western blot

- Data analysis
  - eGFP trial 1 slides
    - ImageJ
    - Excel

August 2, 2016



Gel 1

Lane 1 – Ladder

Lane 2 – Adar 1-1x #1

Lane 3 – Adar 1-1x #2

Lane 4 – Adar 1-1x #3

Lane 5 – Empty

Lane 6 – Adar 1-2x #1

Lane 7 – Adar 1-2x #2

Lane 8 – Adar 1-2x #3

Analysis

They ran too small.



Gel 2

Lane 1 – Ladder

Lane 2 – Empty

Lane 3 – Adar 1-3x #1

Lane 4 – Adar 1-3x #2

Lane 5 – Empty

Lane 6 – Adar 2-1x #1


Lane 7 – Adar 2-1x #2

Lane 8 – Adar 2-1x #3

Analysis

The gels ran too small.

Especially something is wrong in Lane 8.




### Gel 3

Lane 1 – Ladder  
 Lane 2 – Adar 2-2x #1  
 Lane 3 – Adar 2-2x #2  
 Lane 4 – Adar 2-2x #3  
 Lane 5 – Empty  
 Lane 6 – Adar 2-3x #1  
 Lane 7 – Adar 2-3x #2  
 Lane 8 – Adar 2-3x #3

### Analysis

Lanes ran too small.



### Gel 4

Lane 1 – Ladder  
 Lane 2 – Apobec 1x #1  
 Lane 2 – Apobec 1x #2  
 Lane 3 – Apobec 1x #3

### Analysis


Lanes ran too small.



- New cell culture media  
100 mL
- Cell plating  
Cell count: 327; 395 → Average = 361
  - $1.444 \times 10^6$  cells/mL
 Plating:
  - 12-well plate:  $3 \times 10^5$  cells \* 13 wells → 2.7 mL
  - 4 flasks:  $3 \times 10^5$  cells \* 13 wells → 4.16 mL
- Transfection:
  - APOBEC plate
  - Row 1: pRetro empty
  - Row 2: pRetro-APO2X-dCas9
  - Row 3: pRetro-APO3X-dCas9

Transfected HEK293T cells with pcDNA, peGFPN3, ATG+ACG eGFP  
Created cultures for interlab in CAM


August 3, 2016



Gel 1 – Digest with only  
xho1  
Lane 1 – Ladder  
Lane 2 – Empty  
Lane 3 – Neg. Control  
Lane 4 – Empty  
Lane 5 – Positive Control  
(using Apobec clone)  
Lane 6 – Empty  
Lane 7 – Adar 1-3x #1  
Lane 8 – Empty  
Lane 9 – Adar 2-1x #2

#### Analysis

They all ran too small, even  
the positive control, which  
is known to be a clone.  
Something may be wrong  
with the xho1 cut cite or the  
enzyme itself.



Gel 2 – Digest with xho1 and EcoR1

Lane 1 – Ladder

Lane 2 – Empty

Lane 3 – Neg. Control

Lane 4 – Empty

Lane 5 – Positive Control

Lane 6 – Empty

Lane 7 – Adar 1-3x #1


Lane 8 – Empty

Lane 9 – Adar 2-1x #2

#### Analysis

The clones ran correctly for having the dCas9 insert and the editing enzyme insert.


We will run a gel with just EcoR1 to see if something is wrong with xho1.



Gel 3 – Digest with only  
EcoR1  
Lane 1 – Ladder  
Lane 2 – Empty  
Lane 3 – Negative Control  
Lane 4 – Empty  
Lane 5 – Positive Control  
Lane 6 – Adar 1-2x #1  
Lane 7 – Empty  
Lane 8 – Adar 2-1x #3

#### Analysis

The positive control and Adar 1-2x #1 ran correctly (they have dCas9 and the editing enzyme in the vector). However Adar 201x #3 did not run correctly. We can assume from this that there is something wrong with the xho1 enzyme. We will do a test digest with only EcoR1 with the rest of the clones to test them.



Gel 5 – Digest of all Clones  
with EcoR1 only

Lane 1 – Ladder

Lane 2 – Neg. Control

Lane 3 - Adar 1-1x #1

Lane 4 – Adar 1-2x #1

Lane 5 – Adar 1-3x #2

Lane 6 – Adar 2-1x #2

Lane 7 – Adar 2-2x #1

Lane 8 – Adar 2-3x #3

Lane 9 – Apobec 1x #3

- Transfection

eGFP trial 3

1  $\mu$ g pgRNA + 3  $\mu$ g of:

- pEGFP N3 (+)
- pcDNA (-)
- ATG GFP (ATG)
- ACG GFP (ACG)

Adjust to 200  $\mu$ L with 0.1X TE  $\rightarrow$  add 200  $\mu$ L 2X HBS

Add 100  $\mu$ L to corresponding well

- Fixing cells

eGFP trial 2

- New cell culture media

200 mL

- Induction

Doxycycline  $\rightarrow$  APOBEC plate

- Column 1: 0
- Column 2: 10 ng/mL
- Column 3: 100 ng/mL
- Column 4: 1000 ng/mL

Transfection, fixed cells  
Performed Interlab Experiment

## August 4, 2016

- Imaging
  - eGFP trial 2
- Cell plating
  - Cell count: 335; 323 → Average = 329
    - $1.316 \times 10^6$  cells/mL
  - Plating:
    - 12-well plate:  $3 \times 10^5$  cells \* 13 wells → 2.96 mL cells
- Harvest APOBEC plate (for Western blot)
  - 1.5 mL 2xSDS sample buffer + 30  $\mu$ L DTT
  - Heat at 80 °C, 10 min
  - Store in freezer
  - Gibson of globin prep
  - PCR amplified globin
  - Interlab Data Management
  - PCR'd ACG/ATG eGFP using biobrick primers



August 5, 2016


Ran globin on gel  
PCR purified globin  
PCR enzymes for biobrick

Gibson of globin onto Kan plates (3uL vector, 6uL insert)  
Gibson on gRNA to AMP plates  
(1uL vector, 9uL insert)


Key:  
purple ACG sgRNA  
Red globin sgRNA  
brown scramble  
black and purple sgRna neg control  
green WT  
orange PTC  
black globin neg control

PCR order  
ADAR1-1, 1-2, 1-2, 2-1, 2-2, 2-3, APOBEC1, 2, 3, 3x linker

Gel 1:  
Ladder, ADAR1-1, 1-2, 1-2, 2-1, 2-2, 2-3, APOBEC1, 2, 3, 3x



Gel 2:  
Ladder, 3x linker



Nothing there

Tested PCR products for eGFP biobrick parts  
Gibson PCR product into BioBrick plasmids

- Fixing cells  
eGFP trial 3
- Transfection  
eGFP trial 2 (redo)
- Imaging  
Exposure: Green = 26.17; Red = 1432  
46 images total