

Bee T



PCR product clean-up (Spin Column Extraction)

NucleoSpin® PCR clean-up Gel Extraction Kit

1. DNA Purification

- a. Adjust DNA binding condition, For very small sample volumes < 30 µL adjust the volume of the reaction mixture to 50–100 µL with water. It is not necessary to remove mineral oil. Mix 1 volume of sample with 2 volumes of Buffer NTI (e.g., mix 100 µL PCR reaction and 200 µL Buffer NTI).

(Note: For removal of small fragments like primer dimers dilutions of Buffer NTI can be used instead of 100 % Buffer NTI. Please refer to section 2.3. Incubate the mixture at 50°C for 5-10 min or until the gel has completely melted.)

- b. Place a NucleoSpin® PCR clean-up Gel column in a provided 2 ml collection tube. Apply up to 700 µl of the DNA/agarose solution to the NucleoSpin® PCR clean-up Gel column, and centrifuge at 11,000 x g for 30s at room temperature.
- c. Discard liquid and place the NucleoSpin® PCR clean-up Gel column back into the same collection tube. For volumes greater than 700 µl, load the column and centrifuge successively, 700 µl at a time. Each NucleoSpin® PCR clean-up Gel column has a total capacity of 25 µg DNA. If the expected yield is larger, divide the sample into the appropriate number of columns.
- d. Add 700 µl of Buffer NT3 into the NucleoSpin® PCR clean-up Gel column. Centrifuge at 11,000 x g for 30s at room temperature to wash the column. Discard the flow-through and re-use the collection tube.
- e. Repeat step k with another 700µl of Buffer NT3
- f. Discard liquid and centrifuge the empty NucleoSpin® PCR clean-up Gel column for 1 min at 11,000 x g to dry the column matrix. Do not skip this step, it is critical for the removal of ethanol from the NucleoSpin® PCR clean-up Gel column.
- g. Place a NucleoSpin® PCR clean-up Gel column into a clean 1.5 ml eppendorf tube. Add 15-30µl (depending on desired concentration of final product) of Elution Buffer (10 mM Tris-HCl, pH 8.5) directly onto the column matrix and incubate at room temperature for 1 minute. Centrifuge for 1 min at 11,000 x g to elute DNA. This represents approximately 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

