

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU-Eindhoven

SDS-PAGE Analysis

Table of contents

SDS-PAGE Analysis

1	SDS-Page Analysis	3
1.1	Preparation of reagents	3
1.1.1	Materials	3
1.1.2	Protocol	3
1.2	Sample preparation	4
1.2.1	Materials	4
1.2.2	Protocol	4
1.3	Preparation of the gel	5
1.3.1	Materials	5
1.3.2	Protocol	5
1.4	Running of the gel	5
1.4.1	Materials	5
1.4.2	Protocol	5
1.5	Staining of the gel	6
1.5.1	Materials	6
1.5.2	Protocol	6

1 SDS-Page Analysis

1.1 Preparation of reagents

Estimated bench time: 30 minutes

Estimated total time: 30 minutes

Purpose: Preparation of buffers and solutions.

1.1.1 Materials

- Balance
- Bottle (glassware)
- Bromophenol blue
- dH₂O
- ddH₂O
- Glycerol
- Glycine
- DTT
- SDS
- Tris base
- Tris-HCl
- β -mercaptoethanol

1.1.2 Protocol

- For 1 M DTT, the correct amounts are:
 - 2.32 g DTT
 - 15 mL dH₂O
- For 2X sample buffer SDS-PAGE, the correct amounts are:
 - 1 mL 1M Tris-HCl, pH 6.8
 - 10 mL dH₂O
 - 2 mL glycerol
 - 2.5 mL β -mercaptoethanol
 - 4 mL 10% SDS
 - 500 μ L 1% Bromophenol blue
- For 2X Laemmli buffer, the correct amounts are:
 - 450 μ L of 2X sample buffer SDS-PAGE
 - 50 μ L 1M DTT
- For 1X SDS running buffer, the correct amounts are:
 - 1.8 L ddH₂O
 - 2 g SDS
 - 28.8 g glycine
 - 6.04 g Tris base

1.2 Sample preparation

Estimated bench time: 20 minutes

Estimated total time: 45 minutes

Purpose: Preparing samples for loading on the SDS-PAGE gel.

1.2.1 Materials

- 10 μL 2X Laemmili buffer
- 2X Laemmili buffer
- BugBuster™ Protein Extraction Reagent
- Culture tube with cells and expressed proteins
- DTT stock (1M)
- Eppendorf tubes
- Heat/shaking block
- MiniSpin Centrifuge
- Tabletop Centrifuge
- Thermal cycler
- Vortex

1.2.2 Protocol

- Spin down the culture tube with cells and expressed proteins for 15 minutes at 3000 xg.
- Discard supernatant.
- Add BugBuster™ Protein Extraction Reagent. This needs to be 1/50th volume of the original volume of the culture.
- Shake at 250 rpm at room temperature until the pellet is dissolved.
- Spin down for 10 minutes at 13,400 rpm.
- Transfer the supernatant to a new Eppendorf tube.
- Mix 10 μL sample with 10 μL 2X Laemmili buffer in a PCR tube.
- Mix by vortexing and short-spin to get the samples at the bottom of the tubes.
- Denature for 5 minutes at 95 °C in the thermal cycler.

1.3 Preparation of the gel

Estimated bench time: 10 minutes

Estimated total time: 10 minutes

Purpose: Preparing a SDS-PAGE gel (10%) for loading the samples.

1.3.1 Materials

- 1X SDS running buffer
- Precasted SDS-PAGE gel (10%)
- SDS-PAGE kit components:
 - Cell for SDS-PAGE gels
 - SDS-PAGE casting tray
 - SDS-PAGE dummy

1.3.2 Protocol

- Take a cell for SDS-PAGE gels and fill this with 1X SDS running buffer for 50%.
- Take a precasted gel (10%) and take it out of the packing. Remove the strip at the bottom of the gel.
- Take a casting tray and clamp the gel (in the right direction!) and the dummy in this tray.
- Fill the chamber between the gel and dummy completely with 1X SDS running buffer.
- Make sure the chamber is not leaky.

1.4 Running of the gel

Estimated bench time: 45 minutes

Estimated total time: 5 minutes

Purpose: Analysis of the protein product with a SDS-PAGE gel.

1.4.1 Materials

- Extra-long pipet points
- Precision Plus Protein™ Dual color protein ladder
- Prepared loading samples
- Precasted SDS-PAGE gel (10%)

1.4.2 Protocol

- Load 10 μ L of the prepared loading samples on the gel with extra-long pipet points.
- Load 10 μ L of Precision Plus Protein™ Dual color protein ladder on the gel.
- Run the gel for 40 minutes on 200 V.

1.5 Staining of the gel

Estimated bench time: 20 minutes

Estimated total time: 2.5 hours

Purpose: Analysis of the protein product with a SDS-PAGE gel.

It is essential to work with gloves when working with the staining reagents used.

1.5.1 Materials

- Bio-Safe™ Coomassie G-250 stain
- dH₂O
- ImageQuant
- Rocking plate
- Tray (for transporting the SDS-PAGE gel)
- Tin foil

1.5.2 Protocol

- Carefully take the gel out of the casting frame and put it in a tray.
- Wash the gel with a layer of dH₂O for 10 minutes on a rocking plate. Cover the tray with tin foil.
- Stain the gel with a layer of Bio-Safe™ Coomassie G-250 stain for 1 hour on a rocking plate. Again, cover the tray with tin foil.
- Discard the Bio-Safe™ Coomassie G-250 stain in the sink.
- Destain the gel with a layer of dH₂O for 1 hour on a rocking plate.
- Create an image of the gel using the ImageQuant.