

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU_Eindhoven

NanoDrop

Table of contents

NanoDrop

1	NanoDrop
1.1	Materials
1.2	Setup & Protocol

Error! Bookm
Error! Bookm
Error! Bookm

1 NanoDrop

Estimated bench time: -

Estimated total time: 5 minutes start-up, 2 minutes per sample.

Purpose: Determining the concentration of samples.

1.1 Materials

- Autoclaved dH₂O
- dH₂O
- Fiber-free tissues
- NanoDrop spectrophotometer
- Pipettes and tips
- Samples

1.2 Setup & Protocol

- Start the NanoDrop spectrophotometer.
- Select the DNA measurement 'Nucleic Acid' in the NanoDrop menu for examining DNA sample, or;
- Select the protein measurement 'UV-Vis' in the NanoDrop menu for examining proteins by absorption.
- Clean the surface of the NanoDrop with dH₂O and a fiber free-tissue.
- Perform a calibration and blank measurement by entering one drop of 2 µl autoclaved dH₂O.
- Clean the surface again and place 2 µl per sample on the NanoDrop and measure the concentration.