

iGEM TU/e 2016

Biomedical Engineering

Eindhoven University of Technology

Room: Ceres 0.04

Den Dolech 2, 5612 AZ Eindhoven

The Netherlands

Tel. no. +31 50 247 55 59

2016.igem.org/Team:TU_Eindhoven

PCR Purification

Table of contents

PCR Purification	1	PCR Purification	Error! Bookm
	1.1	Materials	Error! Bookm
	1.2	Setup & Protocol	Error! Bookm

1 PCR Purification

Estimated bench time: 30 minutes

Estimated total time: 30 minutes

Purpose: Purifying the product obtained from a PCR reaction.

Important before starting are the following points (retrieved from QIAGEN):

- Add ethanol (96–100%) to Buffer PE before use (see bottle label for volume).
- Add 1:250 volume pH indicator I to Buffer PB. The yellow color of Buffer PB with pH indicator I indicates a pH of 7.5.
- Add pH indicator I to entire buffer contents. Do not add pH indicator I to buffer aliquots.
- If the purified PCR product is to be used in sensitive microarray applications, it may be beneficial to use Buffer PB without the addition of pH indicator I.

1.1 Materials

- Autoclaved Eppendorf tubes
- Autoclaved H₂O
- MiniSpin Centrifuge
- PCR tubes with PCR product
- Pipettes and tips
- QIAquick PCR Purification Kit

1.2 Setup & Protocol

- Mix 1 volume of the PCR sample with 5 volumes of buffer PB.
- Load the sample on a QIAquick spin column which is inserted in a collection tube (with a maximum of 800 µl per run).
- Centrifuge the sample for 1 minute at 13,400 rpm. Weight-balance the sample well.
- Discard the flow-through.
- Wash the sample with 750 µl of PE buffer and centrifuge for 1 minute at 13,400 rpm.
- Discard the flow-through.
- Dry spin the sample for 1 minute at 13,400 rpm.
- Transfer the spin column to a new autoclaved Eppendorf tube.
- Load 42 µl of autoclaved H₂O on the column (pipette drops in the middle of the membrane, do not touch the membrane). Incubate for 1 minute and centrifuge for 1 minute at 13,400 rpm.
- The resulting elution product will contain purified PCR product