

iGEM 2016: Pittsburgh
Week 3 Lab Notebook

Monday, June 6

Lab meeting

[Transformation](#), iGEM protocol

amilCP BBa_K592009 Kit Plate 1, 19E

T7 RBS (strong) BBa_K525998 Kit Plate 1, 3M

Terminator BBa_B0015 Plate 3, 3F

66 uL competent cells per tube (one 68 uL)

5 uL resuspended DNA per tube

To resuspend: pierce foil w/ tip, then add dH₂O

(Ice for > 30 min, b/c water bath not at temp yet)

Incubate 2 PM - 4 PM (~20 C)

Plate 100 uL each (LB + CM)

OUTREACH

Praneeth emailed Steve (tissue engineering camp)

Tuesday, June 7

Many more colonies from transformation

5 mL [liquid cultures](#) in 15 mL centrifuge

5 mL LB + 5 uL CM (1:1000)

Each person makes one (4*3 = 12)

Into room temp shaker 1:40 PM

Wednesday, June 8

[Mini-prep](#) of liquid cultures

Wash solution + 170 mL, 100% ethanol

Transfer cells into 1.5 tube to pellet @ 8000 rpm

1 liquid culture * 3 plasmids = 3 (for glycerol stocks)

From 6/6 plates

Into shaker at 12:15 PM

40% glycerol solution

Final volume 25 mL: 20 mL 50% glycerol and 5 mL H₂O (10 mL glycerol)

Thursday, June 9

[Long-term bacterial stocks](#) from yesterday's liquid cultures -- in -80 box with Top 10 cells

DNA concentrations from yesterday's mini-prep

Blank w/ elution buffer

(unit: ng/uL)	T7-RBS	amilCP	terminator
CC	74.6	123.1	88.6
MLK	108.5	109.0	88.5
PP	76.3	245.2	96.9
AN	64.3	118.2	58.7

Restriction digest

T7-RBS MLK

2 uL SpeI
2 uL 10X buffer
1 ug * 1000 ng/ug * uL/108.5 ng = 9.22 uL DNA
6.78 uL H₂O

20 uL Total

T7-RBS PP

2 uL SpeI
2 uL 10X buffer
13.1 uL DNA
2.9 uL H₂O

20 uL Total

amilCP CC

1.5 uL SpeI
1.5 uL XbaI
3 uL 10X buffer
2 ug = 16.25 uL DNA
7.75 H₂O

30 uL Total

amilCP PP

1.5 uL SpeI
1.5 uL XbaI
3 uL 10X buffer
8.16 uL DNA
15.84 uL H₂O

30 uL Total

Terminator CC

2 uL XbaI
2 uL 10X buffer
11.29 uL DNA
4.71 uL H₂O

20 uL Total

Terminator PP

2 uL XbaI
2 uL 10X Buffer
10.32 uL DNA
5.68. uL H₂O

20 uL Total

General Case (20 uL)

2 uL enzyme
10X buffer → 1X
Plasmid 1-2 ug
Nuclease-free H₂O to volume

Incubate @ 37 C

Incubate @ 80 C for 20 min (deactivate enzymes)

Gel electrophoresis

1% agar in 1X TAE buffer (50 mL, .5 g agarose, microwave)

Add 5 uL ethidium bromide

Pour gel

Lane 1 - 3-4 uL ladder, 4 uL loading dye

Lane 2 - 1 uL T7 MLK control (uncut), 2 uL dye

Lane 3 - 1 uL T7 MLK in buffer

Lane 4 - 1 uL T7 PP in buffer

Lane 5 - 1 uL amilCP control (uncut), 2 uL dye

Lane 6 - all amilCP CC in buffer


Lane 7 - all amilCP PP in buffer

Lane 8 - 1 uL terminator control (uncut), 2 uL dye

Lane 9 - 1 uL terminator CC in buffer

Lane 10 - 1 uL terminator PP in buffer

T7-RBS good, amilCP okay, terminator questionable--XbaI issue?



[Purify](#) linearized plasmids (T7 and terminator) (Maya)

Concentrations: Terminator CC -- 43.9 ng/uL
Terminator PP -- 39.0 ng/uL
T7-RBS MLK -- 58.2 ng/uL
T7-RBS PP -- 45.8 ng/uL

[Cut and purify](#) amilCP bands (Praneeth)

Concentrations: 6.6 ng/uL
9.0 ng/uL

[Sequencing](#) reaction (Pitt's [Genomics Research Core](#)) (Claire)

[Sample submission protocol](#)

15 uL

8 pmoles of primer (10 picomolar) = 100 uL H₂O + .1 uL 10 mM

1 ug DNA -- 10 mM

T7 MLK *Forward* -- CC1

9.22 uL DNA

1 uL 10 pM forward primer

4.78 uL H₂O

Backward -- CC2

9.22 uL DNA

1 uL 10 pM backward primer

4.78 uL H₂O

amilCP PP *Forward* -- CC3

4.08 uL DNA

1 uL 10 pM forward primer

9.92 uL H₂O

Backward -- CC4

4.08 uL DNA

1 uL 10 pM backward primer

9.92 uL H₂O

Terminator PP *Forward* -- CC5

10.32 uL DNA

1 uL 10 pM forward primer

3.68 uL H₂O

Backward -- CC6

10.32 uL DNA

1 uL 10 pM backward primer

3.68 uL H₂O

Friday, June 10

Reattempt [Restriction Digest](#) of Terminator with XbaI-- 30 uL

3 uL XbaI

3 uL 10X buffer

11.30 uL MLK terminator (1 ug)

12.7 uL nuclease-free water

Ligation reaction

Dephosphorylate T7-RBS

DNA = 10 uL

Add 0.5 uL rSAP (recombinant Shrimp Alkaline Phosphatase)

37 C for 30 minutes

65 C for 5 min

Ligation

Control Plasmid: 1 uL with 58.2 ng/uL

(no insert) Insert : 0

T4 buffer: 2 uL

T4 DNA Ligase: 1 uL

H2O up to 20: 16 uL

3:1 Plasmid: 1 uL with 58.2 ng/uL

Insert @ 50 ng: 5.56 uL of 90 ng/uL

Buffer: 2 uL

T4 DNA Ligase: 1 uL

H2O: 10.44 uL

~5:1 Plasmid: 1 uL of 58.2 ng/uL

Insert @ 81 ng: 1.7 mL of 9 ng/uL and 9.1 mL of 6.6 ng/uL

Buffer: 2 uL

Ligase: 1 uL

H2O: 2.46 uL

T4 Ligase should be added last

Room temperature for 10 min, heat inactivate at 65 C for 10 min

Transformation


Use 5 uL of each with Cheryl's Top 10 cells (200 uL aliquot, so 66.66 uL for each)

66 uL of cells in each tube, 5 uL of control, 3:1, and 5:1 solutions

Gel electrophoresis of terminator (MLK)

1 kb ladder (lane 3), uncut (lane 4) and cut (lane 5)


Gel run for 45 minutes (seems to work)



Phosphatase (rSAP) added to MLK terminator

Saturday, June 11

Remove colonies from incubator (Maya)



Sunday, June 12

Liquid cultures of T7 -- amilCP colonies (10 total cultures, 5 each of 3:1 and 5:1)
(Praneeth)