

## Mechanical stimulation

---

### Ultrasound stimulation

For real-time imaging of samples a Leica TCS SP5 laser scanning microscope mounted on a Leica DMI 6000 CS inverted microscope (Leica Microsystems, Germany) with an HCX plan apo 40× oil immersion objective was used. For image analysis we used CaPTURE software, developed by our team.

- Cells were seeded on a 6-well glass-bottom plate (Cellvis) and after 24 h transfected with selected plasmids.
- Cells were loaded with fluorescent calcium indicators and incubated for 30 min at 37°C and 5% CO<sub>2</sub>.
- Medium was exchanged for fresh medium with 4 mM CaCl<sub>2</sub> 10 min before starting the stimulation.
- Cells were stimulated for 10 s and responses were analyzed in CaPTURE.

---

### Stimulation of cells with touch, ultrasound and shaking

For detection of luminescence we used G:BOX (Syngene). For image analysis we used ImageJ (Image Processing and Analysis in Java) software (<http://rsbweb.nih.gov/ij/>).

- Cells were seeded in a petri dish (61 mm, TPP) and after 24 h transfected with selected plasmids.
- 24 h after transfection the medium was exchanged for medium with 1 mM luciferin and 4 mM CaCl<sub>2</sub>.
- Cells were incubated for 30 min at 37°C, 5% CO<sub>2</sub>.
- After the incubation cells were transferred to G:BOX machine and stimulated touching with a glass rod, by ultrasound or shaking.
- Successive images were acquired with 60 s exposure time after stimulation.
- Images were processed using ImageJ software.