

iGEM Notebook

Made with Benchling

Project: Protocols

Authors: Kamila Drezek

Dates: 2016-05-19 to 2016-10-09

THURSDAY, 5/19/16

Preparing Reagents

SATURDAY, 5/21/16

LB Broth (1L)

20 grams per liter

LB Agar (1L)

40 grams per liter

Prepare two liter containers of each- one per antibiotic (Cam and Carb)

Autoclave for 1 hour

Preparing chloramphenicol stock (1000x)

Final concentration of 1x chloramphenicol should be 25 ug/ml

Stock concentration of 1000x chloramphenicol should be 25 mg/ml

Prepare 50 ml of stock to aliquot into microcentrifuge tubes by measuring 1.25 g of chloramphenicol and ~50ml of 200 proof ethanol in a 50 ml Falcon

Filter chloramphenicol through 20 um syringe filter

Preparing Plates

20ml agar per plate

Innoculate Vector 1721

5 ml of LB+carb

Stab cell stock

Grow overnight at 37C shaker

SUNDAY, 5/22/16

Preparing Cell stock of 1721

Box Name- iGEM 2016 Source Box

Freezer Number- 3

Shelf- 3

Vector 1721 is miniprep

See protocol sheet

Vector 1721 is digested (300ul digest)(1-2hrs)

24000 ng digestion of DNA

30 ul buffer (cold room)

9ul of BpsI (-80C, small green box)(incubation in 37C)

101 ul of water

Digestion Clean-Up of Vector 1721

See PCR cleanup protocol

Note- vector was eluted in 50 ul of EB not 30 ul.

Dilution of Oligos (ideal final concentration is 100uM)

Short Cut-

300 ul of TE for values ~25- ~35 nmol

200 ul of TE for value ~15 - ~24 nmol

Note- Oligo 12, 13, 14, 17, 18, and 20 was vortex without liquid initially

Annealing of pGEX oligos

See protocol

Annealing time- 2 minutes

dilute oligos to 1ng/ul

Ligation of 1721 and insert

Note- diluted annealed oligos named- pGEX#A where # is the number of the construct

24 reactions

139.2 ul of water

24 ul of T4 Ligase

24 ul of T4 Ligase Buffer

26.4 ul of Vector

Transformation

See protocol

Plates put in incubator at 7:09 PM- Desired pick up time is 10:09

MONDAY, 5/23/16

pGEX plates were pulled from incubator at 10:15

Growth and density of all samples exceeded that of the negative control

TUESDAY, 5/24/16

Picked two colonies per plate and grew in 4ml of LB+carb in 37C shaker.

Created stab plate for each colony.

WEDNESDAY, 5/25/16

Minipreping colonies

Notes: During the miniprep samples from colony 10 were contaminated

Test Cut

4 ul DNA

1 ul 2.1


.25 of ul each enzyme (Bsp1 and Nco1)


4.5 ul water

Run digestion for at least 15 minutes

Test Cut Results

All colonies appeared with proper bands at 1.1 kb and 4.2kb and no bands at 440 bp. We decided to send only the first colony from each set in for sequencing





THURSDAY, 5/26/16

Sequencing results are in

All sequences with the exceptions of pGEX1-1, pGEX14-1, and pGEX 15-1 were successful.

All sequences remaining sequences were then grown up for cell stock

Note: While pGEX 6-1 sequenced properly the tube containing its reference appeared to have been contaminated and pGEX 6-2 will be sent in for sequencing to reverify

Sequencing sent in

pGEX: 1-2, 6-2, 14-1, 14-2, 15-1, 15-2, 10-1, 10-2

Restriction digest of 2490 and DC16 with MluI and NruI

5ul buffer3.1

6ul DNA

1.5ul NruI


1.5ul MluI

36ul water

Prerequisites to subvectors which will be used to create one plasmid with Cre and both Flps

Subvectors to be made have homology with each other, allowing for easy gibson assembly with no PCR (promoters and ERT2 - modified estrogen receptor - have repeat regions)

Digestion of DC16 yielded 3 bands rather than the expected 2. Culprit may be unreliable cell stock.



We were supposed to have:

at 10723bp, 1172bb.

2490 bands at 10590bp and 227bp.

DC16 bands

Ligate 2490 and DC16 to create 2534 (10ul)

See Ligation Calculator

4.2 ul DI water

1.0 ul DC16 (insert)

2.8 ul 2490 (vector)

1 ul T4 Ligase Buffer

1 ul T4 Ligase

Construct 2534 and a negative control are transformed

Tissue Culture Room

Prepping hood for transfection

Remember to clean the hood down with 70% ethanol

Aspirate all media from the plate to aid in the trypsin function

Kidney cells are stored in T75 flask

After aspiration perform a PBS wash

Phosphate Buffer Solution- a salt water solution (~5ml wash)

Trypsin is used remove kidney cells from the plate.

Add 2.5ml of Trypsin to the T75

Dont shoot directed at the cells, let the trypsin run over them.

Addition of Medium will neutralize trypsin

Add 7.5 ml of media to plate

And tritate (Pipette up and down)

To improve cell count- spin cells down then resuspend them

Spin at 300 RCF for 2-5 minutes

To resuspend- tap the hood with the 50ml falcon tube

Add back 10 ml of media and tritate

Performing a cell count

Add 10 ul of cell solution to a cell counting slides

Divide the cell count concentration by 2

Actual Cell Counts

~2.5 Million cells per ml

Magic Number for transfection:

200,000 cells per ml

For each T75- 20 ml of solution

For desired concentration of cells-

Add 4.8 ml of cell concentrate

Add 55.2 ml of media

Continuing the cell line

Note:

Split/Passage: diluting cells in the same vessel

Split the parent cell line

FRIDAY, 5/27/16

pGEX Cells moved to Stock

1-13, 16-20

Sequencing results

6-2 and 1-2 are perfect

pGEX 14 and 15 failed to sequence again

Transformation for Construct 2534

No significant growth in the positive- less than what appeared in the negative -> failure

Digest a different aliquot of DC16 at a low concentration

Need 4ng of DNA with DNA at 45ng/uL so need 89 uL of DNA; will do 100 ish uL digest

89 uL DNA
1.5 uL MluI
1.5 uL NruI
10 uL Buffer 3.1
No water

Ligation of dc16 insert and 2490 vector

2.7 ul of vector
1.7 ul of insert
3.5 ul of water
1 ul of T4 Ligase
1 ul of T4 Ligase Buffer

Digestion of pSb1A3 and pBS1C3

40 ul of DNA
5 ul of Cutsmart Buffer
1.5 ul of EcoRI
1.5 ul of PstI
2 ul of Water
After digestion treat with Dnpi for 15 minutes
PCR Clean- up and elute in 30 ul of EB

Gibson of pSb1A3 (2533) and pSb1C3 (2531)

CONSTRUCT NAME 2531

Mass Vec =25

Ratio = 3

V A B C D W

SIZE 2000 1200 600 0 0

CONC 5.6 10 10 1 1

VOL 4.5 4.5 2.3 0.0 0.0 -1.2

CONSTRUCT NAME 2533

SIZE 2000 1200 600 0 0

CONC 23 10 10 1 1

VOL 1.1 4.5 2.3 0.0 0.0 2.2

Both the Gibson Assembly Transformation and the Ligation Transformation are over 10kb

See Protocol

Stable Line Integration of 2494-2496

See Protocol for Stable Line Transfection

Construct 900- Piggy back transpose

Cell status preceding transfection

Low cell count ~30-40% confluences

Thaw out 2 tubes of PEI

Use an incubator in a BL2 room

Note: We are making 4x the number of cells than what is required

Mixes

DNA part

.78 ul of transposase

3.2 ul of DNA

.15 M NaCl (3ml per tube)

PEI part

2ml of PEI

.15 M NaCl (10.5 ml)

Add 4 ml of PEI part to DNA part

Transfection


Don't pour directly onto the- shoot away from the cells

SATURDAY, 5/28/16

Create cell stock of pGEX10-1

Mlprep of 14-3/4/5 15-3/4/5 and 2534-1/2

Test Cut Results



pGEX 14-3, pGEX 15-5, 2534-2 => sent to sequencing

Cell inspection preceeding puromycin selection

~90% confluence

Puromycin Selection

4ul of concentrated puromycin - 2ug per ml to each flask

Wait 10 days for total transient plasmids are diluted outs.

SUNDAY, 5/29/16

Sequencing Results

pGEX 15-5 can be moved to cell stock

pGEX 14-3 failed

2534 failed

Sent Sequencing

pGEX 14-4, pGEX 14-5

Digestion

2531

4 uL DNA

4 uL DI

1 uL CutSmart

0.5 uL EcoRI and AgeI

2533

4 uL DNA

4 uL DI

1 uL CutSmart

0.5 uL EcoRI and AgeI

DC16 midiprep

1 uL DNA at 1:10 dilution

7 uL DI


1 uL Buffer 3.1

0.5 uL MluI and NruI

DC16 miniprep
2 uL DNA
6 uL DI
1 uL Buffer 3.1
0.5 uL MluI and NruI

05-27-2016 DC16 DC56 2531 test cut.jpeg

Printed on Friday, May 27, 2016 19:33:43 (PM)



dc16 Primers

BP0165, BP0162 reverse

BP0198, BP0176 forward

dc56 Primers

BP0165, BP0130 reverse

BP0038, BP0188 forward

pGEX 14-4 and pGEX14-5 moved to cell stock pending sequencing

TUESDAY, 5/31/16

Phosphorylating and Annealing Oligos for pGOP

Moved Construct 2531 to cell stock

Grow up more 2531 for ligation

Because of Chloramphenicol's effect on growth, we will grow this up for 24 hours

Oligos are at 1000x concentration- dilute to 1x (~1ng/ul)

Picked more colonies from 2533 and grew up overnight

WEDNESDAY, 6/1/16

Miniprep 2531 (grown again due to low initial concentration) and 2533 (initial sequencing was incorrect)

Digest 2531 with BbsI (50 uL total)

30 uL DNA

5 uL Buffer 2.1

2 uL BbsI

13 uL DI water

Ligate oligos with 2531 to create pGOP1-20 (10uL total)

See protocol

0.4 uL 2531 (BbsI)

1.3 uL annealed oligos

6.4 uL DI water

1 uL T4 Ligase Buffer


1 uL T4 Ligase

All pGOP plasmids were transformed and plated (25 ul) with Control

Test Cut 2533 with EcoRI Age1

06-01-2016 testcut of 2533.jpeg

id on Wednesday, June 01, 2016 18:55:56



2533-3/4/5 were sent in for sequencing

Prepped for creation of comp cells for tomorrow

Mix 25g PEG8000 with 30ml water and dissolve over night

THURSDAY, 6/2/16

Pick 2 colonies for each pGOP

Miniprepped dc16 and dc56


Prepared 200ml TSS Solution

dc 16 test cut with MluI and NruI

expected bands: 10723, 1172

06-02-2016 DC16 testcuts.jpeg

d on Thursday, June 02, 2016 20:04



DC16: 1,2,3,5 were good

dc 56 test cut with AgeIHF and EcoRIHF

expected bands: 6700, 1000

Primer 0170 for 2533 failed and this impeded the sequencing

Preparing Top10 comp cells

First OD (optical density) measurement: 0.2 (time: 2:15 p.m.)

Second OD measurement: 0.34 (time 2:45 p.m.)



Third OD measurement: 0.419 (time 3:05 p.m.)

Fourth OD measurment: 0.486 (time 3:35 p.m.)

FRIDAY, 6/3/16

miniprepped pGOP plasmids

test cut pGOP plasmids with Bbs1 and Spe1



grew up more _____ for midiprep

SATURDAY, 6/4/16

Picked three colonies of pGOP 20 (3, 4, and 5)

Grew up pGOP 20-3/4/5

Miniprep 2533-7/8/9

SUNDAY, 6/5/16


Miniprep 20-3/4/5

Test Cut for 2533 and pGOP 20 addition colonies

06-05-2016 pGOP20 and 2533 test cuts.jpeg



inted on Sunday, June 05, 2016 18:22:05 (PM)



MONDAY, 6/6/16

Miniprep pGOP20-R1/R2/R3

Test cuts for pGOP20-R1/R2/R3

Cell Stock pGOP 1-19

Grew up DC16-1 and DC56-1 for cell stock

Sent pGOP 20-3/4/5 and 2533-7/8/9 for sequencing

Diluted pGEX 1-20 and pGOP 1-19 to 50 ng/uL

Made a liter of LB + Carb

Nanodropped pGEX 1-20 and pGOP 1-19

Inoculated pGEX16-1 due to low concentration

Midiprep for Transfection

Prepared 150 mL of LB + Carb with pCJH2

Prepared 150 mL of LB + Carb with T40

Prepared 150 mL of LB + Carb with pBW363

TUESDAY, 6/7/16

Plan for transfection of pGEXs and pGOPs to be performed on Thursday



MM1:

CJH2 (5 μ L)
T40 (5 μ L) $\times 22(+2)$ } Repeat dispense 15 μ L
363 (5 μ L)

MM2:

CJH2 (5 μ L)
T40 (5 μ L) $\times 22(+2)$ } Repeat dispense 10 μ L

CONSTRUCT	mass g \times 460 (ng)	Volume g \times 460 (μ L)	Vol. 4 \times 460 (vol)
Txn marker: CJH2	625	1.25	5
dCas9-VPE (T40)	"	"	5
"pGEX" or 363	"	"	5
pGOP	"	"	5
TOTAL	250	5	20

DNA 20 μ L
NaCl 30 μ L
50 μ L

Big min of (44 \times 2)
PEI 8 μ L
NaCl 42 μ L
50 μ L

Σ
= 100 μ L

\downarrow
25 μ L/txn. (transfection)

Repeat dispense 25 μ L per well
via yellow pipettor.

Miniprep of pGOP12, pGEX16, DC174 and DC100 (both DCs will be used as positive controls transfection on Thursday)

Midiprep of pBW363 (Blank), T40 (dCas9), pCJH2 (BFP transfection marker), DC16 and DC56

Evaporated supernatant

WEDNESDAY, 6/8/16

Cell stocked pGOP 20

50 μ L Restriction digest of 2490 and DC16 with MluI and NruI to create 2534

5 μ L buffer3.1

6 μ L DNA

1.5 μ L NruI

1.5 μ L MluI

36 μ L water

incubate 1 hr in 37C

Got 3 bands so threw out digest and are now investigating our enzymes to see if there has been contamination

Stocked midiprep DNA with 1x TE

Grew up T40 for further minipreps

Prepared HEK293 plates for tomorrow's transfection

Add 250 μ L of 200000 cells/ml cell solution to each well on 48 well plate

THURSDAY, 6/9/16

Transfected HEK cells

PLATE 1								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP1	pGOP2	pGOP3	pGOP4	pGOP5	pGOP6	pGOP7	pGOP8
E	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
	pGOP9	pGOP10	pGOP11	pGOP12	pGOP13	pGOP14	pGOP15	pGOP16

PLATE 2								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP17	pGOP18	pGOP19	pGOP20	DC174	DC100	pGEX1	pGEX2
E	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
	pGEX3	pGEX4	pGEX5	pGEX6	pGEX7	pGEX8	pGEX9	pGEX10
	pGOP3	pGOP4	pGOP5	pGOP6	pGOP7	pGOP8	pGOP9	pGOP10

PLATE 3								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	pGEX11	pGEX12	pGEX13	pGEX14	pGEX15	pGEX16	pGEX17	pGEX18
D	pGOP11	pGOP12	pGOP13	pGOP14	pGOP15	pGOP16	pGOP17	pGOP18
E	pCJH2	pCJH2	pCJH2	pCJH2	WT			
F	T40	T40	T40	T40	WT			
	pGEX19	pGEX20	pGEX1	pGEX2	WT			
	pGOP19	pGOP20	DC174	DC100	WT			

Miniprep T40

Diluted pCJH2, pBW363, and T40 for transfection

Checked and rediluted pGEXs and pGOPs for transfection

note- pGEX11 is contaminated

pGEX7 is too low of concentration

FRIDAY, 6/10/16

Checked transfection plates under a microscope

all worked; however, Cre showed leakiness

Resent DC16-1 for sequencing due to QuintaraBio sequencing failure

Passaging HEK293 cells

MONDAY, 6/13/16

Miniprep pGOP1-20 (inoculated in Cam last night)

Cut DC16 and 2490 with MluI and NruI

8 uL DC16 or 11uL of 2490

1.5 uL MluI

1.5 uL NruI

5uL buffer 3.1

34 uL water or 31 uL water (2490)

Run Gel Extraction on 2490 and DC16

Bands

DC16


10723

1172 - WANTED BAND

2490

10590 - WANTED BAND

227



Ligation of 2490 and DC16 to make 2534


Ligations were transformed and plated

TUESDAY, 6/14/16

Diluted primers for PCR reaction of pGPX1 (2531)

PCR GFP, BFP and mRuby

Digested pGOPs 3,5,6,8,10,14,15,19 and pGPX1 with Mlu1-HF and Age1-HF in preparation of making pGOPs with GFP, BFP and mRuby



Miniprep 2531

Picked colonies for 2534

Transfected HEK cells with different concentrations of the 8 pGOPs given above, to test whether the concentration of pGOPs in the cell was causing high levels of basal expression)

Transfected transfection marker, T40 (dCas9) pGEXs, Blank plasmid and pGOP. pGOP at 62.5, 30 and 7.5 ng


Screen Shot 2016-08-01 at 12.17.24.png

PLATE 1								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP3	pGOP5	pGOP6	pGOP8	pGOP10	pGOP14	pGOP15	pGOP19
E	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
	pGEX3	pGEX5	pGEX6	pGEX8	pGEX10	pGEX14	pGEX15	pGEX19

PLATE 2								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP3	pGOP5	pGOP6	pGOP8	pGOP10	pGOP14	pGOP15	pGOP19
E	pCJH2	pCJH2	pCJH2	#REF!	pCJH2	pCJH2	pCJH2	pCJH2
F	T40	T40	T40		T40	T40	T40	T40
	pGEX3	pGEX5	pGEX6		pGEX10	pGEX14	pGEX15	pGEX19

PLATE 3								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP3	pGOP5	pGOP6	pGOP8	pGOP10	pGOP14	pGOP15	pGOP19
E	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
	pGEX3	pGEX5	pGEX6	pGEX8	pGEX10	pGEX14	pGEX15	pGEX19

PCR digested GFP, BFP and mRuby with Mlu1-HF and Age1-HF



WEDNESDAY, 6/15

Miniprepred T40 (was running low)

Ligated GFP, BFP and mRuby into pGOP 3,5,6,8,10,14,15,19. Transformed onto Cam plates

THURSDAY, 6/16

Ran FACs on transfection from Tuesday

PCR the silencer out of DC100 (has UAS and pG5) to put into pGOPs and pGPX6 with PrA4 and PrA4r

Digested DC 100 and prepared annealed oligos o83-o122 to prepare _____


FRIDAY, 6/17

Miniprepred Transformed plasmids from 6-15-16

PCR for GFP, BFP, and mRuby


06-17-2016 GFP BFP mRuby.jpeg

Printed on Friday, June 17, 2016 19:20:27 (PM)



06-17-2016 2534t and 2534b.jpeg

Printed on Thursday, June 16, 2016 16:45:05 (PM)





Sent plasmids in for sequencing

Sequence for 2534B-4 is verified as functioning and will be moved forward with cloning

06-17-2016 bfp gfp pgops.jpeg

Printed on Friday, June 17, 2016 22:34:58 (PM)





MONDAY, 6/20

Cell Stocked plasmids for BFP, GFP, and mRuby that sequenced properly

Minipreped plasmids

Picked colonies for mRuby pGOP 10, 14, and 15

TUESDAY, 6/21

Cell stocked BFP pGOP 15

Picked colonies for GFP pGOP 8, GFP pGOP 19, mRuby pGOP 10, mRuby pGOP 15, and mRuby pGOP 19

Minipreped plasmids from 6-20-16 (mRuby pGOP 10-3, 10-4, 14-3, 14-5, 14-6, 15-3)

Sent mRuby pGOP 14-3/4/5 and BFP pGOP 19-2 for sequencing

Passaged new HEK cells to P7

Made 2 new bottles of D5 Media (D5= 5%FBS)

500ml of DMEM base

25ml FBS (alloquated in fridge across hallway from TC room)

5ml L-Glut

5ml P/S


5ml sodium pyruvate

WEDNESDAY, 6/22

Cell stocked picked colonies from 6/21/16 (GFP pGOP 8, GFP pGOP 19, mRuby pGOP 10, mRuby pGOP 15, and mRuby pGOP 19).

THURSDAY, 6/23

Gel extract PrA4/4r (silencer)



Digested DC100 with NheI and dropped in annealed oligos with ligation to make pDOP 1-20 (backbone from DC100 with gRNAs dropped in)

sent pGOP 6-9 back in for sequencing - suspected error in gRNA target site

MONDAY, 6/27

ligate and transform pDOP 1-20 (DC100 and annealed oligos)

inoculate pGPX1 so we can make pGPX2 and 3

TUESDAY, 6/28

Prepped for transfection

Picked pDOP1-20


Gel extracted pGPX1


WEDNESDAY, 6/29

ligate and transform pGPX2 and pGPX3 (pGPX1 MluI AgeI with PrA1/1r and PrA2/2r MluI and AgeI)

miniprep pDOP 1-20

Test Cut pDOP1-20 with NheI-HF and StuI





Transfected pGOP 3,5,6,8,10,14,15,19 of iRFP, BFP, GFP, and mRuby


PLATE 1								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP3	pGOP5	pGOP6	pGOP8	pGOP10	pGOP14	pGOP15	pGOP19
E	IRFP							
F								
PLATE 2								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP23	pGOP25	pGOP26	pGOP28	pGOP30	pGOP34	pGOP35	pGOP39
E	gfp							
F								
PLATE 3								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP43	pGOP45	pGOP46	pGOP48	pGOP50	pGOP54	pGOP55	pGOP59
E	bfp							
F								
PLATE 4								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP63	pGOP65	pGOP66	pGOP68	pGOP70	pGOP74	pGOP75	pGOP79
E	mRuby							
F								

THURSDAY, 6/30

Picked pGPX2 and pGPX3 colonies

Reannealed-phosphorylated- ligated - and transformed pDOPs
recut pDOPs

06-30-2016 recut pDOPs.jpeg



FRIDAY, 7/1

Miniprep of GPX 2-1/2/3 and GPX 3-1/2/3

Test-cut GPX 2 and 3

Picked colonies and grew up pDOP 1-4; 6-11; 14-20


Ran FACS (flow cytometry)

Did not test BFP due to wrong transfection marker

PCR clean up of 2485; 2487; and 2489 with NotI and Cip

Gel extraction of 2484, 2486, 2488, and 2534 digested with AscI and NotI

07-01-2016 pGPX2 3.jpeg



SATURDAY, 7/2

Miniprep and test cut pDOP 1-4; 6-11; 14-20 colonies 1 and 2

digested with NheI and PstI

Vectors pGPX2 and pGPX3 are sequenced verified. pGOP 66,34, and 26 could not be sequenced properly (quintara error)

Reenoculate pGEX and GFP pGOP 3, 5, 6, 8, 10, 14, 15, 19 gh54

SUNDAY, 7/3

Miniprep pGEX and GFP pGOP 3, 5, 6, 8, 10, 14, 15, 19

Gibson and transform pREC 1, 2, 3

MONDAY, 7/4

Passage HEK cells

Inoculate pREC 1, 2, 3

TUESDAY, 7/5

Miniprep pREC 1, 2, 3


Enoculate BW1720

Digested BW1720, and BW1942-43 with BbsI


Digested BW758, and BW762-65 with AgeI and EcoRI

Digested DC 100 with NheI and Cip
Digested pGPX1-4 with XbaI and BbsI


07-05-2016 dc100, pGPX1:2:3:4.jpeg



07-05-2016 BW contstructs.jpeg



07-05-2016 pREC.jpeg



Transfected HEK cells for gRNA orthogonality experiment


PLATE 1								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP23	pGOP25	pGOP28	pGOP30	pGOP35	pGOP39	pGOP23	pGOP25
E	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
	pGEX3	pGEX3	pGEX3	pGEX3	pGEX5	pGEX5	pGEX5	pGEX5
	pGOP28	pGOP30	pGOP35	pGOP39	pGOP23	pGOP25	pGOP28	pGOP30

PLATE 2								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	pGEX5	pGEX5	pGEX8	pGEX8	pGEX8	pGEX8	pGEX8	pGEX8
D	pGOP35	pGOP39	pGOP23	pGOP25	pGOP28	pGOP30	pGOP35	pGOP39
E	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
	pGEX10	pGEX10	pGEX10	pGEX10	pGEX10	pGEX10	pGEX15	pGEX15
	pGOP23	pGOP25	pGOP28	pGOP30	pGOP35	pGOP39	pGOP23	pGOP25


PLATE 3								
	1	2	3	4	5	6	7	8
A	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2	pCJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	pGEX15	pGEX15	pGEX15	pGEX15	pGEX19	pGEX19	pGEX19	pGEX19
D	pGOP28	pGOP30	pGOP35	pGOP39	pGOP23	pGOP25	pGOP28	pGOP30
E	pCJH2	pCJH2						
F	T40	T40						
	pGEX19	pGEX19						
	pGOP35	pGOP39						

WEDNESDAY, 7/6

Miniprep pGEX, GFP pGOP, and mRuby pGOP 6& 14
 Digested 1720 with Agel and EcoRI



Redigested 2485, 2487, and 2489 with NotI and Cip



Ligated and transformed pBEX1-24 pPV1,5-8
 Retransformed pGEX6 and 14 pGOP26,34,66,74

Made LB + Carb plates

THURSDAY, 7/7

Picked colonies for pBEX 1 - 24, pGEX 6 and 14, GFP and mRuby pGOPs 6 and 14, adh pPVs 1, 5, 6, 7, 8


Miniprep pGEX 1-4

Ran FACs

Grew more 2534

FRIDAY, 7/8

IMG_9323 copy.jpg



MONDAY, 7/11

Grew up pPV7 and pBEX 10, 12, 13, 17, 18, 22

Growing more 2534

Gibsoned pREC1, 2, and 3

Passaged HEK cells

TUESDAY, 7/12


Cell Stock pPV7 and pBEX 10, 12, 13, 17, 18, and 22

Digested pGPX2, mRuby2, and pGOP 1, 2, 4, 6, 7, 9, 11, 12 13, 16, 17, 18, and 20

mRuby2 and pGOPs digested with MluI and AgeI

pGPX2 digested with BbsI

Ran on a gel and gel extracted




WEDNESDAY, 7/13

Miniprep pGPX2, 1720, 1942-44

Digested pGPX2, 1720, 1942-44 with _____

Re-digested pGPX2 (bbs1, extract 3kb band) ,mRuby2(mlu1, age1; extract 700bp band) ,and pGOP21 (mlu1,age1, extract 3kb band)

ran on gel and gel extracted



THURSDAY, 7/14


Ligating and transforming

pGOP 21 22 24 27 29 31 32 33 34 36 37 38 40 41 42 44 47 49 51 52 53 56 57 58 60 61 62 64 66 67 69 71 72 73 76 77 78
80

Digesting

1270 758 763 765 with Agel and EcoRI

07-14-2016 758 763 765 1720.jpeg




Circuit Designs Digesting

Part1 Ascl NotI
Part2 Ascl NotI
Part3 Ascl NotI
Part4 Ascl NheI
Dest EcoRI NotI

Circuits (P1-P2-P3-P4)

pBEX1-pBEX6-pBEX15-pBEX24
pBEX1-pBEX10-pBEX19-pBEX20
pBEX5-pBEX18-pBEX19-pBEX16
pBEX5-pBEX2-pBEX15-pBEX8

07-14-2016 pBEX 8 16 20 24.jpeg




07-14-2016 pBEXs.jpeg



Test cut 2578, 2579, 2580 (formely pRECs)


07-14-2016 2578 2589 2580.jpeg



FRIDAY, 7/15

Miniprep 2578, 2579, and 2580 3-7 & testcut

07-15-2016 2578 2579 2580 test cut.jpeg




2578-2 and 2579-1 were sequenced verified and moved to cell stock

Miniprep 1716 pDest

Digested pDest with EcoRI and NotI

07-15-2016 1716.jpeg



Religated pPV1, pPV6, and pPV8
Picked new pGOP colonies for GFP only

SATURDAY, 7/16

Miniprepped new GFP colonies
Miniprepped 2580-8 through 2580-17
Transformed pCir1 pCir2 and pCir3

SUNDAY, 7/17

Test Cuts for 2580 additional colonies
Test Cuts for pPV1, pPV6, and pPV8 additional colonies
Test Cuts for pGOP GFP additional colonies

MONDAY, 7/18


Test Cuts for pCIR1, pCir2, and pCir3
Test Cuts for pGOp1 and 2 under BFP and mRuby
Split HEK cells 1:10
Created two 48 well plates at 200,000 cells/ml for transfection of BFP pGOP constructs
Created two T175plates with 25ml of 100,000cells/ml for stable integration

TUESDAY, 7/19

Midiprep BW361, BW465, BW471, BW 474, 2578, and 2579
Digested pBEX 19 with Ascl and NotI
Gel extracedracted pBEX 19



Digested pCir 1-3, 1-1, 3-2, 3-3 with NotI



Transfected BFP

PLATE 1							
	1	2	3	4	5	6	7
A	BW471 (IRFP)	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40
B	T40	BW363	BW363	BW363	BW363	BW363	BW363
C	BW363	pGOP5	pGOP6	pGOP8	pGOP10	pGOP14	pGOP15
D	BW471	BW471	BW471	BW471	BW471	BW471	BW471
E	T40	T40	T40	T40	T40	T40	T40
F	pGEX3	pGEX5	pGEX6	pGEX8	pGEX10	pGEX14	pGEX15

PLATE 2							
	1	2	3	4	5	6	7
A	BW471	CJH2	BW471	WT blank			
B	BW363	BW363	CJH2				
C			BW363				
D							
E							
F							

PROJECT _____

Notebook No. _____

37

Continued from Page _____

construct	Mais in 2 wells	Vol 3 well	Vol 4 wells	(control) x 10 wells	(experimental) x 10 wells	* controls → special																																																									
BW471 (T40 master)	62.5ng	1.25μl	5μl	50μl	50μl	control																																																									
T40	62.5ng	1.25μl	5μl	50μl	50μl	control																																																									
pGEX 363	62.5ng	1.25μl	5μl	50μl	50μl	control																																																									
pGEX	7.5ng	1.25μl	5μl	50μl	50μl	control																																																									
<p>↓ dispense 19.4μl to each well add 10μl pGEX</p>						<p>control</p> <table> <tr> <th>1</th><th>max. 2 wells</th><th>Vol 3 well</th><th>Vol 4 wells</th></tr> <tr> <td>BW471</td><td>62.5ng</td><td>1.25</td><td>5μl</td></tr> <tr> <td>BW363</td><td>187.5</td><td>3.75</td><td>15μl</td></tr> <tr> <td colspan="4">20 tot</td></tr> </table> <p>2</p> <table> <tr> <th>1</th><th>max. 2 wells</th><th>Vol 3 well</th><th>Vol 4 wells</th></tr> <tr> <td>CJH2</td><td>62.5</td><td>1.25</td><td>5μl</td></tr> <tr> <td>BW363</td><td>187.5</td><td>3.75</td><td>15μl</td></tr> <tr> <td colspan="4">20 tot</td></tr> </table> <p>3</p> <table> <tr> <th>1</th><th>max. 2 wells</th><th>Vol 3 well</th><th>Vol 4 wells</th></tr> <tr> <td>BW471</td><td>62.5</td><td>1.25</td><td>5μl</td></tr> <tr> <td>CJH2</td><td>62.5</td><td>1.25</td><td>5μl</td></tr> <tr> <td>BW363</td><td>125</td><td>2.5</td><td>10μl</td></tr> <tr> <td colspan="4">20 tot</td></tr> </table>						1	max. 2 wells	Vol 3 well	Vol 4 wells	BW471	62.5ng	1.25	5μl	BW363	187.5	3.75	15μl	20 tot				1	max. 2 wells	Vol 3 well	Vol 4 wells	CJH2	62.5	1.25	5μl	BW363	187.5	3.75	15μl	20 tot				1	max. 2 wells	Vol 3 well	Vol 4 wells	BW471	62.5	1.25	5μl	CJH2	62.5	1.25	5μl	BW363	125	2.5	10μl	20 tot			
1	max. 2 wells	Vol 3 well	Vol 4 wells																																																												
BW471	62.5ng	1.25	5μl																																																												
BW363	187.5	3.75	15μl																																																												
20 tot																																																															
1	max. 2 wells	Vol 3 well	Vol 4 wells																																																												
CJH2	62.5	1.25	5μl																																																												
BW363	187.5	3.75	15μl																																																												
20 tot																																																															
1	max. 2 wells	Vol 3 well	Vol 4 wells																																																												
BW471	62.5	1.25	5μl																																																												
CJH2	62.5	1.25	5μl																																																												
BW363	125	2.5	10μl																																																												
20 tot																																																															

MM 3	20 μL
PEI	8 μL
NaCl	42 μL
50 μL	50 μL
100 μL	

MM 3	20 x2 → 40
PEI	320 μL
NaCl	1680 μL

dispense 25μL per 100μL well

WEDNESDAY, 7/20


THURSDAY, 7/21

Performed FACs on the BFP transfected cells. Tried to create the data in MEFL but had trouble with getting the channel names to work with the software.

Split HEK cells to P16. Expanded them into large T175 plate


Transfected the HEK cells in T175 with 2578 and 2579. Brought media volume up to 50ml to do so

IMG_1338.jpeg



FRIDAY, 7/22

07222016 pCir test cut2.jpeg



will did a testcut of pCir1 and pCir3. sent one of each in for sequencing.


SATURDAY, 7/23

added puromycin to the two stable line integrations of 2578 and 2579. 2ug per ml so we added 10ul of puromycin to 50ml of media in the flask

many more cells today than yesterday

Sent in a sample of the pGOP41 used in the previous transfection to troubleshoot what went wrong and why that particular gRNA was not working with BFP

07232016 BFP PCR.jpeg



SUNDAY, 7/24

Split the stable line integration HEK cells of 2578 and 2579 each into two new T175 flasks because they were too crowded. Added puromycin to the fresh media.

Transformed (GFP) pGOP 26, 27, 32, 33, 38, 40

(BFP): pGOP 41, 42, 43, 44, 46, 47, 49, 51, 52, 53, 56, 57, 58, 60


(mRuby): pGOP 61, 62, 64, 66, 67, 69, 71, 72, 73, 76, 77, 78, 80

MONDAY, 7/25

Picked colonies from all transformed GFP pGOPs and less successful BFPs and mRubys

Prepped overnight ligation for the BFPs and mRubys

due to miscommunication, the ligation mix was made incorrectly but we went ahead anyway



TUESDAY, 7/26

Digested pBEX 1,2,5,6,16,19 (Ascl and NotI except 16 which was Ascl and NheI) and 1716 (EcoRI and NotI) and pGOP68 and 35 (MluI and AgeI)

did not extract the pGOPs or 1716


extracted the rest



Inoculated pBEX 1,6,16,19, and 1716 and pGPX2 (so we could drop in the new oligos)

Transformed DNA from overnight ligation


made Cam plates




WEDNESDAY, 7/27

Transformation did not work

Digested 1716 (EcoRI and NotI) and pGOP68 and pGOP35 (MluI and AgeI)
gel extracted



annealed multimerized oligos
digested more pGPX



Diluted all the DNA for the transfection

Transfected HEK cells with the "mini circuit" (1902/1903, pGOP1, pGOP22, BW390 (constitutive cre) and 2196

also transfected pGOP22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38, 39 to test all of the GFPs. we did not have pGOP 21 and we did not have pGOP 40. We did not have pGEX11 to transfect so pGOP31 should not work.
organized our DNA boxes

PLATE 1								
	1	2	3	4	5	6	7	8
A	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	IRFP	IRFP
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP1	pGOP1	pGOP1	pGOP1	pGOP1	pGOP1	2196	2196
E	pGOP22	pGOP22	pGOP22	pGOP22	pGOP22	pGOP22		
F	IRFP	IRFP						
G	T40	T40						
H	1902	1903						
I	BW390	BW390						
J	2196	2196						

PLATE 2								
	1	2	3	4	5	6	7	8
A	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
D	pGOP22	pGOP23	pGOP24	pGOP25	pGOP26	pGOP27	pGOP28	pGOP29
E	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
G	BW363	BW363	BW363	BW363	BW363	BW363	BW363	BW363
H	pGOP30	pGOP31	pGOP32	pGOP33	pGOP34	pGOP35	pGOP36	pGOP37

PLATE 3								
	1	2	3	4	5	6	7	8
A	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2
B	T40	T40	T40	T40	T40	T40	T40	T40
C	BW363	BW363	pGEX2	pGEX3	pGEX4	pGEX5	pGEX6	pGEX7
D	pGOP38	pGOP39	pGOP22	pGOP23	pGOP24	pGOP25	pGOP26	pGOP27
E	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2
F	T40	T40	T40	T40	T40	T40	T40	T40
G	pGEX8	pGEX9	pGEX10	pGEX11	pGEX12	pGEX13	pGEX14	pGEX15
H	pGOP28	pGOP29	pGOP30	pGOP31	pGOP32	pGOP33	pGOP34	pGOP35

PLATE 4								
	1	2	3	4	5	6	7	8
A	CJH2	CJH2	CJH2	CJH2	GFP	BFP	IRFP	GFP
B	T40	T40	T40	T40				BFP
C	pGEX16	pGEX17	pGEX18	pGEX19				
D	pGOP36	pGOP37	pGOP38	pGOP39				
E	GFP	BFP	GFP	WT				
F	IRFP	IRFP	IRFP					

THURSDAY, 7/28

FRIDAY, 7/29


Digested more BFP and mRUBY

26ul pGOP55 30ul pGOP64

1.5 of MLU1-HF & Age1-HF 1.5 of Mlu1-HF & Age1-HF


5 cutsmart 5 cutsmart

20 h2o 16 h2o



ligated to create pGOP41, 51, 82, 83 to see what wasnt working with our cloning by testing our T4 ligase vs Wong labs t4 ligase
Testcut pCir1-6, -7, -8, -9, -10, -11, -12, -13, -14, -15, pCir3-4 with Acs1 Not1-HF

07292016 pCir test cut.jpeg



Cell stocked pGOP38-4, pGOP42-1, pGOP43-1, pGOP44-1, pGOP26-5, pGOP27-4, pGOP32-4, pGOP32-4, pGOP33-4, pGOP60-1, pGOP64-1

performed FACs

Split Deboki's cells for her and took a t175 plate at 1:10. (passage again on monday)

Moved the stable line integrations to T50 plates. added 2ul of puromycin to 10ml of media

Picked colonies for GOP 82, and 41, 51 (no colonies for 83), pCIR1s that Ben made (picked 10 colonies)

Met with Leidy to discuss what we are doing wrong with our ligation and transformation protocol.

Re did digesection of BFP and mRuby from pGOP 43 and 64, respectively. (age1 mlu1)

Redigested pgpx2 to drop in the oligos.

Religating fluorescence into regular pGOPs

SATURDAY, 7/30

Picked colonies of pGOP 66-1*(failed),-2, 67-1*(failed),-2, 71-1*(failed),-2, 73-1*(failed),-2, 76-1*(worked),-2, 46-1*(worked),-2, 47-1*(failed),-2., 51-1*(failed),-2, 52-1*(failed), 53-1*(worked),-2, 56-1(failed)*,-2, 49-1*(failed),-2, 69-1*(failed sequencing reaction),-2, 41-1*(worked),-2,-3 -4,-5, 51-1,-2,-3,-4,-5,-6, 40-1,-5

pGOP 81-1,-2, 82-1*, -2, 83-1,-2,84-1,-2, 85-1,-2, 82-0* (-0 means grew up day before), 86-1,-2,87-1,-2,88-1*, -2, 89-1,-2, 90-1,-2, 83-0

pCIR1 -16,-17,-18, -19,-20*, -21,-22,-23, -24,-25*, -26*, -27*, Ben's1, Ben's2*


(starred**** ones sent for sequencing)

SUNDAY, 7/31


miniprep of all of the colonies picked yesterday

Ran test cuts of all of the colonie picked


07312016 pCIR test cuts.jpeg




07312016 pGOP80 test cuts.jpeg



07312016 pGOP test cuts more.jpeg



07312016 pGOP cuts again.jpeg



MONDAY, 8/1

split HEK cells 1:10 from large flask into small flask.

let stable line of cells grow for another day.

Grew up.... pGEX11, pGEX12, pGEX13, pGEX17, pGEX18, pGEX20, pGOP6(x2), pGOP21, pGOP22, pGOP23, pGOP24, pGOP25, pGOP26, pGOP27, pGOP28, pGOP29, pGOP30, pGOP31, pGOP32, pGOP33, pGOP34, pGOP35, pGOP36, pGOP37, pGOP38, pGOP39, pGOP42, pGOP43, pGOP44(x2), pGOP45(x2), pGOP50, pGOP54, pGOP55(x2), pGOP59, pGOP68, pGPX2

TUESDAY, 8/2

(Will and Marisa) minipreped pGEX11, pGEX12, pGEX13, pGEX17, pGEX18, pGEX20, pGOP6(x2), pGOP21, pGOP22, pGOP23, pGOP24, pGOP25, pGOP26, pGOP27, pGOP28, pGOP29, pGOP30, pGOP31, pGOP32, pGOP33, pGOP34, pGOP35, pGOP36, pGOP37, pGOP38, pGOP39, pGOP42, pGOP43, pGOP44(x2), pGOP45(x2), pGOP50, pGOP54, pGOP55(x2), pGOP59, pGOP68, pGPX2

Grew up pGOP41, pGOP46, pGOP53, pGOP71, pGOP76 for cell stock


Split the stable line HEK cells into new T75 flasks (p18)

WEDNESDAY, 8/3

suspended the crRNA gBlocks oligos in 1xTE

Digested pGPX2 with BBs1

20160803 jeffery's stuff.jpeg



Digested pGOP76-1 and pGOP 41-1 for mRuby and BFP

Transformed pGOPs 57, 58, 61, 62, 66, 77, 78, 81, 84, 85, 86, 87, 88, 89 using Leidys transformation protocol

- use Ben's Top Ten cells, let cells sit on ice

- add 10ul of 5x KCM to ligation mix

- add 30ul of DI H₂O to ligation mix

- add all 50 ul to cells

- plated 95ul of cells onto 1/2 plate

Sent sequencing pGOP 47-1(worked), 49-1(worked), 51-2(worked), 52-1(worked), 56-2 (failed), 67-1(worked), 69-1(worked), 71-1(FAIL), 72-1(worked), 73-1(resendfor seuquencing), 83-0(resend for sequencnig)

Cell stocked pGOP41, pGOP46, pGOP53(need to throw awayb/c miscommunication), pGOP76 (throw away b/c miscommunication)

THURSDAY, 8/4

(Marisa) Picked 2 colonies for pGOPs 57, 58, 61, 62, 66, 77, 78, 81, 84, 85, 86, 87, 88, 89

(Kami) Split wildtype cells 1:20 to P19

(Jeffery) dimerized guides 1, 3, 8, 13 in various configurations with each other (1-3, 3-1, 8-13, 13-8)


- this involved PCR and a gel extraction

- Designed decoders

20160805 dec2.jpeg



20160805 decoder parts.jpeg



(Will) Grew up pGOP 47, 49, 51, and 82-0 for cell stock

Due to a miscommunication error, we had to retransform and plate pGOP52-1, pGOP53-1, pGOP67-1, pGOP69-1, pGOP71-2, pGOP72-1, pGOP73-1. We did this with our typical transformation protocol.

FRIDAY, 8/5

(Marisa and Will) Miniprepped pGOPs 57, 58, 61, 62, 66, 77, 78, 81, 85, 86, 87, 88, and 89

Due to a miscommunication error, miniprepped pGOP 47, 49, 51, and 82-0

(Will) Gibsoned

pCir 2 using 1716, A, B, C, D

pCir 2 using 1716 Part1/2, Part 3/4

pCir 2 using 1716, Decoder 2

pCir 4 using 1716, lambda, Part 3/4

pCir 1 using 1716, alpha, Bex19


pCir 1 using 1716, Part 1/2, Bex19, Bex16

pCir 1 using 11716 Bex 1, Bex6, Bex19, Bex16

pCir 3 using 1716, Bex5, Bex2, Bex19, Bex16

(Kami) Split stabel line (2578 and 2579) 1:10 in T75

(Rachel) Test-cut pGOPs 57, 58, 61, 62, 66, 77, 78, 81, 85, 86, 87, 88, and 89





SATURDAY, 8/6

(Rachel) digested pGPX2 (5mg) with 3ul of BbsI. Double digested - too low of concentration to proceed with ligation

SUNDAY, 8/7

(Will) grew up

MONDAY, 8/8

(Jeffery) Digested pgpX2 with bbsi

(kami) made WT Hek plates with Yash and split WT and 2578 and 2579 integrated cells 1:10

(will) cell stocked pGPX2

(Will) Test cut the pCir gibsons with AscI and NotI

(Will) transformed pGOP 84, 86, 87, 88, 89, 90 with a ligation using Leidy's protocol
transformed pGOP47, 49, 51, 82-0 without a ligation using Leidy's protocol

TUESDAY, 8/9

(Rachel) Growing pCir4_2_4 from stab plate as sequencing indicated it worked!

(Rachel) Picked colonies from pGOP 81,84,86,87,88,89,90

(Rachel) made more LB + Cam broth

(Kami) autoclaved tips and agar

(Marisa and Kami) made carb plates

(Jeffrey) performed overhang extension PCRs on G-blocks (that make decoder) to build Circuit 2

WEDNESDAY, 8/10

(Marisa and Will) miniprep pGOP 81,84,86,87,88,89,90
test cut with SpeI and BbsI

(Marisa) Cell Stock pGOP 47, 49, 51, and 82; pCir 4
miniprep pCir 4

(Jeffrey) Digested Bex 19 and 16

THURSDAY, 8/11

(Kami and Marisa) grew up BW 363, 390, 391; pSB1C3; and T40 for midiprep

(Kami) split WT, 2578, 2579 HEK cells into T175

(Jeffrey) Picked Colonies for 87 and 88

FRIDAY, 8/12

(Marisa) Midiprep BW 363, 390, 391, T40, pSB1C3

(Jeffrey) Mlinprepped and test 87 and 88... picked colonies for 90. Transformed interlab study plasmids

SATURDAY, 8/13

(Marisa) Midiprep BW 363, 390, 391, T40, pSB1C3

Cell stock pGOP 87

(Marisa and Rachel) Gel extracted

(Jeffrey) Mini prep and test cut 90... Retransformed interlab device 1

SUNDAY, 8/14

(Marisa and Rachel) Cell stock Interlab study: Pos 1, Neg 1, Part 2-1, Part 2-2, Pos 2, Neg 2, Part 3-1, Part 3-2

Miniprep Interlab Study and pGOP 56-3/4/5, pGOP 83-5/6/7/8/9

(Rachel) Digested pGOPs 56 and 83 with SpeI and XmaI

(Marisa) Test Cut pGOPs 56 and 83

(Rachel) Ligate

MONDAY, 8/15

(Jeffrey and Marisa) picked and grew up colonies: pGOP 57-1, 61-1/2/3, 62-1/2/3, 66-1/2/3, 77-1/2/3, 78-1/2/3, 80-1/2/3, 88-1/2/3/4/5, 90-1/2/3/4/5

(Kami and Marisa) Diluted and nanodropped

(Kami) made 6 plates for transfection. 3 with WT, 2 with stable line integrations of 2578 and 2579 each.

TUESDAY, 8/16

(Marisa) miniprepped pGOP 57-1, 61-1/2/3, 62-1/2/3, 66-1/2/3, 77-1, 78-1/2/3, 80-1/2/3, 88-1/2/3/4/5, 90-1/2/3/4/5

(Will) Digested Marisa's minipreps

(Kami and Rachel) 7:30 pm transfected HEK cells to test the multimerization and the circuit

for the circuit. Transfected 1/8 of 250ng of each component. (2 pgops made up 1/8 of the ng transfection and the extra was filled with Luis's pSB2C3 blank)

(Kami And rachel) 2 hours after transfection, we added 4OHT/Absistic acid/Rapalog to appropriate wells

WEDNESDAY, 8/17

THURSDAY, 8/18

(Kami and Rachel) prepared the transfected cells for FACS, and ran facs

(Kami) Cells stocked stabel lines of HEK 2578(m150) and 2579 (m151)

trypsonize cel

neutralize with media

spin cells down

remove media and break up pellet.

resuspend pellet in 90%FBS 10%DMSO (1ml per cryotube)

one T75 makess 2-3 cell stocks

(Marisa) Made LB Broth

inoculated T40 and BW363 in falcon tubes

turn LB Broth into LB + Carb broth

FRIDAY, 8/19

(Kami) made cam agar and broth

miniprep pGOP 2, 6, 83-11/12/13/14/15

(Marisa) made Cam plates

Cell stocked pGOP 57, 61, 80, 88, and 90

inoculated 1mL of each T40 and BW363 sampled into flasks with 150 mL LB + Carb

Digested pGOPs 83 with SpeI and BbsI for test cut

SATURDAY, 8/20

(Marisa) Midiprep three of both T40 and BW363

MONDAY, 8/22

(Marisa) Grew up pGOP 82
 (Jeffrey) Grew up pGOP 83-14
 (Jeffrey) took OD 600 of the interlab and diluted cells down to .02
 (Will) PCR'ed the pCMV fro the iGEM registry
 (Will) Gibsoned and transformed the CMV PCR product into pGEX
 (Marisa) made 96 well plate for InterLab Study
 250 uL DPBS and 2uL cells per sample
 (Kami Rachel) prepared HEK plates for transfection at 200,000 cells

TUESDAY, 8/23

(Marisa) Cell stock pGOP 83
 Grew up BW 1945, 1946, 1947, and 1948
 (Will) Miniprep pGOP 82 ad 83
 (Jeffrey) Picked colonies from GFP and BFP under igem CMV
 (Jeffrey) Transformed and plated interlab bacteria
 (Kami Rachel Marisa) made trasfection DNA mixes and transfected the HEK cells we plated yesterday

IMG_20160823_145414238.jpg

Multimerization				X13	
construct	mass in 2 well	vol in 4 well	vol in 4 cells	MM1 (control)	MM2 (exp)
C3A2	62.5	1.25	5	65	65
T40	62.5	1.25	5	65	65
pGOP+363	71.5+55	0.15+1.1	0.6+4.4		57.2 363
pGEX/363	62.5	1.25	5	60.2 363	65 pGEX15
				19.4 uL each well + 0.6 pGOP	19.4 uL each well + 0.6 pGOP

IMG_20160823_145437000_HDR.jpg

0123116	3 induc Cre+Flp	1 induc Cre+Flp
Const Cre+Flp	+ 842 LT	+ 842 LT
MM3	MM4	MM5
15 uL 474	15 uL 474	15 uL 474
45 uL 363	15 uL T40	15 uL 2578
10 uL each well	15 uL DC56	60 uL 363
+ 2.5 uL 390/363	15 2286	15 uL each well
+ 2.5 uL 391/363	15 2287	+ 5 uL 842/363
+ 5 uL 842/363	15 363	
	15 uL each well	3 induc Cre+Flp
	+ 5 uL 842/363	+ cr4 comb0
		LT
1 induc Cre+Flp	const Cre+Flp	MM8
+ 842 LT	+ cr4 comb0	
MM6	WT	
15 uL 474	MM7	15 uL 474
15 2579	15 uL 474	15 uL T40
60 363	15 uL T40	15 uL DC56
15 uL each well	15 uL 363	15 uL 2286
+ 5 uL 842/363	3.6 pGOP45	15 uL 2287
	3.6 pGOP23	3.6 pGOP45
	3.6 pGOP10	3.6 pGOP23
	3.6 pGOP75	3.6 pGOP10
	15.6 Luis's blank	3.6 pGOP75
		15.6 Luis's blank
	12.5 uL each well	
	+ 2.5 uL 390/363	17.5 uL each well
	+ 2.5 uL 391/363	+ 2.5 uL cr4/363
	+ 2.5 uL cr4/363	



<p>Induc Cre+Flp +Cir 4 comb 0 UT MM 9</p> <p>15 μL 474 15 μL T40 15 μL 2578 30 μL 363 3.6 μL p60P 45 3.6 μL p60P 25 3.6 μL p60P 10 3.6 μL p60P 75 15.6 μL Luv's blank 17.5 μL each well +2.5 μL Cir 4/363</p>	<p>Induc Cre+Flp +Cir 4 comb 0 UT MM 9</p> <p>15 μL 474 15 μL T40 15 μL 2578 30 μL 363 3.6 μL p60P 45 3.6 μL p60P 25 3.6 μL p60P 10 3.6 μL p60P 75 15.6 μL Luv's blank 17.5 μL each well +2.5 μL Cir 4/363</p>	<p>Const Cre+Flp +Cir 4 comb 1 UT MM 11</p> <p>15 μL 474 15 μL T40 3.6 μL p60P 25 3.6 μL p60P 5 3.6 μL p60P 45 3.6 μL p60P 43 3.6 μL p60P 3 3.6 μL p60P 63 3.6 μL p60P 50 3.6 μL p60P 30 3.6 μL p60P 70 3.6 μL p60P 55 3.6 μL p60P 35 3.6 μL p60P 15 1.8 μL Luv's blank 12.5 μL each well +2.5 μL 390/363 +2.5 μL 391/363 +2.5 μL Cir 4/363</p>
<p>Const Cre+Flp +Cir 4 comb 2 UT MM 12</p> <p>15 μL 474 15 μL T40 15 μL 363 3.6 μL p60P 55 3.6 μL p60P 30 3.6 μL p60P 3 3.6 μL p60P 45 15.6 μL Luv's blank 12.5 μL each well +2.5 μL 390/363 +2.5 μL 391/363 +2.5 μL Cir 4/363</p>	<p>Const Cre+Flp +Cir 4 comb 3 UT MM 13</p> <p>15 μL 474 15 μL T40 3.6 μL p60P 45 3.6 μL p60P 45 3.6 μL p60P 23 3.6 μL p60P 50 3.6 μL p60P 30 3.6 μL p60P 10 3.6 μL p60P 55 3.6 μL p60P 35 3.6 μL p60P 15 3.6 μL p60P 75 9 μL Luv's blank 12.5 μL each well +2.5 μL 390/363 +2.5 μL 391/363 +2.5 μL Cir 4/363</p>	



<p>Const Cre+Flp +Cir 4 comb 4 1/2 UT MM 14</p> <p>15 μL 474 15 μL T40 15 μL 363 3.6 μL p60P 23 3.6 μL p60P 30 3.6 μL p60P 75 19.2 μL Luv's blank 12.5 μL each well +2.5 μL 390/363 +2.5 μL 391/363 +2.5 μL Cir 4/363</p>	<p>Const Cre+Flp +Cir 4 comb 5 UT MM 15</p> <p>15 μL 474 15 μL T40 15 μL 363 3.6 μL p60P 23 3.6 μL p60P 30 3.6 μL p60P 70 19.2 μL Luv's blank 12.5 μL each well +2.5 μL 390/363 +2.5 μL 391/363 +2.5 μL Cir 4/363</p>	<p>Induc Cre+Flp +84d Luv 162 MM 16</p> <p>30 μL 474 150 μL 363 15 μL each well +5 μL 84d/363</p>
<p>Induc Cre+Flp +Cir 4 comb 0 Luv 1+2 MM 17</p> <p>30 μL 474 30 μL T40 90 μL 363 7.2 μL p60P 45 7.2 μL p60P 23 7.2 μL p60P 10 7.2 μL p60P 75 36.8 μL Luv's blank 17.5 μL each well +2.5 μL Cir 4/363</p>		

PLATE 1							
1	2	3	4	5	6	7	8
CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2
T40	T40	T40	T40	T40	T40	T40	T40
pGFP35	pGFP81	pGFP82	pGFP83	pGFP84	pGFP85	pGFP86	pGFP87
blank	blank	blank	blank	blank	blank	blank	blank
CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	CJH2
T40	T40	T40	T40	T40	T40	T40	T40
pGFP88	pGFP89	pGFP90	pGFP35	pGFP81	pGFP82	pGFP83	pGFP84
blank	blank	blank	GEX15	GEX15	GEX15	GEX15	GEX15

PLATE 2							
1	2	3	4	5	6	7	8
CJH2	CJH2	CJH2	CJH2	CJH2	CJH2	363	GFP
T40	T40	T40	T40	T40	T40		
pGFP85	pGFP86	pGFP87	pGFP88	pGFP89	pGFP90		
GEX15	GEX15	GEX15	GEX15	GEX15	GEX15		
BFP	mRuby	IRFP	Orange	Orange	Orange	Orange	Orange
				363	363	363	363
				842	842	842	842
				363	390	391	390+391

PLATE 3							
1	2	3	4	5	6	7	8
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
T40	T40	T40	T40	T40	T40	T40	T40
363	842	842	842	842	363	363	363
390	DC56	DC56	DC56	DC56	DC56	DC56	DC56
	2286 & 2287	2286 & 2287	2286 & 2287	2286 & 2287	2286 & 2287	2286 & 2287	2286 & 2287
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
363	363	363	363	363	363	363	363
842	842	842	842	842	842	842	842
2578	2578	blank	2578	2578	2578	2578	blank

PLATE 4							
1	2	3	4	5	6	7	8
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
T40	T40	T40	T40	T40	T40	T40	T40
363	390	391	390+391	363	pCIR4	pCIR4	pCIR4
DC56	DC56	DC56	DC56	DC56	DC56	DC56	DC56
2286&2287	2286&2287	2286&2287	2286&2287	2286&2287	2286&2287	2286&2287	2286&2287
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
T40	T40	T40	T40	T40	T40	T40	T40
pCIR4	pCIR4	pCIR4	pCIR4	pCIR4	pCIR4	pCIR4	pCIR4
DC56	DC56	DC56	DC56	DC56	DC56	DC56	DC56
2286&2287	2286&2287	2286&2287	2286&2287	2286&2287	2286&2287	2286&2287	2286&2287

PLATE 5							
1	2	3	4	5	6	7	8
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
T40	T40	T40	T40	T40	T40	T40	T40
363	363	363	363	363	363	363	363
2578	2578	2578	2578	2578	2578	2578	2578
390	390	390	390	390	390	390	390
391	391	391	391	391	391	391	391
390+391	390+391	390+391	390+391	390+391	390+391	390+391	390+391

PLATE 6							
1	2	3	4	5	6	7	8
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
T40	T40	T40	T40	T40	T40	T40	T40
363	363	363	363	363	363	363	363
842	842	842	842	842	842	842	842
pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3
363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)
GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23
GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75

PLATE 7 2578							
1	2	3	4	5	6	7	8
Orange	Orange	Orange	Orange	Orange	Orange	Orange	Orange
T40	T40	T40	T40	T40	T40	T40	T40
363	363	363	363	363	363	363	363
842	842	842	842	842	842	842	842
pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3	pSB1C3
363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)	363(for Cre Flip)
GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23	GOP25 GOP23
GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75	GOP10 GOP75

Added drugs to inducible cells 2hrs after transfection

Drugs are located in pink box labeled "Small molecules"

Rapalog= AC heterodimerizer (1000x)

Absistic Acid (1000x)- vortex before use to get rid of precipitant

4OHT (1000x)

WEDNESDAY, 8/24

(Marisa) Miniprep and nanodropped pCIR 4, pBex 16 pBEX 19, CMV BFP-1, CMV GFP-1/2, 1945, 1946, 1947, 1948
Grew up Bex 16 and Bex 19
(Jeffrey) Digested 1945, 1946, 1947, and 1948 with BbsI
Grew up pPV 1/5/6/7/8, 1942, 1943, 1944, and 1716
(Marisa) Grew up pCIR 4 (x2), DC 56, BW 390, and BW 391
Made LB + Carb

THURSDAY, 8/25

(Rachel and Kami) Ran FACS
(Kami) split HEK293 and stable line cells into big flasks
(Marisa) miniprep pPV 1, 5, 6, 7, 8, pBEX16, 19(x3), BW1942, BW1943, BW1944, BW1716
midiprep DC56, BW390, BW391, pCIR4(x2)
(Will) cut for gel extraction 1916, 1945, 1947, 1948 with BbsI
gel extracted those as well
(Will) grew up iGEM CMV for cell stock
(Rachel) annealed the mutated oligos

FRIDAY, 8/26

(Jeffrey and Marisa) Midiprep BW 2286, BW 2287, and pCIR 4 (x2)
Cell stock CMV GFP-1
(Rachel) picked two colonies for pGOP 101-150
(Will) Gel extracted pBEX 16 and 19
pBEX 16 digested with Ascl and NheI
pBEX 19 digested with Ascl and NotI
PCR overhang extension with sequences 37 and 56

SATURDAY, 8/27

(Kami, Will, Rachel) Performed minipreps for Mismatches and additional multimerized operators pGOP 101-150

MONDAY, 8/29

(Will) Test Cut for Mismatches
(Jeffrey) Attempted overhang extension PCR ... Failed

TUESDAY, 8/30

(Jeffrey) Attempted overhang extension PCR ... Failed
(Jeffrey) Completed Interlab Study
(Marisa) PCR GEX plasmids for CAM transfer
(Will) Gel extracted PCR

WEDNESDAY, 8/31

(Marisa) Digested pGOP 109, 110, 112, 114, 134, 135, 136, 138
(Marisa) Grew up pGEX 7, pGEX 8, pGPX1,
(Jeffrey) Attempted overhang extension PCR ... Failed
(Jeffrey) Test Cut 109, 110, 112, 114, 134, 135, 136, 138

THURSDAY, 9/1

(will) pcr mruby and bfp to drop into pgops
(Jeffrey) Miniprep pgor 80-90 for multiplecolors , bexs for circuits, gexs f
(Marisa) Passaged HEK cells
Made 2L of LB Broth, 2L of LB+Carb broth, and 2L of LB + Cam agar for plates
Ligated and transformed pGEX 3, 5, 8, 10, 15, and 19 into Cam (with BW 1721 -n330 backbone)
Grew up T4o (x2), BW 363, GEX 5, GEX 15, pCIR 4 (x2), and GPX 2
(Rachel) autoclaved epi tubes and tips
Ligated and transformed BW 1945-1949

FRIDAY, 9/2

- (Marisa) Made LB + Cam plates
- (Jeffrey) Resuspended cells grown for midiprep and stored in the cold room
- (Rachel) Prepared and ran cells for FACS
- (Will) Annealed oligos for multimerized pGOP 3

SATURDAY, 9/3

- (Jeffrey) transformed mismatched pGOP 101-150
- (Marisa and Jeffrey) grew up GEX 3, 5, 8, 10, 15, 19, BW 1721, 1945, 1946, 1947, 1948, and GOP 29

SUNDAY, 9/4

- (Jeffrey) picked and grew up colonies from pGOP101-150
 - miniprep of GEX 3, 5, 8, 10, 15, 19, BW 1721, 1945, 1946, 1947, 1948, and GOP 29
 - Digested BW1945, 1946, 1947, 1948
- (Marisa) Cell stock GOP 110, 112, and 114
 - Digested GOP 29(MOP) with EcoRI and PstI
 - PCR GEX 3, 5, 8, 10, 15, 19 with PrA6 and PrA6r
 - PCR clean up BW 1945, 1946, 1947, 1948 and MOP
 - Made and ran gel for GEX 3, 5, 8, 10, 15, 19 and BW1721
 - (with Kami) gel extracted GEX 3, 5, 8, 10, 15, 19, and BW 1721
 - Digested GEX 3, 5, 8, 10, 15, 19 and BW 1721 with EcoRI and PstI
 - PCR clean up GEX 3, 5, 8, 10, 15, 19, and BW 1721

MONDAY, 9/5

- (Marisa) Cell stock GOP 101-103; 105-108; 111; 115-134; and 146
 - Miniprep pGOP 101-103; 105;108; 111; and 115-122
 - Ligate and transformed pGEX 3, 5, 8, 10, 15, and 19 into pSB1C3 backbone
- (Will) Miniprep pGOP 123-134 and 146
- (Kami Rachel) seeded 8 48 well plates at 100,000cells for transfection on wednesday

TUESDAY, 9/6

- (Marisa) picked two colonies of GEX 3, 5, 8, 10, 15, and 19 (now termed 103, 105, 108, 110, 115, and 119) from transformation 9/5
- (Jeffrey) Performed overhang extension pcr to build the decoder for circuit 2
- (Jeffrey) transformed the 18T entry vectors for circuit 5

WEDNESDAY, 9/7

- (Marisa) miniprep GEXs 103, 105, 108, 110, 115, and 119
 - Test cut with AseI
 - picked two colonies for 18T-terminator circuits
- (Rachel, Marisa and Kami) transfected HEKcells

THURSDAY, 9/8

- (kami Marisa) passaged HEK cells
- (jeffrey) performed a gibson to build circuit 2
- (jeffrey) test cuts for 18T entry vecotrs for circuit 5 indicated low digestion efficeincy and the circuits needs to be retransformed

FRIDAY 9/9

- (Will) Ligated and transformed GEX 21-40 and GOP 151-168
- (Will and Marisa) Grew up 3 colonies for GEX 21-40 and GOP 151-168

Saturday 9/10

- (Will Jeffrey Marisa) Miniprepped GEX 21-40 and GOP 151-168
- (Marisa) PCR GEX 3, 5, 8, 10, 15, 19, and BW 1721 with PrA6 and PrA6r (previous GEX 103, 105 108, 110, 115, 119 were not successful)
- (Marisa and Will) gel extracted GEX 3, 5, 8, 10, 15, and 19 and BW 1721 (~400); will be insert for GEX 103, 105, 108, 110, 115, and 119

SUNDAY, 9/11

(Marisa) Digested GEX 103, 105, 108, 110, 115, and 119 and BW 1721 inserts with EcoRI-HF and PstI-HF

(Jeffrey) Picked colonies for 18T entry vectors and miniprep circuit 2 for test cuts

(Will) Test cut GEX 21-40 and GOP 151-168

MONDAY, 9/12

(Marisa) picked four colonies for 18T g1, g17, g3, g8, and g13

Digested and test-cut 18T g1, g17, g3, g8, and g13 colonies 1 and 2 with BbsI and SpeI

(Marisa) Ligated and transformed inserts for GEX 3, 5, 8, 10, 15, and 19 into pSB1C3 backbone to create GEX 103, 105, 108, 110, 115, and 119

(Kami) made 8 48 well plates with HEK cells at 100,000 cells/ml

(kami) retransformed pGOP7 because it looked wierd when we transfected it. so going to pick one colony from this transformation to try to start fresh.

(Marisa) Made LB + Cam broth

Grew up pSB1C3, GOP 35, and GPX2 for midipreps

Grew up pSB1C3 and GOP 45 (x2) for minipreps

TUESDAY, 9/13

(Marisa) picked two colonies for GEX 103, 105, 108, 110, 115, and 119

Midiprep pSB1C3, pGPX2, pGOP35

Miniprep GOP45 (x2), SB1C3, and g17 18T-1/2/3/4

made master mixes for transfection to take place on 9/14

picked and grew up a colony from GOP7 re-transformation

Added Carbacillin to plain LB Broth to make LB+Carb broth

Grew up Cir4 (x2), SB1C3 (x2), GEX5, GEX 15, GPX2, T40 (x2), and BW471

WEDNESDAY, 9/14

(Marisa) miniprepped pGEX 103-1/2, 105-1/2, 108-1/2, 110-1/2, 115-1/2, and 119-1/2

Test cut GEXs with EcoRI-HF and PstI-HF

Test Cut g17 18T with BbsI and SpeI

(Marisa and Jeffrey) made 2L of plain LB broth

(Jeffrey) midiprep pSB1C3 (x2), T40(x2), and pCir4(x2)

(Kami Rachel) transfected Hek cells. with pGOPs diluted to 5ng/ul (not 50ng/ul) and add a bigger volume to the transfection mix.

Screen Shot 2016-09-15 at 01.43.39.png

PLATE 1								
	1	2	3	4	5	6	7	8
A	BW 471	BW 471	BW 471	BW 471	BW 471	BW 471	BW 471	BW 471
B	pGOP 21	pGOP 22	pGOP 23	pGOP 24	pGOP 25	pGOP 26	pGOP 27	pGOP 28
C	BW 363	BW 363	BW 363	BW 363	BW 363	BW 363	BW 363	BW 363
D	T40	T40	T40	T40	T40	T40	T40	T40
E	BW 471	BW 471	BW 471	BW 471	BW 471	BW 471	BW 471	BW 471
F	pGOP 29	pGOP 30	pGOP 31	pGOP 32	pGOP 33	pGOP 34	pGOP 35	pGOP 36
	BW 363	BW 363	BW 363	BW 363	BW 363	BW 363	BW 363	BW 363
	T40	T40	T40	T40	T40	T40	T40	T40


PLATE 2								
	1	2	3	4	5	6	7	8
A	BW 471 pGOP 37	BW 471 pGOP 38	BW 471 pGOP 39	BW 471 pGOP 21	BW 471 pGOP 22	BW 471 pGOP 23	BW 471 pGOP 24	BW 471 pGOP 25
B	BW 363 T40	BW 363 T40	BW 363 T40	pGEX 1 T40	pGEX 2 T40	pGEX 3 T40	pGEX 4 T40	pGEX 5 T40
C								
D	BW 471 pGOP 26	BW 471 pGOP 27	BW 471 pGOP 28	BW 471 pGOP 29	BW 471 pGOP 30	BW 471 pGOP 31	BW 471 pGOP 32	BW 471 pGOP 33
E	pGEX 6 T40	pGEX 7 T40	pGEX 8 T40	pGEX 9 T40	pGEX 10 T40	pGEX 11 T40	pGEX 12 T40	pGEX 13 T40
F								

PLATE 4								
	1	2	3	4	5	6	7	8
A	CJH2 pGOP 15	CJH2 pGOP 19	BW 471 pGOP 43	BW 471 pGOP 45	BW 471 pGOP 48	BW 471 pGOP 50	BW 471 pGOP 55	BW 471 pGOP 59
B	pGEX 15 T40	pGEX 19 T40	BW 363 T40	BW 363 T40	BW 363 T40	BW 363 T40	BW 363 T40	BW 363 T40
C								
D	BW 471 pGOP 43	BW 471 pGOP 45	BW 471 pGOP 48	BW 471 pGOP 50	BW 471 pGOP 55	BW 471 pGOP 59	BW 471 pGOP 63	BW 471 pGOP 65
E	pGEX 3 T40	pGEX 5 T40	pGEX 8 T40	pGEX 10 T40	pGEX 15 T40	pGEX 19 T40	BW 363 T40	BW 363 T40
F								

PLATE 5								
	1	2	3	4	5	6	7	8
A	BW 471 pGOP 68	BW 471 pGOP 70	BW 471 pGOP 75	BW 471 pGOP 79	BW 471 pGOP 83	BW 471 pGOP 85	BW 471 pGOP 88	BW 471 pGOP 90
B	BW 363 T40	BW 363 T40	BW 363 T40	BW 363 T40	pGEX 3 T40	pGEX 5 T40	pGEX 8 T40	pGEX 10 T40
C								
D	BW 471 pGOP 75	BW 471 pGOP 79	BW 471 pGOP 83	BW 471 pGOP 85	BW 471 pGOP 88	BW 471 pGOP 90	BW 471 pGOP 93	BW 471 pGOP 95
E	pGEX 15 T40	pGEX 19 T40	pGEX 23 T40	pGEX 27 T40	pGEX 31 T40	pGEX 35 T40	pGEX 39 T40	pGEX 43 T40
F								

PLATE 6								
	1	2	3	4	5	6	7	8
A	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40
B	pCir4 390	pCir4 391	pCir4 390 + 391	pCir4 390	pCir4 391	pCir4 390	pCir4 391	pCir4 390 + 391
C	pGOP 25, 3, 70, 55	pGOP 25, 3, 70, 55	pGOP 25, 3, 70, 55	pGOP 25, 3, 70, 55	pGOP 25, 3, 70, 55	pGOP 25, 3, 70, 55	pGOP 25, 3, 70, 55	pGOP 25, 3, 70, 55
D	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40	BW 474 T40
E	BW 363 pGOP 5,	BW 363 pGOP 5,	BW 363 pGOP 5,	BW 363 pGOP 5,	BW 363 pGOP 5,	BW 363 pGOP 5,	BW 363 pGOP 5,	BW 363 pGOP 5,
F	63, 50, 35	43, 30, 15	43, 30, 15	43, 30, 15	43, 30, 15	43, 30, 15	43, 30, 15	43, 30, 15

PLATE 7								
	1	2	3	4	5	6	7	8
A	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40
B	GEX44 pGOP35	GEX45 pGOP35	GEX46 pGOP35	GEX47 pGOP35	GEX48 pGOP35	GEX49 pGOP35	GEX50 pGOP35	GEX51 pGOP35
C								
D	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40	BW471 T40
E	GEX52 pGOP35	GEX53 pGOP35	GEX54 pGOP35	GEX55 pGOP35	GEX56 pGOP35	GEX57 pGOP35	GEX58 pGOP35	GEX59 pGOP35
F								



THURSDAY, 9/15

(Marisa) Miniprep GEX 103-3/4/5/6, 105-3/4/5/6, 108-3/4/5/6, 110-3/4/5/6, 115-3/4/5/6, 119-3/4/5/6

Digest and test cut GEXs with EcoRI-HF and PstI-HF

(Marisa and Kami) passaged stable line HEK cells (Wild Type, 2578, and 2579) 1:20

(Jeffrey) Minipreped 1945-1948

FRIDAY, 9/16

(Kami) Miniprep pGOP 159-1/2/3, 160-1/2/3, 169-1/2/3, and 170-1/2/3

(Will) Test cut pGOP 159, 160, 169, and 170

SATURDAY, 9/17

(jeffrey) ligated and transformed BEX34-38

SUNDAY, 9/18

(Marisa and Jeffrey) picked three colonies of GOP 151-158 and 161-168

(Jeffrey) Picked colonies for BEX37 Redigested 1945,1946,1948

(Jeffrey) Religated and transformed BEX 34,35,36,38

MONDAY, 9/19

(Marisa) Miniprep colonies 1,2, and 3 of GOP 151-158 and 161-168

(Jeffrey) Picked colonies for BEX 34-38

(Kami) Made 6 hek plates at 100,000 (NOT USED)

(Will) Test cut the minipreped GOPs 151-168

TUESDAY, 9/20

(Marisa) Cell stock BEX 30, GOP 159, 160, 169

(Marisa) made 2L of plain LB broth

(Marisa) Diluted GOP 81-90 to 25 ng/ul

(jeffrey) Miniprepred bex 34-38 and test cut

(Jeffrey) Digested BEX 1,6,15,20, 30,31,32,33

(Jeffrey) digested parts vectors with bbsi for 3 input circuit

(Will) Sent in colony 1 of GOPs 151-168 that test cut correctly

WEDNESDAY, 9/21

(Marisa and Rachel) Made LB agar (x2)

(Marisa) Cell stock pGOP 159, 160, and 169

(Marisa and Kami) Ligated GOP 153, 154, 157, 163, 165, 166, 168, 170, 162, and 167 with respective oigos

(rachel) transformed the pGOPs ligated above

(Kami) Made LB + Cam plates

(Marisa) Made LB + Carb plates

(Jeffrey) Gibson Cir2 and Cir 5 using BEX parts vectors

THURSDAY, 9/22

(Kami) picked 4 colonies for GOP 153, 154, 157, 162, 163, 165, 166, 167, 168, 170

(Jeffrey) Cell stock GOP 151, 152, 155, 156, 158, 161, 164

Miniprep colonies 5-8 of GOP 153, 154, 157, 162, 163, 165, 166, 167, 168, 170

Test cut GOPs with BbsI and SpeI

Digested BEX 34-37 with Ascl and NotI

Digested BEX 38 with Ascl and NheI

(Kami) passaged HEK cells

(Will) Sent in colony two of GOPs 151-168 that did not sequence properly

FRIDAY, 9/23

Cell stocked pGOP 154,

SUNDAY, 9/25

(Rachel) Ligated and transformed 1721 insert into MOP (GOP with biobrick cut out)

(Jeffrey) Gibson Circuit 2 and 5

MONDAY, 9/26

(Marisa) miniprep colonies of pGOP 104, 109, 135, 137, 144, 149, and 150

Test cut GOP (x3) of 91-100, GOP (x4) 104, 109, 135, 137, 144, 149, 150, 153, 154, 157, 162, 163, 165, 166, 167, 168, 170, and GOP (x3) 113, 135, 136, 138, 142, with BbsI and SpeI

(Marisa) Transformed Circuit 2 and Circuit 5

(Marisa) picked four colonies of MOP + 1721

(Will) Sent in GOP151 and 163 for sequencing

(Kami) made HEK plates

TUESDAY, 9/27

(Marisa) PCR pGEX 3, 8, 15, and 19 and BW 1721

(Marisa and Will) Gel extracted pGEX 3, 8, 15, 19 and BW 1721

(Will) Digested GEX 3, 8, 15, 19 and BW 1721 with EcoRI-HF and PstI-HF

PCR cleanup of GEXs and BW1721

(Rachel) Ligated GEX 3, 8, 15, 19, and BW 1721 with MOP

(Marisa) Transformed GEX 3, 8, 15, 19 and BW 1721 (to GEX 103, 108, 115, and 119)

(Will) Cell stocked pGOP 151 and 163

WEDNESDAY, 9/28

(Marisa) Test cut GOP 162, 163, and 167 colonies 5-8 with BbsI and SpeI

(Marisa) Cell stock 91-100, 104, 109, 113, 136, 137, 138, 142, 144, 149, 150, 153, and 157

(Marisa) miniprep GOP 112, 114, 135,

(Marisa) Replated Circuit 2 and 5

THURSDAY, 9/29

(Marisa) Cell stock pGOP 151 and 163

(Marisa) Replated Circuit 2 and Circuit 5

(Marisa and Rachel) ligated and transformed pGOP 142, 162, and 167

FRIDAY, 9/30

(Marisa) picked 5 colonies of pGOP 167

(Kami) made 11 hek plates at 100,000cells/ul

(kami)Diluted pGOPs to 25ng/ul

SATURDAY, 10/1

(Marisa) Miniprep and test cut pGOP 167-1/2/3/4/5 with BbsI and SpeI

SUNDAY, 10/2

(Kami Rachel) Diluted DNA for transfection.

(Kami Rachel) Transfected hek cells with all the multimerized and mutated pGOPs



PLATE 1								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP25 363	471 T40 pGOP91 363	471 T40 pGOP92 363	471 T40 pGOP93 363	471 T40 pGOP94 363	471 T40 pGOP95 363	471 T40 pGOP96 363	471 T40 pGOP97 363
B								
C								
D	471 T40 pGOP98 363	471 T40 pGOP99 363	471 T40 pGOP100 363	471 T40 pGOP25 GEX5 (g3)	471 T40 pGOP91 GEX5	471 T40 pGOP92 GEX5	471 T40 pGOP93 GEX5	471 T40 pGOP94 GEX5
E								
F								

PLATE 2								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP95 GEX5	471 T40 pGOP96 GEX5	471 T40 pGOP97 GEX5	471 T40 pGOP98 GEX5	471 T40 pGOP99 GEX5	471 T40 pGOP100 GEX5	471 T40 pGOP30 363	471 T40 pGOP151 363
B								
C								
D	471 T40 pGOP152 363	471 T40 pGOP153 363	471 T40 pGOP154 363	471 T40 pGOP155 363	471 T40 pGOP156 363	471 T40 pGOP157 363	471 T40 pGOP158 363	471 T40 pGOP159 363
E								
F								

PLATE 3								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP160 363	471 T40 pGOP30 GEX10 (g8)	471 T40 pGOP151 GEX10	471 T40 pGOP152 GEX10	471 T40 pGOP153 GEX10	471 T40 pGOP154 GEX10	471 T40 pGOP155 GEX10	471 T40 pGOP156 GEX10
B								
C								
D	471 T40 pGOP157 GEX10	471 T40 pGOP158 GEX10	471 T40 pGOP159 GEX10	471 T40 pGOP160 GEX10	471 T40 pGOP23 363	471 T40 pGOP161 363	471 T40 pGOP162 363	471 T40 pGOP163 363
E								
F								



PLATE 4								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP164 363	471 T40 pGOP165 363	471 T40 pGOP166 363	471 T40 pGOP167- 1 363	471 T40 pGOP168 363	471 T40 pGOP169 363	471 T40 pGOP170 363	471 T40 pGOP23 GEX3 (gt1)
B								
C								
D	471 T40 pGOP161 GEX3	471 T40 pGOP162 GEX3	471 T40 pGOP163 GEX3	471 T40 pGOP164 GEX3	471 T40 pGOP165 GEX3	471 T40 pGOP166 GEX3	471 T40 pGOP167- 1 GEX3	471 T40 pGOP168 GEX3
E								
F								

PLATE 5								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP169 GEX3	471 T40 pGOP170 GEX3	471 T40 pGOP101 363	471 T40 pGOP102 363	471 T40 pGOP103 363	471 T40 pGOP104 363	471 T40 pGOP105 363	471 T40 pGOP106 363
B								
C								
D	471 T40 pGOP107 363	471 T40 pGOP108 363	471 T40 pGOP109 363	471 T40 pGOP110 363	471 T40 pGOP111 363	471 T40 pGOP112 363	471 T40 pGOP113 363	471 T40 pGOP114 363
E								
F								

PLATE 6								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP115 363	471 T40 pGOP116 363	471 T40 pGOP117 363	471 T40 pGOP118 363	471 T40 pGOP119 363	471 T40 pGOP120 363	471 T40 pGOP101 GEX15 (g13)	471 T40 pGOP102 GEX15
B								
C								
D	471 T40 pGOP103 GEX15	471 T40 pGOP104 GEX15	471 T40 pGOP105 GEX15	471 T40 pGOP106 GEX15	471 T40 pGOP107 GEX15	471 T40 pGOP108 GEX15	471 T40 pGOP109 GEX15	471 T40 pGOP110 GEX15
E								
F								

PLATE 7								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP111 GEX15	471 T40 pGOP112 GEX15	471 T40 pGOP113 GEX15	471 T40 pGOP114 GEX15	471 T40 pGOP115 GEX15	471 T40 pGOP116 GEX15	471 T40 pGOP117 GEX15	471 T40 pGOP118 GEX15
B								
C								
D	471 T40 pGOP119 GEX15	471 T40 pGOP120 GEX15	471 T40 pGOP141 363	471 T40 pGOP142 363	471 T40 pGOP143 363	471 T40 pGOP144 363	471 T40 pGOP145 363	471 T40 pGOP141 GEX15 (g13)
E								
F								

PLATE 8								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP142 GEX15	471 T40 pGOP143 GEX15	471 T40 pGOP144 GEX15	471 T40 pGOP145 GEX15	471 T40 pGOP35 363	471 T40 pGOP81 363	471 T40 pGOP82 363	471 T40 pGOP83 363
B								
C								
D	471 T40 pGOP84 363	471 T40 pGOP85 363	471 T40 pGOP86 363	471 T40 pGOP87 363	471 T40 pGOP88 363	471 T40 pGOP89 363	471 T40 pGOP90 363	471 T40 pGOP35 GEX15 (g13)
E								
F								

PLATE 9								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP81 GEX15	471 T40 pGOP82 GEX15	471 T40 pGOP83 GEX15	471 T40 pGOP84 GEX15	471 T40 pGOP85 GEX15	471 T40 pGOP86 GEX15	471 T40 pGOP87 GEX15	471 T40 pGOP88 GEX15
B								
C								
D	471 T40 pGOP89 GEX15	471 T40 pGOP90 GEX15	363					
E								
F								

PLATE 10								
	1	2	3	4	5	6	7	8
A	471 T40 pGOP167- 2 363	471 T40 pGOP167- 3 363	471 T40 pGOP167- 4 363	471 T40 pGOP167- 5 363	471 T40 pGOP167- 2 GEX3	471 T40 pGOP167- 3 GEX3	471 T40 pGOP167- 4 GEX3	471 T40 pGOP167- 5 GEX3
B								
C								
D								
E								
F								

FRIDAY, 10/7

(kami) made 7 hek plates at 100,000 for transfection on sunday and 1 hek plate at 200,000 cells/ml for transfection on saturday with WPI

(Rachel) prepared the transfection DNA mix for WPI

PLATE 1								
	1	2	3	4	5	6	7	8
A	mRUBY1000 marker dCASS-VPR Recombination circuit	mRUBY1000 marker dCASS-VPR Recombination circuit	mRUBY1000 marker dCASS-VPR Recombination circuit	mRUBY1000 marker dCASS-VPR Recombination circuit				
B	Blank pGOF25 (OFF) pGOF55 (BFP)	Cm pGOF25 (OFF) pGOF55 (BFP)	Pig pGOF25 (OFF) pGOF55 (BFP)	Cm + Pig pGOF25 (OFF) pGOF55 (BFP)	Wildtype			
C								
D								
E								
F								

SATURDAY, 10/8

(Kami + Rachel) transfected cells for WPI testing the and/nor circuit

(Kami + Rachel) prepared all the DNA mixes for transfection tomorrow

SUNDAY, 10/9

(Kami) Transfected HEK cells with screen of mutated p GOPs, CMV compared to single and triple multimerized, and 4 analog circuits

Screen Shot 2016-10-10 at 11.32.13.png

PLATE 1								
	1	2	3	4	5	6	7	8
A	BW 471 T40 363 pCir4 pGOP 170 pGOP160 pGOP90	BW 471 T40 390 pCir4 pGOP 170 pGOP160 pGOP90	BW 471 T40 391 pCir4 pGOP 170 pGOP160 pGOP90	BW 471 T40 390 + 391 pCir4 pGOP 170 pGOP160 pGOP90	BW 471 T40 363 pCir4 pGOP23 pGOP163 pGOP90	BW 471 T40 390 pCir4 pGOP23 pGOP163 pGOP90	BW 471 T40 391 pCir4 pGOP23 pGOP163 pGOP90	BW 471 T40 390 + 391 pCir4 pGOP23 pGOP163 pGOP90
B								
C								
D	BW 471 T40 363 pCir4 pGOP100 pGOP23 pGOP160 pGOP110	BW 471 T40 390 pCir4 pGOP100 pGOP23 pGOP160 pGOP110	BW 471 T40 391 pCir4 pGOP100 pGOP23 pGOP160 pGOP110	BW 471 T40 390 + 391 pCir4 pGOP100 pGOP23 pGOP160 pGOP110	BW 471 T40 363 pCir4 pGOP94 pGOP169 pGOP160 pGOP90	BW 471 T40 390 pCir4 pGOP94 pGOP169 pGOP160 pGOP90	BW 471 T40 391 pCir4 pGOP94 pGOP169 pGOP160 pGOP90	BW 471 T40 390 + 391 pCir4 pGOP94 pGOP169 pGOP160 pGOP90
E								
F								

PLATE 2								
	1	2	3	4	5	6	7	8
A	BW 471 T40 CMV	BW 471 T40 pGOP 23 353	BW 471 T40 pGOP 25 363	BW 471 T40 pGOP 30 363	BW 471 T40 pGOP 35 363	BW 471 T40 pGOP90 363	BW 471 T40 pGOP100 363	BW 471 T40 pGOP160 363
B								
C								
D	BW 471 T40 pGOP170 363	BW 471 T40 pGOP 23 GEX 3 (g1)	BW 471 T40 pGOP 25 GEX 5 (g3)	BW 471 T40 pGOP 30 GEX10 (g8)	BW 471 T40 pGOP 35 GEX15 (g13)	BW 471 T40 pGOP90 GEX 15 (g13)	BW 471 T40 pGOP100 GEX 5 (g3)	BW 471 T40 pGOP160 GEX10 (g8)
E								
F								

PLATE 3								
	1	2	3	4	5	6	7	8
A	BW 471 T40 pGOP170 GEX 3 (g1)	471 T40 pGOP101 363	471 T40 pGOP102 363	471 T40 pGOP103 363	471 T40 pGOP104 363	471 T40 pGOP105 363	471 T40 pGOP106 363	471 T40 pGOP107 363
B								
C								
D	471 T40 pGOP108 363	471 T40 pGOP109 363	471 T40 pGOP110 363	471 T40 pGOP111 363	471 T40 pGOP112 363	471 T40 pGOP113 363	471 T40 pGOP114 363	471 T40 pGOP115 363
E								
F								

PLATE 4								
	1	2	3	4	5	6	7	8
A	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40
B	pGOP116 363	pGOP117 363	pGOP118 363	pGOP119 363	pGOP120 363	pGOP101 GEX15 (g13)	pGOP102 GEX15	pGOP103 GEX15
C								
D	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40
E	pGOP104 GEX15	pGOP105 GEX15	pGOP106 GEX15	pGOP107 GEX15	pGOP108 GEX15	pGOP109 GEX15	pGOP110 GEX15	pGOP111 GEX15
F								

PLATE 5								
	1	2	3	4	5	6	7	8
A	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40
B	pGOP112 GEX15	pGOP113 GEX15	pGOP114 GEX15	pGOP115 GEX15	pGOP116 GEX15	pGOP117 GEX15	pGOP118 GEX15	pGOP119 GEX15
C								
D	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40	471 T40
E	pGOP120 GEX15	pGOP141 363	pGOP142 363	pGOP143 363	pGOP144 363	pGOP145 363	pGOP141 GEX15 (g13)	pGOP142 GEX15
F								

PLATE 6								
	1	2	3	4	5	6	7	8
A	471 T40	471 T40	471 T40	363				
B	pGOP143 GEX15	pGOP144 GEX15	pGOP145 GEX15					
C								
D								
E								
F								