

## Team Bax

Week	Date	Jobs
35	29.08.16	Digest of cloning vector PMZ333 with NotI and XbaI and electrophoresis to isolate backbone from insert
	30.08.16	PCR Gel Clean up from Phusion PCR product PCR Gel Clean up from cloning vector PMZ333
	31.08.16	Another Phusion PCR from same vector to generate mCherry Baxshort, because PCR was not successful Electrophoresis of PCR product, which were cut out of gel
	01.09.16	PCR Clean up from Phusion PCR product