

Team Mammal

Week	Date	Jobs
34	2016-08-22	<p>2x test digestion of lov (19.8) and pdz (18.8) midi preps.</p> <p>Overnight cultures: 10 colonies from pJET-pdz plasmid and 5 colonies from iGEM-pdz</p> <p>Retransformation of pKM022 & pMZ333 (Note: pellet retransformation has to many colonies. Use 25 µl & 50 µl in the future)</p>
	2016-08-23	Mini prep of 8 pJET-pdz and 5 iGEM-pdz overnight cultures
	2016-08-24	<p>Mini prep & test digest + electrophoresis of 3 pmz333 and 5 pkm022 colonies</p> <p>Phusion PCR of LOV and PDZ</p>
	2016-08-25	<p>Mini prep of 10 iGEM-pdz-9 colonies + concentration via evaporazation (-> c = 310.9ng / µl)</p> <p>Phusion PCR (x2) of IRES (failed) and backbone (failed)</p>
	2016-08-26	<p>Phusion PCR of IRES and backbone (failed)</p> <p>Phusion PCR (x8) of backbone</p> <p>Electrophoresis of IRES</p>