

TITLE "Good to Great"

PROJECT

Continued from page

Goals for Gold

Wiki — clean, easy to navigate, subpages, nice color scheme, thematic, graphics, animation

Presentation — charismatic, diagrams/visuals > words, limited interaction

Human Practices — multidisciplinary, multi-layered, interaction with other teams, visit schools, distribute pamphlets @ library, iGEM-sponsored tutoring, summer camp, ASF, NGF, local charities, create curriculum, social media

Poster — thematic, subtle graphics, color scheme, flow, well-organized

Science — mathematical & graphical models, repeatability, functional prototype, results, larger impacts

Meeting

- compared William & Mary and SAS Taipei
 - ↳ rubric for each part of the project
- made our emails for Gold

HW — Pre-Lab for a Transformation using Puc19 (due Thursday morning)

- Title
- Purpose
- Background → 7-9 sentences
- Materials
- Methods
- Pre-Lab ? → What would a successful transformation look like

SIGNATURE

DATE

Continued to page

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

12/10/2015

Continued from page

Purpose: The objective of a transformation is to make multiple copies of DNA, once they are manipulated. It will then be transferred into bacteria, specifically E. coli.

Background Information: In a transformation, circular pieces of DNA called plasmids are manipulated. In biotechnology, plasmids can be easily engineered and transferred into bacteria to see its effects. Sterile techniques are used in this lab to prevent contamination of the lab space and the data. Specifically for this transformation, pUC19 will serve as a plasmid cloning vector in E. coli, and this is how DNA will make several copies. After these copies are made, the cells will be grown on agar plates. The agar serves as food to fuel the cells with energy so that they can divide. LB broth helps the cells to recover and grow. The cells are incubated at 37°C because that is their optimal temperature.

Materials: NEB 10-beta E. coli cells	room temperature SOC
ice	centrifuge
micropipettes, tips	Amp plates
waste beaker	incubator
transformation tubes	
plasmid DNA	

Protocol:

- 1) Thaw a tube of NEB 10-beta E. coli cells ~~thaw~~ for 10 minutes or until the ice crystals disappear. Mix gently and pipette 50 μL into a transformation tube on ice.
- 2) Add 1-5 μL of plasmid DNA containing 1 pg-100 ng to the cell mixture. Flick 4-5 times but do not vortex. ^{↳ PUC19} Add 1 μL of dH₂O to control (no plasmid).
- 3) Place mixture on ice for 30 minutes.
- 4) Heat shock at 42°C for 30 seconds.
- 5) Place on ice for 5 minutes.
- 6) Pipette 950 μL of room temperature SOC into the mixture.
- 7) Place at 37°C for 40 minutes. Shake at 250 rpm or rotate, using a ^{centrifuge} ~~centrifuge~~.

SIGNATURE

incubator

DATE

12/17/2015

Continued to page 3

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE NEB 10-Beta High Efficiency Transformation

PROJECT

Continued from page 2

8) Warm selection plates to 37°C.

9) Flick the tube and invert to mix the cells. Make 10-fold dilutions in SOC.

10) Spread 50-100 μ L of each dilution onto an Amp plate. Incubate overnight at 37°C or incubate at 30°C for 24-36 hours or 25°C for 48 hours.

Pre-Lab Question:

1) What would a successful transformation look like?

Multiple colonies of about the same size will form. They will be distributed relatively evenly throughout the plate. When the DNA is put into gel electrophoresis, results will show that the specific sequence is present in the plasmid. **Control?**

Plating Labeling: EC 5 0

LH 5 +

NC 5 ⊕

0 - no plasmid

+ - plasmid - Plain

⊕ - plasmid - Amp

Continued to page

SIGNATURE

DATE

12/17/2015


DISCLOSED TO AND UNDERSTOOD BY

DATE


PROPRIETARY INFORMATION

TITLE Transformation Results


Continued from page



no plasmid/Amp
negative control
no growth



NEB 10
plain plate with plasmid
a lot of colonies → both with and without plasmids



NEB 10
Amp plate with plasmid
less colonies → can calculate transformation efficiency

Future Experiments

- find a single colony → inoculate it (collect part of the colony) → streak → liquid colony → etc.

- decrease/dilute to have less colonies

Continued to page

DATE

1/7/2016

READ TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE Potential Projects - Brainstorm/ PROJECT QUESTIONS

Continued from page

Zinc Biosensor Questions:

- 1) Does the biosensor detect outside of 8-15 μM ?
- 2) What do you envision about our iGEM team?
- 5 3) Would it be more helpful to engineer a pathway with standard systems rather than using a zinc biosensor?
↳ proof of concept for precision vs. dynamic engineering?
- 4) Can you provide an analogy to clarify the differences between precision and dynamic engineering?
- 10 5) Could you clarify which colors correlate to which levels of Zinc? There was an error in the press release.
- 6) Is there a kill switch?
- 7) How exclusive is the biosensor to Zinc?

15 Notes - Skype Interview

- #6 → no killswitch for the system, but precision metabolic engineering can be used to create one
- degradation tags → adding deg tags affect how quickly you can move from one state to another
- 20 • developing a biosensor →
 - need something that can sense whatever is needed to be identified
ex/ transcriptional regulatory mechanism/quorum sensing
 - challenges:
 - 25 ① sensing mechanism must be exclusive to only what you want to detect
 - ② time (time require ~24 hrs to obtain desired result)
↳ possible solution - use of catalysts, especially pigment producing enzymes
 - ideal catalyst would be one that can expedite transcription
 - need a well-designed, elegant construct
- 30 • inducer → a chemical that binds to a molecule and invokes a certain behavior
↳ standard inducible promoters → chromoproteins → ex/ TS Purple
ex/ LacI

SIGNATURE

Continued to page

DATE

1/21/2016

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

Continued from page

• Reporters → [in syn bio] mechanism to display information

ex/ If you design a biosensor system that is sensitive to substrate A, then it would be beneficial to have a way of displaying results rather than performing additional assays

• fluorescent proteins → cells fluoresce when excited by light of a certain wavelength

• luciferases → enzymes that cause a cell to catalyze a reaction that produces light

• Types of reporters →

*** ① chromoproteins

② fluorescent proteins (i.e. GFP/RFP)

③ luciferase proteins

④ enzymes that produce colored substrates

• Different ways of measuring/quantifying reporters

↳ depends on particular reporter and type of characterization data desired

① microscopy → useful for obtaining spatial and temporal information, especially at a single cell level

② flow cytometers → measuring distribution in reporter activity across a large population of cells

③ plate readers → population average measurements of many different samples over time

• Chromoproteins → responsible for coloration in corals/sea anemones

• expressed pigment can be seen with the naked eye

• collection of reporter genes

• function → quantify the level of induction and amount of protein produced
ex/ TS Purple, SOrange, ASPink

• RBS = Ribosomal Binding Site → RNA sequence found in mRNA which ribosomes can bind and initiate translation

• translation in bacteria must be initiated by both an RBS sequence and a start codon

↳ protein coding sequence begins with start codon in iGen registry

↳ if you want to build a BioBrick system that produces a protein, you need to select an RBS port and put it upstream of the protein coding sequence to be translated

SIGNATURE

Continued to page 7

DATE

2/4/2014

DISCLOSED TO AND UNDERSTOOD BY


DATE

PROPRIETARY INFORMATION

TITLE Chromoproteins/RBS/Degradation TAGS Research PROJECT

Continued from page 6

RBS Image



- Degradation Tags → short polypeptide sequence that marks a protein for degradation by the cell's protein recycling machinery
- effectively decreases the protein's half life or time that a protein will exist in a cell

Constructs

- Three main repressible/inducible promoter systems (rep-prom):

- ① LacI-plac
 - ② TetR-pTet
 - ③ AraC-pBAD
- } all found in the JGEM registry

- Two different ways to approach construct designs:

- ① Production fully expressed in "off" state then turned "on" to a specific level in addition to an inducer

- have a chromoprotein under control of an inducible protein, study/measure the steady state induced level of production, measure how protein is produced

- changing RBS from strong to weak > ↓ amount of protein, ↑ time required to produce

- adding degradation tags

- main thing to characterize → how different modifications (RBS) and protein deg tags change levels of induction

- ② Production is "on" initially then turned "off" upon addition of an inducer

- use 2 different inducible promoters to create a switch

- little more "tricky"

SIGNATURE

Continued to page

DATE

2/4/2014

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE Gene Expression

Continued from page

• Gene Expression in Prokaryotes (ex. lac operon)

• how genes are turned ON or OFF

• prokaryotes have simpler gene expression than eukaryotes → few controls

• Operons: sections of DNA that consist of ≥ 1 genes and their controlling elements- middle of an operon are structural genes → sections that code for ≥ 1 mRNA molecules [that code for proteins]

- promoter region → "turns on" the prokaryotic gene; in the beginning of the operon


• RNA pol must attach itself to the promoter to begin mRNA synthesis

- operator region → "turns off" the whole operon if a regulatory molecule (i.e. repressor) attaches at the operator

• RNA pol is blocked from continuing down the strand

* selective gene expression - blocking/unblocking the operator - is how bacterial cells produce only certain proteins at certain times

• Operon Model:



• Gene Expression in Eukaryotes

• no operons, but instead eukaryotes have enhancers and silencers

- silencer → sections of DNA that decrease the expression of a gene

- enhancer → sections of DNA that increase a gene's expression

• Eukaryotic genes are always "turned on", but expression can be significantly increased or decreased with enhancers and silencers.

• Prokaryotic Gene Regulation of Metabolism

① enzyme inhibition → slows down enzymatic processes and provides effective control of the metabolic pathway

② regulatory protein → protein that binds to the operator [region of the operon] to control their metabolic pathways over time

ex. repressors

SIGNATURE

Continued to page

DATE

2/18/2016

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE Gel Electrophoresis Lab PROJECT

Continued from page

Purpose: The objective of completing gel electrophoresis is to analyze DNA samples by separating bands of DNA according to their size.

Background Information: The agarose gel contains wells that the DNA/RNA samples are placed in using a micropipette. The gel has pores that allow the DNA/RNA to move towards the opposite side. The gel is put in an electric box with [TAE] buffer. The wells are on the negative (black) side, while the other end is at the positive side (red). Because DNA is negative, it will "run to red" since it is attracted to the positive charge. Longer strands take more time to maneuver through the pores, so shorter bands will travel further down the gel.

<u>Materials</u> :	agarose gel	staining tray
	methylene blue staining solution (10 mL)	thermometer
	TAE buffer	biological reagents/samples
	dH ₂ O	metric ruler
	2 beakers (100 mL)	resealable bag
	digital micropipets	white surface/paper
	electrophoresis chamber/power supply	lightbox
	marker	disposable pipet tips
	paper towels	

Protocol:

- 1) Gently slide the gel into the chamber with the wells on the negative side.
- 2) Pour the TAE buffer on top of the gel until it is completely submerged.
- 3) Gently flick the tubes containing the DNA samples to mix its contents.
- 4) Using a micropipette, gently transfer about 20-30 μ L of sample into a well. Be careful not to puncture the well.
- 5) Repeat step 4, and make sure a new, clean tip is used for each sample.
- 6) Place the lid on the chamber and turn it on by plugging the

Continued to page 10

SIGNATURE

DATE

2/26/2016

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

Continued from page 9

power source to an electrical source. Make sure the wires are properly matched.


7) Run it for about 30 minutes to 2 hours. The first bands should be about 1cm from the end of the gel. Turn the power off when done.

8) Carefully take out the gel and put it on a lightbox to analyze the bands.
* If necessary, add 40 ml of methylene blue staining solution onto the gel. Allow for staining to occur for about 5 minutes before putting it on the lightbox.

TITLE Genetic Circuits


PROJECT

Continued from page



• The Toggle Switch → circuit consisting of two repressors

- ① Tetr
 - ② LacI
- } both mutually repress each other



- The LacI repressor is inhibiting the P-Tet-2 promoter → repressor reporter
- The Tetr repressor is inhibiting the P₂ TetO1 promoter → promoter reporter
- The ATC inducer inhibits the expression of Tetr repressor, therefore promoting the expression of P₂ TetO1 promoter and the LacI repressor.
- The IPTG1 inducer inhibits the expression of LacI repressor, therefore promoting the expression of P-Tet-2 promoter and the gene of interest.

Continued to page

SIGNATURE

DATE

3/1/2016

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE Gatech Notes on Assemblies PROJECT



Continued from page

Types of antibiotics to use in plasmids: [4]

- ① Kanamycin (Kan)
- ② Chloramphenicol (Cam)
- ③ Tetracycline (Tet)
- ④ Ampicillin/Carbenicillin (Amp/Carb)

Types of Plasmids: [3]

- ① IC3 → Cam resistant
- ② IAK3 → Amp and Kan resistant
- ③ IAT3 → Amp and Tet resistant



Golden Rule of
3A Assemblies:

- ① ES
- ② XP
- ③ EP

ATURE

Continued to page

DATE

3/4/2016

USED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

7/11/2016 GIT WPLAC SATS PROJECT

Continued from page

1) Mini-Prep

- 2500 xg for 5 min instead of 10,000 xg for 1 min
- Once solution II is put in, do not let it sit for > 5 min
- After soln III microcentrifuge for $\geq 13,000$ mg or 17900 rcf for 10 min
- Centrifuge wash buffers for 1 min.
- Step 14 - ~~centrifuge~~ w/ elution buffer for 1 min (not 60 s)
- not necessary to repeat last wash buffer
- * Step 12 \rightarrow centrifuge empty mini column @ max speed for 2 min. to remove ethanol
- elution buffer - 50 μ L

* Vertex before nanodrop

2) Nanodrop

73.4 ng/ μ L

A260

1.469

A280

0.833

260/280

1.76

(1.8)

260/230

1.91

(2.1)

73.4 ng

μ L

1 μ g

1000 ng

=

0.0734 μ g/ μ L

13.624 μ L = 1 μ g DNA

13.624 μ L

avg DNA

2 μ L restriction enzyme (E)

1 μ L standard enzyme (S)

+ 3.38 dH₂O

20 μ L total volume

3) Digest

- ① water

② DNA

③ enzyme

order

* 200 band did not show up

* re-mix before non-re mix

- 1 hr @ 37°C, 20 min at 80°C

4) Ligation

3:1 molar ratios

vector: insert

	vector	IC3	insert length	vector length	total
			2000	2200	3200
digest 1 μ g	IAT3 PART 1 (BW P1AC)		200	3450	3650
	IAT3 PART 2 (C1PXP)		2000	3200	5200

~ [DNA]
ng/ μ L \cdot μ L

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

Continued to page

1/total length (kb)

0.31

0.27

0.19

TITLE 1/11/2016 GIT

PROJECT

Continued from page

vector

$$1/\text{total } \lambda = 0.31 \rightarrow \frac{.31}{.21} \times 3$$

μL

1

Part 1

$$= 0.27$$

3.4

Part 2

$$= 0.19 \rightarrow \frac{.31}{.19} \times 3$$

4.9

* actual measurements

+ 2 μL ligase buf

+ 1 μL ligase

+ 7.7 μL H_2O

20 μL

EXPERIENCE PRACTICE

length insert

length vector

total

1/total λ

μL

vector

375

1000

3200

4200

0.24

1

Part 1

PKR lac 5 GFP DAs 1800

3400

5200

0.19

3.8

3.789

Part 2

BR11 ClpXP-cl 375 3000

2900

5900

0.17

4.2

4.235

+ 2 ligase buf

+ 1 ligase

\Rightarrow 8 μL H_2O 7.976

20 μL

= phosphatase - 1 hr, 20 min deactivation

= last = enzyme

- 1 hr

- no deactivation required

- buffer \rightarrow has ATP, aliquate to smaller amts

5) TRANSFORMATION

- thaw on ice for 5 min

- 10 μL ligation mixture into 100 μL competent cells

- flick to mix

- ice for 30 min

- 200 μL of LB

- incubator @ 37°C for one hr

6) Plate

150 μL of cells on plate

SIGNATURE

Continued to page

DISCLOSED TO AND UNDERSTOOD BY

DATE

DATE

PROPRIETARY INFORMATION

TITLE 7/12/2016 GT

PROJECT

Continued from page

① Digest of weak PLAC (and chromoprotein)

$$\frac{73.4 \text{ ng}}{1 \text{ mL}} \times \frac{1 \text{ mg}}{1000 \text{ ng}} = 0.0734 \text{ mg/mL} \quad 13.62 \text{ mL} = 1 \text{ mg DNA}$$

13.62 mL DNA
2 mL RE (E)
1 mL SE (S)
+ 3.38 mL dH₂O
20 mL total

reminder: ① water ② DNA ③ enzyme
*never use the last bit of DNA

* Well 1 (edge) - LH digest weak PLAC
* Well 2 (next to ladder) - digest XP blue chromoprotein

200bp band still did not show up

② Ligation of pAR and the chromoprotein

		length	length	total	~[DNA]	mL
vector	AK3	1000	3200 3200	4200	0.229	1
Part 1	pAR lacI 3TG	1113	3200	4313	0.232	2.91
Part 2	chromoprotein IC3 2100	669	2200 2100	2869	0.341	1.93

1 mL vector
2.91 mL Part 1
1.93 mL Part 2
2 mL ligase buffer
1 mL ligase
+ 11.16 dH₂O
20 mL total

ACTUAL MEASUREMENTS

		length	length	total	~[DNA]	mL
vector	AK3	1000	3200	4200	0.2381	1
Part 1	pAR lacI 3TG	1113	3200	4313	0.2319	3.13
Part 2	blue chromo IC3	669	2100	2769	0.3611	1.93

2 lig buf
1 ligase

10.87 mL dH₂O

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE WORKFLOW WEEK 1 (8/15-8/19) PROJECT

Continued from page

① Liquid culture (overnight) → 8/16/2016

1a and 1b: pAR LACI TS PURP DAS DH10 8/16/2016

2a and 2b: pAR LACI TS PURP AK3 DH10 8/16/2016

3a and 3b: pAR LACI TS PURP LAA 8/16/2016

② Miniprep (with Qiagen) and Nanodrop → 8/17/2016

1a = A260: 1.691

2a = A260: 2.257

3a = A260: 1.952

A260/A280 = 1.81 1.8

A260/A280: 1.78 1.8

A260/A280: 1.72 1.8

82.5 ng/μL

112.8 ng/μL

97.6 ng/μL

1b = A260: 1.434

2b = A260: 2.226

3b = A260: 2.285

A260/A280: 1.74 1.8

A260/A280: 1.79 1.8

A260/A280: 1.83 1.8

71.7 ng/μL

111.3 ng/μL

114.3 ng/μL

③ Digest (① water, ② enzymes, ③ DNA) → 8/17/16

1a = 12.20 μL DNA

2a = 8.87 μL DNA

3a = 10.25 μL DNA

2 μL RE

2 μL RE

2 μL RE

1 μL SE

1 μL SE

1 μL SE

4.80 μL dH₂O4.13 μL dH₂O6.75 μL dH₂O

1b = 13.95 μL DNA

2b = 8.98 μL DNA

3b = 8.75 μL DNA

2 μL RE

2 μL RE

2 μL RE

1 μL SE

1 μL SE

1 μL SE

3.05 μL dH₂O8.02 μL dH₂O8.25 μL dH₂O

RESULTS from gel (8/18/16) → bands matched up at correct band lengths

④ LIGATION → 8/18/16 (only LAA & DAS)

	Insert μL	vector μL	Total μL	~[DNA]	μL
Part 1: pAR RBS LACI TS purple AK3	690	3200	3890	0.257	3.8
Part 2: RII ClXP CI AK3	2974	3200	6174	0.162	6
vector: IC3	1000	2100	3100	0.323	1

2 ligase buffer

1 ligase

+ 6.25 μL dH₂O

20 μL rxn

SIGNATURE

DATE

Continued to page

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

PROJECT

TITLE

Continued from page

AK3

	Insert μ	Vector μ	Total	~[DNA]	μ L
Part 1: pAR RBS RBS LACI TS PURP DAS	921	3200	4121	0.243	4
Part 2: RII ClpXP CI AK3	2974	3200	6174	0.162	4
Vector: IC3	1000	2100	3100	0.323	1

2 ligase buff
1 ligase
+ 6 ~~20~~ μ L dH₂O
20 μ L rxn

	Insert μ	Vector μ	Total	~[DNA]	μ L
Part 1: pAR RBS LACI TS PURP LAA AK3	920	3200	4120	0.243	4
Part 2: RII ClpXP CI AK3	2974	3200	6174	0.162	4
Vector: IC3	1000	2100	3100	0.323	1

2 ligase buff
1 ligase
+ 6 ~~20~~ μ L dH₂O
20 μ L rxn

- ⑤ Transformation \rightarrow LAA & DAS \rightarrow 8/18/2014 (includes plates)
- 1) 8/18/2014 pAR LACI TS PURP DAS AK3 DH10 iGEM
 \rightarrow 8/19/14 \rightarrow Results: ~~lots of small, purple colonies~~ ~~dull purple colonies, <100~~ \rightarrow a few white colonies
 - 2) 8/18/14 pAR LACI TS PURP LAA AK3 DH10 iGEM
 \rightarrow 8/19/14 \rightarrow Results: ~~<100 colonies, purple & more dull than~~ many small, purple colonies
 - 3) 8/18/14 R ClpXP CI in AK3
 \rightarrow digested (2.2 μ L DNA, 2 μ L X, 1 μ L P, 14.8 dH₂O)
 \rightarrow 8/19/14 \rightarrow Results: lots of colonies, various sizes, some are more clustered than others
 - 4) 8/18/14 pUC19 control
 \rightarrow 8/19/14 \rightarrow Results: no growth, part of agar gone?

TITLE 10/13/2016 Inoculation

PROJECT

Continued from page

- ① PAR LACI GFP LAA ClpXPCI
 - ② PAR LACI GFP Ø
 - ③ PAR LACI GFP DAS
 - ④ PAR LACI GFP LAA
- } CONSTRUCTS

cells { DH10
K12 Wild
K12 ClpP

inducing IPTG levels { NO IPTG
10 µM IPTG
100 µM IPTG
1 mM IPTG

not enough cells w/ constructs

10/15/2016 Inoculation
DH10 cells

K12 Wild

K12 ClpP

PAR LACI
GFP Ø

x 3

x 3

x 3

PAR LACI
GFP DAS

x 3

x 3

x 3

PAR LACI
GFP LAA
ClpXPCI

x 3

x 3

x 3

in triplicates,
induced
with
0 µM IPTG
and 100
µM IPTG

CONTROLS

x 1

x 1

x 1

x 1 Plain LB

RESULTS: cells grew, but did not fluoresce

SIGNATURE

DATE

Continued to page

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE 10/11/2016 Inoculations

PROJECT

Continued from page

CELL TYPE	CONSTRUCT	# TUBES	IPTG
DH10	PAR-LACI TS Purple CIPXPCI	3	0 μ M
	PAR-LACI TS Purple DAS CIPXPCI	3	0 μ M
	PAR-LACI TS Purple LAA CIPXPCI	3	0 μ M
	control	1	0 μ M
K-12 wild	PAR-LACI TS Purple CIPXPCI	12	0 μ M, 10 μ M, 100 μ M, 1 mM ↓
	PAR-LACI TS Purple DAS CIPXPCI	12	
	PAR-LACI TS Purple LAA CIPXPCI	12	
	control	1	
K-12 CIP knockout	PAR-LACI TS Purple CIPXPCI	12	0 μ M, 10 μ M, 100 μ M, 1 mM ↓
	PAR-LACI TS Purple DAS CIPXPCI	12	
	PAR-LACI TS Purple LAA CIPXPCI	12	
	control	1	
no cells	plain LB	+ 1 38 85 tubes	0 μ M

4 mL ~~plain~~ ^{KAN} LB into each tube
 140 mL → KAN LB with NO IPTG
 100 mL → KAN LB with 10 μ M IPTG
 100 mL → KAN LB with 100 μ M IPTG
 100 mL → KAN LB with 1 mM IPTG

3 μ L DNA into DNA
 ↳ take smear of colony
 from plate, dilute into
 50 μ L water

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION