

iGEM 2016: Team Pittsburgh

Week 17 Lab Notebook

Monday, September 12

Mini-prep codon optimized BFP, GFP, RFP, YFP, and OFP (Maya)

Enzyme cleaned linearized PT3 and linearized PT7-RBS (Maya)

T3 Amplification from RBS-T3 (Maya)

50 uL reaction

10 uL 5X Phusion HF Buffer

1 uL 10mM dNTP

2.5uL 10 uM PstI-T3 Reverse 2

2.5uL 10 uM XbaI-T3 Forward 2

0.4 uL DNA (26.3 ng/uL) RBS-T3

0.5 uL Phusion

28.6 uL Water

50 uL reaction

10 uL 5X Phusion HF Buffer

1 uL 10mM dNTP

2.5uL 10 uM Xba-T3 Reverse 1

2.5uL 10 uM EcoRI-T3 Forward 1

0.4 uL DNA (26.3 ng/uL) RBS-T3

0.5 uL Phusion

28.6 uL Water

PCR: 98C for 30sec

[98C for 10 sec

52C for 30 sec

72C for 2:30] 30X

72C 10 min

4C hold

Transform G switch

2 uL 42.1 ng/uL switch into 50 uL cells

500 LB (SOC contaminated)

Cell-free DNAzyme concentration (Claire)

9/12

Monday, Sep. 12.

Titration of DNAzyme amt (lead) cell-free

no trigger

0.52 nM

1 nM

3.74 nM ~~5 nM~~

10 nM

~~25 nM~~

50 nM

100 nM

1 μ M

DNAzyme and DNA

		no lead	cleavage	cleavage	DNA
I 11-12	no trigger	0.7 μ L H ₂ O	0.26 μ L 10 nM hp		
I 13-18	0.52 nM	0.7 μ L 3.74 nM	0.28 μ L 66.9 μ M hp	0.7 μ L 3.74 nM	
			0.417 μ L 24 μ M Pb		
			0.033 μ L Pb buff.		
	1 nM	0.5 μ L 10 nM hp	0.25 μ L 200 nM hp	0.5 μ L 10 nM	
		0.2 μ L H ₂ O	0.417 μ L 24 μ M Pb	0.2 μ L H ₂ O	
			0.033 μ L Pb buff.		
J 19-24	3.74 nM	0.28 μ L 66.9 μ M hp	0.28 μ L 66.9 μ M hp	0.374 μ L 50 nM	
		0.42 μ L H ₂ O	0.417 μ L 24 μ M Pb	0.326 μ L H ₂ O	
			0.033 μ L Pb buffer		
J 11-16	10 nM	0.67 μ L 74.8 μ M hp	0.25 μ L 200 nM hp	0.5 μ L 100 nM	
		0.03 μ L H ₂ O	0.417 μ L 24 μ M Pb	0.2 μ L H ₂ O	
			0.033 μ L Pb buffer		
J 17-22	50 nM	0.25 μ L 1 μ M hp	0.25 μ L 1 μ M hp	0.25 μ L 1 μ M	
		0.45 μ L H ₂ O	0.417 μ L 24 μ M Pb	0.45 μ L H ₂ O	
			0.033 μ L Pb buffer		
J 23-24	100 nM	0.25 μ L 2 μ M hp	0.25 μ L 2 μ M hp	0.5 μ L 1 μ M	
K 11-14		0.45 μ L H ₂ O	0.417 μ L 24 μ M Pb	0.2 μ L H ₂ O	
			0.033 μ L Pb buffer		

9/12 99

K 15-20 $1 \mu M$ $0.25 \mu L$ $20 \mu M$ hp $0.25 \mu L$ $20 \mu M$ hp $0.5 \mu L$ $10 \mu M$
 $0.45 \mu L$ H_2O $0.417 \mu L$ $24 \mu M$ Pb $0.2 \mu L$ H_2O
 $0.033 \mu L$ Pb buffer

K 21-22 $5 \mu L$ leftover MM $5 \mu L$ H_2O

$20 \mu M$ hp: 1 part $74.8 \mu M$: 2.74 parts water
 5 : 13.7

$10 \mu M$ hp from \uparrow (1:1) \rightarrow 2.5 $74.8 \mu M$: $13.7 H_2O$
 nm

$20 \mu M$ hp: 5 parts $74.8 \mu M$: $5.48 H_2O$

$2 \mu M$ hp from \uparrow (1:9)

$1 \mu M$ hp from \uparrow (1:1)

38 reactions (make 40):

Sol'n A: 2 80

B: 1.5 60

RNAse : 0.25 $\times 40 =$ 10

Sub : 0.25 10

~~Sub~~ 0.3 12 @ $72.1 \mu g/\mu L$

[NEED TO RE-RUN]

2nd time: $10 \mu L$ @ 72.1
 $2 \mu L$ @ 61.9

Heat hairpin: $5 \mu L$ of $74.8 \mu M$.

95° 3 min

90° 2 min

85° 2 min

80° 2 min

70° 4 min

60° 4 min


40° 8 min

20° 8 min

4° ∞

re-run - every 5 deg. 2 min. for

Condition	hairpin				lead				DNA			
	1	2	average	st dev	1	2	average	st dev	1	2	average	st dev
no trigger	0.0936	0.1289	0.11125	0.024961								
0.52 nM	0.0699	0.0739	0.0719	0.002828	0.0689	0.0964	0.08265	0.019445	0.085	0.1003	0.09265	0.010819
3.74 nM	0.0512	0.065	0.0581	0.009758	0.0646	0.076	0.0703	0.008061	0.0565	0.0754	0.06595	0.013364
10 nM	0.0628	0.0668	0.0648	0.002828	0.0385	0.0381	0.0383	0.000283	0.1581	0.0969	0.1275	0.043275
50 nM	0.0687	0.0733	0.071	0.003253	0.1298	0.1662	0.148	0.025739	0.2234	0.2052	0.2143	0.012869
100 nM	0.3162	0.356	0.3361	0.028143	0.1149	0.293	0.20395	0.125936	0.9697	1.3128	1.14125	0.242608
1 uM	0.4054	0.5161	0.46075	0.078277	0.3923	0.5609	0.4766	0.119218	1.0448	1.316	1.1804	0.191767



Tuesday, September 13

Cell-free concentrations (Praneeth)

No trigger

0.52 nM hairpin

0.52 nM hairpin + 2 uM lead + water


0.52 nM DNA

0.52 nM hairpin

0.52 nM hairpin + 2 uM lead + water

0.52 nM DNA

Condition	hairpin				lead				DNA			
	1	2	average	st dev	1	2	average	st dev	1	2	average	st dev
no trigger	0.3593	0.3474	0.35335	0.008415								
0.52 nM	0.6082	0.2351	0.42165	0.263822	0.5069	0.3827	0.4448	0.087823	0.7023	0.8368	0.76955	0.095106
3.74 nM	0.5345	0.5009	0.5177	0.023759	0.9972	0.7883	0.89275	0.147715	0.3262	0.3056	0.3159	0.014566



Wednesday, September 14

Cell-free lead hairpin reaction

Solution A	2 uL	24 uL
Solution B	1.5 uL	18 uL
RNAse inhibitor	0.25 uL x12 reactions =	3 uL
Substrate	0.25 uL	3 uL
Switch	0.40 uL at 61.9 ng/uL = 25 ng	4.84 uL
Trigger	0.60 uL	

Triggers:


Plate	Condition	
C19-20	No trigger	0.60 uL water
C21-22	3.47 uM Pb Hairpin	0.25 uL 74.8 uM renatured hairpin 0.35 uL water
C23-24	3.47 uM Pb Hairpin + 2 uM lead + buffer	0.25 uL 74.8 uM renatured hairpin 0.25 uL Pb at 40 uM = 2 uM Pb 0.10 uL Pb buffer
D19-20	3.47 uM Pb Hairpin + 2 uM lead + buffer	0.25 uL 74.8 uM renatured hairpin 0.25 uL Pb at 40 uM = 2 uM Pb 0.10 uL water

D21-22	3.47 nM DNA trigger	0.374 uL 50 nM DNA trigger 0.226 uL water
--------	---------------------	--

PROCEDURE:

1. Make master mix
2. Add 8.8 uL master mix to 8 tubes
3. Add 2x what is listed in "Triggers" to each tube
4. Separate each tube into 2 tubes of 5 uL each
5. Incubate 2 hours at 37 degrees

Condition	1	2	average	st dev
no trigger	0.0868	0.0717	0.07925	0.010677
hairpin	0.3563	0.5855	0.4709	0.162069
lead with buffer	0.8271	1.134	0.98055	0.217011
lead with water	0.4966	1.4514	0.974	0.675146
DNA	0.8189	0.8812	0.85005	0.044053



Repeat with correct concentrations of trigger:


Triggers:

Plate	Condition	
D23-24	3.47 nM Pb Hairpin	0.25 uL 74.8 nM renatured hairpin 0.35 uL water
E19-20	3.47 nM Pb Hairpin + 2 uM	0.25 uL 74.8 nM renatured hairpin

	lead + buffer	0.25 uL Pb at 40 uM = 2 uM Pb 0.10 uL Pb buffer
E21-22	3.47 nM Pb Hairpin + 2 uM lead + buffer	0.25 uL 74.8 nM renatured hairpin 0.25 uL Pb at 40 uM = 2 uM Pb 0.10 uL water
E23-24	3.47 nM DNA trigger	0.374 uL 50 nM DNA trigger 0.226 uL water

Condition	1	2	average	st dev
no trigger	0.0868	0.0717	0.07925	0.010677
hairpin	0.1095	0.0881	0.0988	0.015132
lead with buffer	0.2032	0.1429	0.17305	0.042639
lead with water	0.5805	0.3999	0.4902	0.127703
DNA	0.6848	0.6271	0.65595	0.0408

*no trigger data from previous run



C

D

E

24 23 22 21 20 19
Agarose gel (Claire)
30 mL TAE buffer
0.3 g agarose
3 uL EtBr
Run at 100 V, 30 min
Lane 1: 1 kb ladder
Lane 2: lacZ 1
Lane 3: lacZ 2
Lane 4: T3 amp 1
Lane 5: T3 amp 2

Clear bands for T3 amp, but don't see anything for lacZ

Friday, September 16

Cell-free lead hairpin concentration variation

Friday, September 16.

Redo concentration curve:

3.74 nM, 50 nM, 100 nM 20 rxns

$$25 \text{ ng switch} @ \overset{50}{57.4} \text{ ng} / \overset{0.5}{\mu\text{L}} = 0.430 \mu\text{L} / \text{rxn}$$

$$= 0.57 \mu\text{L trigger} / \text{rxn}$$

O 13-14 no trigger: 0.50 μL H₂O.O 15-16 3.74 nM hp: 0.25 μL hp @ 74.8 nM
0.25 μL H₂OO 12-18 3.74 nM hp+: 0.25 μL hp @ 74.8 nM
2 μM Pb 0.25 μL Pb²⁺ @ 40 μM
0.07 μL H₂OO 19-20 3.74 nM DNA: 0.5 μL DNA @ 34.7 nMO 21-22 50 nM hp: 0.25 μL hp @ 1 μM
0.25 μL H₂OO 23-24 50 nM hp+: 0.25 μL hp @ 1 μM
2 μM Pb 0.25 μL Pb²⁺ @ 40 μM
8.7 μL P 17-18 50 nM DNA: 0.25 μL DNA @ 1 μM
0.25 μL H₂O

A short on volume

P 19-20 100 nM hp: 0.25 μL 2 μM hp
0.25 μL H₂OP 21-22 100 nM hp+: 0.25 μL 2 μM hp
2 μM Pb 0.25 μL 40 μM PbP 23-24 100 nM DNA: 0.5 μL 1 μM DNA.
$$0.25 \mu\text{L} \rightarrow 0.3 \mu\text{L}$$

$$0.5 \mu\text{L} \rightarrow 0.6 \mu\text{L}$$

9/16

103

Cell extract MM:

Sol'n A: 2 μ L	2.4 μ L	52.8 μ L	14.4
B: 1.5 μ L	1.8 μ L	39.6	10.8
RNAse: 0.25 μ L	0.3 μ L	6.6	1.5
Sub: 0.25 μ L	0.3 μ L	6.6	1.5
switch: 0.43 μ L	0.52 μ L	11.44	2.58
H ₂ O: 0.07 μ L	0.08 μ L	1.76	0.42
trigger	1 μ L		
	5 μ L	6 μ L	10.8 μ L/condition

Trigger MM:

no trigger: 0.6 μ L $\times 2 =$ 1.2 μ L H₂O

3.74 nM

hp: 0.3 μ L hp@74.8 \rightarrow 1 μ L hp@74.8 \Rightarrow 1.2 μ L
 0.3 μ L H₂O \rightarrow 1 μ L H₂O

hp+Pb: 0.3 μ L hp@74.8 \rightarrow 1 μ L hp@74.8 \Rightarrow 1.2 μ L
 0.3 μ L Pb@40 \rightarrow 1 μ L Pb@40

DNA:

3.74
 1.2 μ L @ 3.74 nM
 1 μ L 50 = 3 μ L
 0.34 μ L H₂O = 1.02

50 nM

hp: 0.3 μ L hp@1 μ M \rightarrow 1 μ L hp@1 μ M \Rightarrow 1.2 μ L
 0.3 μ L hp@H₂O \rightarrow 1 μ L H₂O

hp+Pb: \rightarrow 1 μ L hp@1 μ M \Rightarrow 1.2 μ L
 1 μ L Pb@40 μ M

DNA:

\rightarrow 1 μ L DNA@1 μ M \Rightarrow 1.2 μ L
 1 μ L H₂O

diff
 MM
 pipetted sep
 100 nM
 hp:
 hp+Pb
 DNA:


\rightarrow 1 μ L hp@2 μ M \Rightarrow 1.2 μ L 1 μ M hp
 1 μ L H₂O

\rightarrow 1 μ L hp@3 μ M \Rightarrow 1.2 μ L
 1 μ L Pb@40 μ M

1.2 μ L @ 1 μ M



Condition	hairpin			lead			DNA					
	1	2	average	st dev	1	2	average	st dev	1	2	average	st dev
no trigger	0.0557	0.0617	0.0587	0.004243								
3.74 nM	0.058	0.0454	0.0517	0.00891	0.1972	0.2863	0.24175	0.063003	0.4999	0.394	0.44695	0.074883
50 nM	0.6621	0.5589	0.6105	0.072973	0.2301	0.138	0.18405	0.065125	0.3075	0.3693	0.3384	0.043699
100 nM	0.3963	0.4312	0.41375	0.024678	0.2831	0.524	0.40355	0.170342	0.3158	0.6438	0.4798	0.231931



Saturday, September 17

Repeat above for 50 nM and 100 nM

9/17

Saturday, September 17.

Redo concentrations @ 50 nM & 100 nM

Cell extract mM

Sol'n A	2.4 μ L	$\times 4$	36 μ L	9.6
B	1.8 μ L		27 μ L	7.2
RNAse	0.3 μ L	$\times 15$ rxns =	4.5 μ L	1.2
Sub	0.3 μ L		4.5 μ L	1.2
Switch	0.45	@ 55.6 μ L =	6.75 μ L	1.8
1 μL H_2O	0.05		0.75 μ L	0.6
mg	0.6		2.25 μ L	0.6
	1.2 μ L			
	6 μ L		10.8 / rxn	10.24 rxn

remainder + 1.8 μ L
 H_2O no trig

Triggers

F21-22

50 nM hp

1.2 μ L 500 nM hp

F23-24

50 nM hp + Pb

1 μ L hp @ 1 μ M \Rightarrow 1.2 μ L
1 μ L Pb @ 40 μ M

G22-23

50 nM DNA

1.2 μ L 500 nM DNA

G24, H22

100 nM hp

1.2 μ L 1 μ M hp

H23-24

100 nM hp + Pb

1 μ L hp @ 2 μ M \Rightarrow 1.2 μ L
1 μ L Pb @ 40 μ M

diff

K23 & 24

100 nM DNA

1.2 μ L 1 μ M DNA2.5 μ L
RNAse +
0.35 μ L H_2O
added alone

no trig A24 F19-20

Condition	hairpin	1	2 average	st dev	lead	1	2 average	st dev	DNA	1	2 average	st dev
no trigger		0.2048	0.1547	0.17975	0.035426							
50 nM		0.8517	0.3305	0.5911	0.368544	0.6468	0.4565	0.55165	0.134562	2.8397	2.6602	2.74995
100 nM		0.4733	0.2644	0.36885	0.147715	0.3505	0.2162	0.28335	0.094964	0.7288	0.8967	0.81275

