

Lab Journal

Researcher Initials Key

AC = Alan Chen
AG = Amy Guo
AP = Amy Pai
JL = Jared Lamar
OE = Omar Ebrahim

June 27th

AC/OE

Inoculation of strains obtained from Dr. Christopher Hayes at UC Santa Barbara

- 10mL LB; 37C for 16h shaking
 - EPI100 w/ EC869 cosmid
 - 10uL ampicillin
 - X90 w/ EC93 cosmid
 - 10uL ampicillin
 - *Enterobacter aerogenes* (ATCC 13048)
-

June 28th

AC/OE

Mini-prep of EC93 and EC869 containing cosmids (x2 each)

-ZR Plasmid Miniprep Kit from Zymo

- EC869
 - 145.42 ng/μL
 - 198.85 ng/μL
- EC93
 - 90.81 ng/μL
 - 45.95 ng/μL

Isolation of *Enterobacter aerogenes* gDNA

-Using QIAGEN Progene Yeast/Bacteria KitB

- 34.87 ng/μL
 - 1.85 260/280
 - Made glycerol stocks (500uL culture, 500uL 50% glycerol)
 - Gibson primers designed for CDI assemblies in pSC101 AND pSB1C3
-

June 29th

- 12 new primers ordered from IDT
 - 4 for ATCC 13048 in pSB1C3
 - retain endogenous promoter instead of T7 promoter since membrane proteins are detrimental when overexpressed
 - 8 for Aerogenes pSC101

- pSC101 backbone replacement for ColE1 (too high copy)

AC

- Amplified Gibson EBL2 fragments 2 (Second Half of CDI-A) and 3 (CDI accessory factors)
 - Diluted stock primer solution 10 fold (100μM to 10μM)
 - Diluted gDNA to 1ng/μL
 - Original was 35.8ng/μL
 - Added 34μL water to 1μL of stock gDNA
 - Used Q5 2x protocol for 50μL reaction
 - 2.5μL forward and reverse primers
 - 1μL of gDNA as template
 - 25μL Q5 2x Master Mix
 - 19μL water
 - PCR Cycle Protocols
 - 98C----->30s
 - Following done for 28 cycles
 - 98C----->10s
 - 66C----->30s
 - 72C----->140s
 - 72C----->2min
 - 14C----->Infinite

June 30th

- Yesterday's PCR failed and showed no band other than primer dimers
 - Could be issue with gDNA concentration (we assumed 1ng would be enough)
 - Troubleshoot with increased gDNA concentrations
 - Also test with different polymerase, use KapaHiFi instead

Troubleshoot: PCR EBL2 Gibson frag 2 & 3 w/ varying gDNA template amount

4 PCR Rxns: (amplifying 3kb frag 3)

- Q5 2x protocol for 25μL reaction
 - 1.25μL forward and reverse primers
 - 1) 1 ng 2) 10 ng 3) 34.87 ng 4) 69.74 ng template gDNA
 - 12.5μL Q5 2x Master Mix
 - 9μL water
- PCR Cycle Protocols
 - 98C----->30s
 - 28 cycles
 - 98C----->10s
 - 65C----->30s
 - 72C----->88s
 - 72C----->2min
 - 14C----->Infinite


Expecting: 3 kb Fragment

Run separate rxn using KapaHiFi instead of Q5 in a 25μL rxn

- KapaHiFi protocol:
 - 12.5μL KapaHiFi Master mix
 - 0.75μL forward and reverse primers
 - ~10ng of gDNA
 - Fill with water until reach 25μL
 - PCR Cycle Protocols (Amplifying 7.8kb frag 2)
 - 98C----->30s
 - 28 cycles
 - 98C----->10s

- 65C----->30s
- 72C----->2min 40s
- 72C----->2min
- 14C----->Infinite

Expecting: 7.8kb fragment



July 1st

- ordered CPCR and Gibson primers for toxin switch of *Enterobacter aerogenes* CDI loci 1 and 2


AC

PCR Protocol

- 10ng of gDNA
- 12.5uL of Q5 Master Mix
- 1.25uL of forward and reverse primers
- Fill with water up to 25uL

Run 25 uL reaction

- Loci 2 Fragment 2: 65C (3 min 54s extension time on 28x cycle)
- Loci 1 Fragment 2: 66C (3min 30s extension time on 28x cycle) Fragment 3:64C(2min 30 extension time on 28x cycle)
- Everything else was the same as the standard Q5 protocol



Ladder/ EBL2 frag 3/ EBL1 frag 2/ EBL1 frag 3

AG

- gel purified all EBL2 frag 3 bands (x3), EBL2 frag 2, and EBL1 frag 3

Nanodrop

EBL2 frag 3

10ng: 66.99ng/uL ; 1.84

34ng: 65.85ng/uL ; 1.89

68ng: 59.41ng/uL ; 1.98

EBL1 frag 3

43.40ng/uL; 2.03

EBL2 frag 2

19.13; 1.86

July 5th

AC

EBL1 frag 2:

Redoing failed PCR, modified protocol by having elongation time be extended to be 30s/kb, annealing temperature set to 72C, ran 2 rxns (25uL and 50uL)

Using standard Q5 protocol for 25uL

Using standard Q5 protocol for 50uL

- 98C----->30s
- Repeat with 28 cycles
 - 98C----->10s
 - 66C----->30s
 - 72C----->3min 30s
- 72C----->2min
- 14C----->Infinite

EBL2 frag 4:


Using standard Q5 protocol for 25uL

- 98C----->30s
- Repeat with 28 cycles
 - 98C----->10s
 - 64C----->30s
 - 72C----->2min
- 72C----->14C

EBL2 frag 1:

Using standard Q5 protocol for 25uL

- 98C----->30s
- Repeat with 28 cycles
 - 98C----->10s
 - 65C----->30s
 - 72C----->3min
- 72C----->14C




ladder/ EBL2 frag 1/ EBL2 frag 4/ EBL1 frag 2 (25uL rxn)/ EBL1 frag 2 (50uL rxn)

July 6th

AC/OE

- PCR fragments that didn't work
 - EBL2
 - Fragment 4 (pSB1C3), fragment 1 (gDNA)
 - EBL1
 - Fragment 2
- Changes to PCR
 - Performed two 25uL reactions (doubled and quadrupled original template DNA amount)
 - 20ng gDNA/40ng gDNA
 - 2ng pSB1C3/4ng pSB1C3
 - Instead of original 28 cycles, did 32 cycles of PCR
 - Other times and temperatures of individual steps remained the same as yesterday



Ladder/ 2x EBL1 frag 2/ 4x EBL1 frag 2/ 2x EBL2 frag 1/ 4x EBL2 frag 1/ 2x EBL2 frag 4/ 4x EBL2 frag 4

- EBL1 fragment 2 failed again -----> most likely faulty primers
- EBL2 frag 2 and 4 worked by increasing template DNA
- Eluted all fragments into 10uL of elution buffer
 - 2x EBL2 frag 4
 - 17.51 ng/uL
 - 4x EBL2 frag 4
 - 21.49 ng/uL
 - 22.83 ng/uL

- 2x EBL2 frag 1
 - 24.94 ng/uL
 - 20.63 ng/uL
- 4x EBL2 frag 1
 - 14.51 ng/uL

July 7th

AG/OE


1. transform Kan^R into DH5-alpha (target) and transform EC93 cosmid into DH5-alpha (inhibitor)

July 8th

AG

PCR of EBL1 frag 2

- remade 10uM primer dilutions using same stock of 2F and 2R
- tried 10 and 20 ng amt of template DNA using both old and new aliquots
- 94 C for 30s
- 32 CYCLES:
 - 94 C for 10 s
 - 66C for 30 s
 - 72 C for 3:30 min
- 70 C for 20 s (final ext)
- 14C infinite hold



ladder/ old 10ng/ old 20ng/ new 10ng/ new 20ng

- observed growth on both 1:1 and 1:100 plates for DH5-alpha transformed w pSK33/EC93 cosmid
 - make 10uL starter cultures

July 9th

AG

- miniprep 5mL of each 10 mL starter culture
 - pSK33 insufficient yield
- test digest of EC93 with Hayes' strain as a control

X90 w/ EC93

- 500 ng DNA (14uL)
- 5uL NEB 2.1
- 1 uL EcoRI
- H2O to 50uL

DH5-alpha EC93

- 500 ng DNA (10.1uL)
- 5uL NEB 2.1
- 1 uL EcoRI
- H2O to 50uL

July 12th

- ordered: EBL1, EBL2 toxin switch CPCR and Gibson primers





AG/OE







- inoculate for autoaggregation test tomorrow morning
 - *Enterobacter aerogenes*
 - DH5a w/ EC93 cosmid
 - DH5a w/ pSK33 (neg control)
 - X90 w/ EC93
 - EPI100 w/ EC869

July 13th


AG/OE

Aggregation test results:

Strain/Cosmid	t = 0h	t = 2h
ATCC 13048		
DH5-alpha w/ EC93		

DH5-alpha w/ pSK33		
X90 w/ EC93 ** expected aggregation**		
EPI100 w/ EC869		

- X90 w/ EC93 did not show any aggregation → mini-prepped (291.22 ng/uL) and ran test digest (1hr at 37 C, 20 min at 80C w EcoRI)



ladder

test digest

- bands at 10kb and 8kb as expected (third band may be supercoiled cosmid)

cPCR of DH5-alpha w/ pSK33

- selected 8 colonies
- 1uL of 10uM F primer
1uL of 10uM R primer

1 uL colony solution (pipette tip to select one colony and resuspend in 100 uL water)

25 uL 2X APEX MM

water to 50 uL

95C --- 5 min

95C --- 30s

49C --- 30s (annealing)

72C --- 10s (ext)

72C --- 3min



ladder/ colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8

- nothing showed up on the gel
 - used incorrect CPCR primers
- reinoculated X90 w/ EC93 to redo auto-aggregation assay

OE

- Categorized CDI primers that had already been ordered into usable and unusable groups bc many had been designed improperly
- Worked on designing all new gibson primers for CDI constructs

AC


Prepared materials needed Lambda Red CdiA knockout

- Obtained pTKRed plasmid
 - Spectinomycin resistance
- Prepared spectinomycin agar plates
 - 12.5g LB
 - 7.5g Bacto-Agar
 - 500mL dH2O
 - 500uL Spectinomycin
 - Poured into miniplates
- Prepared 800mL of autoclaved Milli-Q water
 - 2 bottles of 400mL
 - Autoclaved in liquid cycle
- Prepared autoclaved 250mL flask
- Ordered primers needed to make KanR genes
 - 50bp homology + KanR + 50bp homology
 - Homology regions come at beginning of CdiA and end of CdiA
 - Doesn't include CT or Cdi-I regions
 - Ordered for EBL1 and EBL2

July 14th


AG

- no aggregation observed in EC93 culture → prepare for competition assay to test for proper CDI complex expression



X90 w/ EC93; t = 2h

- redoing test digest → made glycerol stock (LABEL: X90 EC93 7/14/16)



- replated DH5-alpha pSK33 on new kan plate for single colonies
- plated DH5a w/ Loci 2 pSB1C3 (1:1)
- plated DH5-alpha w/ EC869 (1:1, 1:100, 1:1000)

AC

- gibson assembly of pSB1C3 loci 2
 - Fragment concentrations
 - EBL2 frag 2: 19.3ng/uL
 - EBL2 frag 3: 65.85ng.uL
 - 2 fold dilution for 32.925ng/uL
 - EBL2 frag 4 (pSB1C3 Backbone): 22.83ng.uL
 - EBL2 frag 1: 24.94 ng/uL
 - Add all fragments in equimolar amounts and fill up to 5uL of DNA mix
 - EBL2 frag 1: 0.775uL
 - EBL2 frag 2: 1.3uL
 - EBL2 frag 3 (diluted): 0.147uL
 - EBL2 frag 4: 0.282uL
 - dH2O: 2.35uL
 - Add in Mastermix and incubate at 50C for 1hr

- 5uL DNA mix + 5uL GA HiFi 1-Step Mastermix (2x) on ice
 - Incubate at 50C for 1 hr
- PCR purification (elute into water)

OE

- transform DH5-alpha with pSB1C3 loci 2
- transform DH5-alpha with EC869 cosmid

July 15th

AG/OE

- ordered cPCR primers for pSK33 (the other primers used in the cPCR were for EC93 cosmid accidentally)
- no autoaggregation in Enterobacter observed

EBL1 Toxin Switch (TS)

Q5 protocol

1.25uL of 10uM F primer

1.25uL of 10uM R primer

10ng/20ng gDNA

12.5uL Q5 2X MM

water to 25uL

initial denaturation: 98 C --- 30s

30 cycles: 98 C --- 10s

Tm --- 30s

72 C --- 20-30s/kb

final ext: 72C --- 2 min

Hold: 10 C --- infinite

Fragment 1 (8kb)

- anneal at 64 C
- 72 C for 4 min

Fragment 2 (pSK33)

- digest with EcoR1 and BamH1
 - couldn't do bc not enough pSK33

Fragment 3 (5.8kb)

- anneal at 63 C
- 72 C for 2 min 54 s


Fragment 4 (5.6kb)

- anneal at 65 C
- 72 C for 2 min 48 s

Fragment 5 (900bp)

- anneal at 67 C
- 72 C for 28 s

- run all fragments on gel (except frag 2-- pSK33) and gel purify successful products



lanes: ladder/ frag 1 (10ng)/ 1 (20ng)/ 3 (10ng)/ 3 (20ng)/ 4 (10ng)/ 4 (20ng)/ 5 (1ng);/5 (5ng)

- frag 1, 3, and 4 are from Enterobacter gDNA; frag 5 from EC93 cosmid
- smearing in band 4 probably due to nonspecific primer binding to other regions of gDNA
- gel purified all fragments and combined bands of the same fragment
 - frag 1: 63.11 ng/uL
 - frag 3: 77.48 ng/uL
 - frag 4: 63.87 ng/uL
 - frag 5: 110.34 ng/uL
 - all 260:280 between 1.8 and 1.9
- inoculation for test digests, making electrocompetent cells and cloning more pSK33:
 - 5mL DH5a w/ EC869
 - 5mL pSB1C3 w/ Loci 2
 - 5mL Enterobacter
 - 10 mL DH5a w/ pSK33

AC

- Over-grew Enterobacter (went from 0.116 to 1 in 2.5h)
 - Initially added 0.5mL overnight culture to 50mL of LB in 37C shaker

July 16th

AG

- observed growth in DH5a w/ pSB1C3 EBL2
- observed growth in DH5a w/ EC869
 - made glycerol stocks
- test digest:
 - pSB1C3 EBL2 (8kb, 7.6kb, 4kb)
 - XmaI: 1uL (37 for 1h, 65 for 20min)
 - BsaI: 1uL (37 for 1h, 65 for 20min)
 - H2O to 50 uL (37.6uL)
 - 500 ng DNA = 5.44 uL sol'n (91.82ng/uL)
 - 5uL cutsmart buffer

EC869 cosmid (10kb, 8kb)

- EcoR1: 1uL
- H2O to 50 uL (42.6)
- 500 ng DNA = 1.45 uL sol'n (344 ng/uL)
- 5uL NEBuffer 2.1
- mini-prepped DH5a-pSK33 (66.45 ng/uL)
- test digest pSK33
 - 1uL EcoR1 (37 1h, 65 20min)
 - 1uL BamH1-HF (37 1h)
 - H2O to 50uL (33uL)
 - 10 uL DNA
 - 5uL NEBuffer 2.1
- PCR of flag tag fragments (4 and 5)

Q5 protocol

1.25uL of 10uM F primer

1.25uL of 10uM R primer

10ng gDNA or 1ng of EC93 cosmid

12.5uL Q5 2X MM

water to 25uL

initial denaturation: 98 C --- 30s

25 cycles: 98 C --- 10s

Tm --- 30s

72 C --- 20-30s/kb

final ext: 72C --- 2 min


Hold: 10 C --- infinite

Fragment 4 (5.6kb)

- anneal at 65 C
- 72 C for 2 min 48 s

Fragment 5 (900bp)

- anneal at 67 C
- 72 C for 27 s



4FT 5FT pSK33 EC869 pSB1C3 w Loci 2

- gel purified pSK33 (17.93 ng/uL bad 260:280); frag 4 FT (58.81ng/uL); frag 5 FT(106.3 ng/uL)
 - pSB1C3 wrong size
- inoculate more Enterobacter/DH5-a with pSK33 (2x 10mL)

AC/AG


- made electrocompetent Enterobacter cells
 - Inoculate overnight culture
 - Add 0.5mL of overnight culture to 50mL of LB in a 250mL flask

- Incubate in shaker at 37C
- Take OD600 after 1 hour
 - Continue taking at around 15min intervals until reach an OD of 0.4-0.6
- Remove from shaker and transfer into two 50mL plastic Falcon tubes (25mL of liquid culture in each)
- Immediately place tubes into ice-water slurry and mix tubes and ice for 4-5min
 - Move to cold room to perform other steps
 - Cells must be kept at 4C or below until transformation
- Centrifuge Falcon tubes at 3000rpm (~3000xg) for 3min at 4C
 - Use large centrifuge next to incubator
- In cold room, remove LB without disturbing pellet
- Add in 50mL of autoclaved Milli-Q Water per tube
 - Repeat above step 2 **MORE TIMES**
- Resuspend cells in 5mL of Milli-Q water
- Aliquot into 30uL volumes for transformation
- transformed with pTKred
 - no cells (arc time 6.0 and no pellet from centrifugation)
 - most likely cells were too dilute
 - Perhaps try final resuspension in 100uL of water instead of 5mL of water

July 17th

AG

- mini-prepped two 10mL cultures of DH5a w/ pSK33
 - 80.25 ng/uL and 66.34 ng/uL
- digest of pSK33 to linearize it for gibson
 - 1uL EcoR1 (37 1h, 65 20min)
 - 1uL BamH1-HF (37 1h)
 - H2O to 50uL (28uL; 30.5uL)
 - 1.5 ug DNA
 - 5uL NEBuffer 2.1
 - 2 hr at 37 C → immediately gel purify (BamH1 cannot be heat inactivated)
 - expect band at 4kb



ladder/ pSK33

- PCR PURIFY RIGHT AWAY (to get rid of BamH1 enzyme)
 - don't nanodrop
- gel purify fragments (15uL sample, 5uL dye, 10 uL H2O)
 - pSK33 = 43.55 ng/uL
- Gibson Assembly (no FT)
 - Frag 1: 0.40 uL (25 ng)
 - Frag 2: 0.26 uL (11.21 ng)

- Frag 3: 0.23 uL (17.8 ng)
- Frag 4: 0.27 uL (17.2 ng)
- Frag 5: 0.13 uL (13.8 ng)
- H2O: 3.82 uL
- Gibson Assembly (FT)
 - Frag 1: 0.40 uL (25 ng)
 - Frag 2: 0.26 uL (11.21 ng)
 - Frag 3: 0.23 uL (17.8 ng)
 - Frag 4: 0.27 uL (17.2 ng)
 - Frag 5: 0.13 uL (13.8 ng)
 - H2O: 3.82 uL
- 1. thaw GA MM on ice
- 2. put fragments and H2O in PCR tube
- 3. combine 5uL fragments and 5uL MM → pipette 3 times to mix
- 4. incubate at 50 C for 1 hr
- 5. PCR purify into water (10 uL)
- 6. use 2uL to transform cells
 - a. plated 1:1 and 1:100 of both pSK33-TSL1 with and without flag tag (on kan plates)

OE

- the psk33 w/ Aerogenes Loci 2 [EC93 Toxin] and Flag Tag Benchling files did not have the attached primers after the CDI-I region that was necessary for proper amplification. Luckily some of the old psc101 primers work.
 - 1F_EBL2-p33-TS_GP
 - 1R_EBL2-p1c3_GP
 - 2F_EBL2-p1c3_GP
 - 2R_EBL2-p33-TS_GP
 - 3F_EBL2-p33-TS_GP
 - 3R_EBL2-p101-TS_GP
 - 3R_EBL2-p33-TS_GP
 - 4F_EBL2-p101-TS_GP

1.25μL forward primer (to get 500nM)

1.25μL reverse primer

10ng of DNA insert

12.5 μL of Q5 enzyme (for 2X reaction)

H2O up to 25μL

PCR fragments:

Fragment 1: 1F_EBL2-p33-TS_GP/1R_EBL2-p1c3_GP

Anneal at 65°C and extend for 3 minutes and 22 seconds

Fragment 2: 2F_EBL2-p1c3_GP/2R_EBL2-p33-TS_GP

Anneal at 66°C and extend for 3 minutes

Fragment 3: 3F_EBL2-p33-TS_GP/3R_EBL2-p101-TS_GP

Anneal at 66°C and extend for 32 seconds

Fragment 4: 4F_EBL2-p101-TS_GP/3R_EBL2-p33-TS_GP

Anneal at 64°C and extend for 2 minutes and 42 seconds

initial denaturation: 98 C --- 30s

25 cycles: 98 C --- 10s

Tm --- 30s

72 C --- 20-30s/kb

final ext: 72C --- 2 min

Hold: 14 C --- infinite

July 18th

AG

- colonies grew on both 1:1 and 1:100 plates for both DH5a w/ EBL1 TS w/o flag tag and with flag tag
 - cPCR (6 colonies, 3 CPCR each)

1uL of 10uM F primer

1uL of 10uM R primer

1 uL colony solution (pipette tip to select one colony and resuspend in 100 uL water)

25 uL 2X APEX MM

water to 50 uL (22uL)

95C 5 min

95C 30s


Tm 30s (annealing)

72C 10s (ext)

72C 3min

12C hold

- Region 1 (199bp) red
 - 51 C annealing temp
- Region 2 (193bp) blue
 - 52 C
- Region 3 (140bp) green
 - 50 C annealing temp
- CPCR failed (no bands in the right place → try test digest of 4x 10mL cultures)




ladder (1kb)-- accident, C1R1, C1R2, C1R3, C2R1, C2R2, C2R3, C3R1, C3R2, C3R3, C4R1, C4R2, C4R3, C5R1, C5R2, C5R3, C6R1, C6R2, C6R3, ladder (50bp)

- inoculated colonies 1, 3, 4 and 5 (5 had the best cPCR result)

OE

- Performed PCR twice, both failed... Do not use bio-rad pcr tubes in ProFlex because the tube caps will come off
- Redoing PCR with original tubes



ladder/ EBL2 frag 1/ frag 2/ frag 3/ frag 4


- Gibson of EBL2
 - F1 PCR: 77.99 ng/μL 2.41 260/280
 - F2 PCR: 42.73 ng/μL 3.41 260/280
 - F3 PCR: 50.57 ng/μL 2.47 260/280
 - F4 PCR: 24.91 ng/μL 4.07 260/280
- Use 25ng for biggest frag (F1): 0.33 μL
- Use 22.2ng for 2nd largest frag (F2): 0.52 μL
- Use 20.05ng for 3rd largest frag (F4): 0.81 μL
- Use 14.9ng for 4th largest frag (pSK33): 0.35μL
- Use 18.64 ng for 5th largest frag (toxin/cdil): 0.37μL
 - 2.62 H₂O
- Run at 50°C for an hour
- Preparation of electrocompetent cells (E. Aero)
 - Inoculate culture
 - Add 0.5 mL of overnight culture to 50 mL of the LB in a 250 mL flask
 - Incubate @37 C in shaker
 - Take OD600 after 1 hour
 - Continue taking @ around 15 min intervals until you reach an OD of 0.4-0.6
 - Remove from shaker and transfer into two 50 mL plastic falcon tubes (25 mL of liquid culture in each)
 - Immediately place tubes into ice-water slurry and mix tubes and ice for 4-5 minutes
 - Move to cold room
 - Cells must now be kept @ 4 C or lower for transformation
 - Centrifuge falcon tubes at 3000 rpm (~3000g) for 3 minutes @ 4 C
 - In cold room, remove LB without disturbing pellet
 - Add in 50 mL of autoclaved Milli-Q water per tube
 - Repeat last two steps two more times
 - Resuspend cells in 100 μL of Milli-Q water
 - Aliquot into 30 μL volumes for transformations
- Transform lambda-red

July 19th

AG

- mini-prepped colonies 1, 3, 4, and 5

- test digest (expected frag: 11kb, 5.6kb, 4kb, 2kb, 1.5kb)
 - 1uL XmaI
 - 5uL cutsmart buffer
 - 1ng DNA
 - H2O up to 50uL
- digest for 1.5 hr and run the gel at 95V for 1hr



ladder, colony 1, 3, 4, 5


AC/OE

- Re tried making competent cells
 - Use overnight culture and take OD 600
 - Add 0.5mL culture to 50mL LB and incubate at 37C
 - 0.416 OD after 1hr and 50min
 - Transfer into 2 falcon tubes of ~25mL each
 - Swirl falcon tubes in ice-water slurry for 5min
 - Spin at 3000g for 5min (no pellet collected)
 - Spin again at 3000g for 5min (no visible pellet)
 - Removed supernatant
 - Replaced with 25mL of Milli-Q Water
 - Spin at 3000g for 5min
 - No pellet collected again so assumed failure

July 21st

AG

- cPCR region 1 on 32 colonies
 - H2O to 25 uL (11uL)
 - 0.5uL forward primer
 - 0.5uL reverse primer
 - 1uL colony suspension
 - 12.5uL APEX master mix
- master mix cPCR:
 - 460uL H2O
 - 20uL forward primer
 - 20uL reverse primer
 - 500uL APEX
 - put 24 uL of the MM into each PCR tube and 1 uL of the colony suspension



ladder/ colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8/ 9/ 10/11/ 12/ 13/ 14/ 15/ 16
ladder/ colony 17/ 18/ 19/ 20/ 21/ 22/ 23/ 24/ 25/ 26/ 27/ 28/ 29/ 30/ 31/ 32


OE/ AC

- Remade enterobacter aerogenes
- Changes:
 - Spun @ 8000 rpm for 10 minutes
 - Last pellet was very soft and first two were very hard
 - Resuspending in 100μL was difficult because pellet was disturbed too easily
 - When transforming, the cells looked like they weren't very concentrated
 - Arc time was 5.9 and we accidentally put the heat sensitive lambda red transformations in the shaker @ 37 C for 10 minutes as opposed to the 30 C
 - We also ran out of normal sized petri dishes so we had plated onto tiny Kan-R plates with ¼ the surface area of a normal petri dish


July 22nd

AG

- CPCR of regions 2 and 3 with colonies 1-16
 - Region 2 (193bp) blue
 - 52 C
 - Region 3 (140bp) green
 - 50 C annealing temp
 - ladder (100bp)
 - 14uL water
 - 1uL ladder
 - 3uL dye
 - ran at 100V for 30 min (last time had poor resolution)
 - no positive results (probably primer dimers)



ladder/ cpcr region 2 colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8/ 9/ 10/ 11/ 12/ 13/ 14/ 15/ 16



ladder/ cpcr region 3 colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8/ 9/ 10/ 11/ 12/ 13/ 14/ 15/ 16

- Gibson re-do of EBL1 TS pSK33
 - .004794 pmol
 - frag 1: 25 ng = 0.396 uL x
 - frag 2: 11.199 ng = 0.257 uL x
 - frag 3: 17.827 ng = 0.230 uL x
 - frag 4: 17.175 ng = 0.269 uL
 - frag 5: 13.825 ng = 0.125 uL
 - H2O = 3.723 uL x

OE

- Nothing grew on the transformed lambda-red plates
 - Potential problems:
 - Arc time
 - Cells were too dilute
 - Washed away too many cells
 - Bad incubation
 - Going to redo the whole process but after second spin of milli-q, put the solution in 1.5 mL tubes and spin @ max speed
 - PCR Protocol:
 - Take colony and resuspend in 100µL of dH2O
 - Prepare 16 colonies
 - Add
 - H2O to 25 µL
 - 0.5 µL forward primer
 - 0.5 µL reverse primer
 - 1 µL colony suspension
 - 12.5 µL APEX master mix
 - Cycle
 - 95 C for 5 minutes
 - 25 x
 - 95 C for 30 seconds
 - Annealing Temp for 30 seconds
 - 72 C for extention time
 - 72 C for 1 minute
 - 12 C infinite hold
 - F1: length 360 bp anneal @ 50 C extention time 14.4 seconds
 - F2: length 297 bp anneal @ 49 C extention time 11.88 seconds
 - F3: length 332 bp anneal @ 50 C extention time 13.28 seconds
-

July 23rd


OE

- Enterobacter still has not grown
 - We believe that the issue with the digestions performed on our backbone vectors were done incorrectly because we had used both high fidelity enzyme (BamI) and a non high fidelity enzyme (EcoRI). But these require different buffers to operate whereas we had used both enzymes in the buffer that is suitable only for HiFi enzymes.
 - Performed digestions over again
-

July 24th

AG

- inoculated 4x 10mL cultures of DH5a w/ pSK33 (7/23/16)
 - miniprepd all of them
- digest of pSK33 probably didn't work (used BamHI-HF instead of BamHI which has different buffer conditions → needed cutsmart but used NEBuffer 2.1 so BamHI-HF only had 50% activity)
- double digest
 - 5uL of Cutsmart
 - 1uL EcoRI-HF
 - 1uL BamHI-HF
 - 2ug DNA (30uL of pSK33 conc 67.60 ng/uL)
 - H2O up to 50uL (13uL)
 - digest for 2 hours and PCR purify after
- controls:
 - EcoRI single cutter
 - 5uL NEBuffer 2.1
 - 1uL EcoRI
 - 1ug DNA (18uL of pSK33 conc 56 ng/uL)
 - H2O up to 50uL (26uL)
 - digest 2hr and heat inactivate at 65 for 20 min
 - uncut plasmid



ladder/ uncut/ single cut/ double cut

- Gibson Assembly (no FT)
 - .004794 pmol
 - frag 1: 25 ng = 0.396 uL
 - frag 2: 11.199 ng = 0.195 uL
 - frag 3: 17.827 ng = 0.230 uL
 - frag 4: 17.175 ng = 0.269 uL

- frag 5: 13.825 ng = 0.125 uL
- 3.785uL H2O
- 5uL Gibson MM
- 50 C for 1 hr
- inoculate: DH5a w/ ampR; EC93 in EPI100; EC869 in X90; DH5a w/ pSK33; Enterobacter

OE

- PCR of Kan marker for transposon insertion doesn't seem to be working, Jared had asked me to try it on my own
 - Tried the original annealing temperature of 65 C as well as \pm 3C
 - Used 23.88 seconds as the extension time
 - 1.25 μ L forward primer
 - 1.25 μ L reverse primer
 - 1 μ L of DNA insert
 - 12.5 μ L Q5 2X Master Mix
 - H2O up to 25 μ L
- PCR failed again, checked primer design and noticed a presence of both secondary structures and multiple annealing regions
- Redoing Gibson Assembly with the properly digested fragments
 - Used 25ng for the largest fragment (F1)
 - 22.2 ng for F2
 - 20.05 ng for F4
 - 14.9 ng for pSK33
 - 18.64 ng for the toxin fragment

July 25th

AP


- cpcr of omar's EBL2 w/ TS no FT colonies (regions 1, 2 and 3) - 16 colonies each

AC


- redoing competent cells
 - used DH5-alpha control

AG

- cpcr of EBL1 w/ TS no FT (regions 2 and 3) - 16 colonies each
 - same results as before the fixed digestion
- cpcr of kan marker control
 - MM: 105 uL H2O, 5uL F primer, 5uL R primer, 125 uL APEX
 - colonies 1-8



ladder/ cpcr region 2 colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8/ 9/ 10/ 11/ 12/ 13/ 14/ 15/ 16



ladder/ cpcr region 3 colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8/ 9/ 10/ 11/ 12/ 13/ 14/ 15/ 16


JL

- annealing temp gradient for transposon insert

July 26th

AC

- make a stock of the EPI100 w/ EC869 (500uL stock and 500uL glycerol)
- miniprep the rest of the overnight culture (4.5mL) and do a test digest:
EC869 cosmid (expected bands:10kb, 8kb)
 - EcoR1: 1uL
 - H2O to 50 uL (42.3uL)
 - 500 ng DNA (1.7uL of a 300ng/uL stock)
 - 5uL NEBuffer 2.1
 - digest for 1hr at 37C and 65C for 20 min



ladder/ test digest


OE

- transform DH5-alpha with linearized pSK33 to check for background
 - dilute to 1ng/uL (final volume 25uL)
- transform DH5-alpha with puc19
- inoculate EPI100 w/ EC93 (also streak new plate from the stock in the -80C with a V on it for verified)
- inoculate X90 w/ EC869 (use rainin tip)
 - streak new plate if alan's test digest works
- **make stock of DH5-alpha ampR (negative control)**
- inoculate DH5-alpha w/ pSK33
- inoculate DH5-alpha w/ EBL1 TS and EBL2 TS
 - for sequencing

July 27th

AG/OE

- ran cPCR controls
 - used Danny's cpcr as a positive control to make sure our technique was ok (came back positive)
 - miniprep transformed DH5a w/ EBL1 and EBL2 colonies and used CPCR region 1 primers → no results



ladder/ positive control colony 1/ 2/ 3/ 4/ miniprep colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8/ positive control colony 5/ 6/ 7/ 8

July 28th

AC

- PCR Kan R for transposon insert with Jared's primers

AG/OE

- made more kan plates
- made sequencing primers for pSK33 variable region with cut sites

Competition assay

Purpose: determine whether the EC93 cosmid is expressed in EPI100 since it failed the autoaggregation test

Test inhibitor: EPI100 w/ CDI EC93 ampR

Positive control: Hayes strain (X90 w/ CDI EC869) ampR

Negative control: DH5-alpha ampR

Target (CDI-): DH5-alpha with empty pSK33 plasmid kanR

- all cells are in the shaking incubator in culture tubes with labels
- 1. take 0.5mL of the overnight culture (using p1000) and resuspend in 50 mL fresh LB w/ correct antibiotic inside autoclaved 250mL flasks already labeled with tape
- 2. place back in 37 C shaking incubator w/ tinfoil on
 - a. **prepare the target strain 20 minutes before the controls/inhibitor**
- 3. **after 1 hour, take OD600 of all cells**
 - a. should be all around 0.1-0.2
- 4. **after 10 min take another OD600 (target OD600 = 0.35)**
 - a. if cells are growing at very different rates, back dilute faster growing culture to get it down to a lower OD
 - i. resuspend appropriate vol of cells in 50 mL LB w/ appropriate antibiotic
- 5. mix inhibitor (or control) and target cell at a ratio of 10:1 (determine using growth phase/conc relation) in LB (w/o antibiotic)
 - a. add target cell culture to the inhibitor culture

$$vol\ in\ inhibitor\ flask * (x\ cells/mL) = number\ of\ inhibitor\ cells\ in\ the\ culture = n$$

$$n / 10 = number\ of\ target\ cells\ needed = m$$

$$m * (1mL/x\ cells) = vol\ of\ target\ needed\ to\ be\ added$$

****the cells per mL will be different in each equation; depends on the specific OD of the culture**



6. take a before sample (10-fold serial dilutions from 1:1 to 1:1000000)
7. incubate at 37 C shaking at 225 rpm
8. before and after sample ten-fold serially diluted in SOC (before competition assay, every hour until 4th hour, after 20 hours)
 - a. plate 100uL of samples onto selective agar
 - b. count CFUs per mL


July 29th

AG/OE

Results of competition assay

Strain	CFUs/mL	
	t=0h	t=3h
(+) control X90 w/ EC869	1.6 E 6	3.78 E 4
(-) control DH5a ampR	1.53 E 6	1.07 E 8
test EPI100 w/ EC93	1.26 E 6	7.11 E 3





- inoculate for new competition assay
 - DH5-a w/ EC869
 - DH5-a w/ EC93
 - DH5-a w/ ampR

○

Notes:


- for t=0h, use 1:1000 and 1:10000 and 1:100000
- for t=0h, use 1:100, 1:1000 for inhibitors and 1:100000 and 1:1000000 for neg control

- pSK33 map was incorrect due to a 4 base pair frameshift
- order new primers for EBL1 construct w/ updated pSK33 file
- sent in pSK33 for sequencing
 - 2.5 uL of primer, 8.16uL miniprepmed plasmid, 4.34uL water
 - do each primer separately
- redo kan marker cpcr on higher copy plasmid that has kan
 - found out kan marker primers oriented outward and can't be used together

August 1st

AP/JL

- Transposase Activity Assay
 - 1.35 μ l pSB1C3 (148 ng/ μ l)
 - Transposon at 92.45 ng/ μ l
- Gel electrophoresis
 - Control: unreacted pSB1C3




ladder/ reacted/ unreacted

AE/OE

Competition assay results

	CFUs/mL	
	t=0h	t=3h
DH5a w/ EC869	7.33 E 6	1.29 E 4
DH5a w/ EC93	4.44 E 6	2.55 E 2
DH5a w/ ampR	3.28 E 6	1.3 E 8



- inoculate same strains for another trial


August 2nd

OE/AG

Competition Assay results

Strain	CFUs/mL	
	t=0h	t=3h
DH5a w/ EC869	1.6 E 6	3.78 E 4
DH5a w/ EC93	1.53 E 6	1.07 E 8
DH5a w/ ampR	1.26 E 6	7.11 E 3

- PCR of pSK33 backbone and frag 1 and 3 of EBL1-TS and frag 1 and 4 of EBL2-TS
 - KAPA Hifi
 - 25 uL KAPA Hifi MM
 - 1.5 uL F primer
 - 1.5 uL R primer
 - 1 ng template DNA
 - 21.0 uL H2O
 - 95 C 3 min
 - 98 C 30s
 - Tm 15s
 - 72 C 40s/kb
 - 72 C 1min/kb
 - 12 C 1hr
 - AE frag 1: 63 C, 8kb → 5 min 20s ext time
 - AE frag 3: 62 C, 5.8 kb → 3 min 52s ext time
 - pSK33: 66C, 3.3 kb → 2 min 13s ext time



Ladder pSK33 EBL1 1 EBL1 3 EBL2 1 EBL2 4

- Secondary more preferential product being formed in the PCR of pSK33 → one primer was very nonspecific (15bp homology downstream)
- EBL2 frag 4 appeared very strange -- all other bands the correct size

August 3rd

AP/ JL

Transposase Activity Assay 2

1. Prepare the transposon insertion reaction mixture by adding the following, in order:

- 1 µl EZ-Tn5 10X Reaction Buffer
- 0.22 µl Target plasmid (900 ng/µl)
- 0.67 µl Transposon (Kan R at 92.45 ng/µl)
- 7.11 µl sterile water (to 9 µl volume)

- 1 µl Tn5 Transposase
- 10 µl Total reaction volume

2. Incubate the reaction mixture for 2 hours at 37°C.

3. Stop the reaction by adding 1 µl of EZ-Tn5 10X Stop Solution. Mix and heat for 10 minutes at 70°C.

Kit: <http://www.epibio.com/docs/default-source/protocols/ez-tn5-custom-transposome-construction-kits.pdf?sfvrsn=8>


4. Gel electrophoresis-----> Check for transposase activity (See if transposon got inserted in plasmid)

Controls:

- Reaction without transposon
- Target reacted without transposase
- Reaction without target
- Transposon only
- Unreacted target

Lanes:


Ladder/ Reaction/ Reaction without transposase/ Reaction without transposon/ Reaction without target/ pUC19/ Transposon only



*Need to run for better resolution

AC

- PCR with positive controls
 - Tested Miniprep 1, 2 and used stock pTKRED as a positive control



ladder/ control/ miniprep 1/ 2

- Using self-designed primers
 - 1F_LBR + 1R_LBR (Random sequence on pTKRED)
 - 2F_LBR-pTKRED_cPCR + 2R_LBR-pTKRED_cPCR (On SpectR gene)
 - Mix contidions
 - 0.5uL Fwd 10nM Primer
 - 0.5uL Rev 10nM Primer
 - 1.5uL Template (133ng/uL Stock; 2ng/uL Miniprep 1; 3ng/uL Miniprep #2)
 - 12.5 OneTaq 2x Mastermix
 - 10uL Water
 - Thermocycling conditions
 - 94C--->30s
 - 28 cycles of
 - 94C--->30s
 - 52C--->30s (52C for 1F/1R; 54C for 2F/2R)
 - 68C--->30s (1min/kb)
 - 68C--->15min
 - 14C--->Infinite Hold
- Results:
 - Nothing showed up including the positive control
 - Most likely primer design problems

August 4th

JL

- PCR more transposon
- Transposase activity assay Diagnostics

AP

- Formed transposomes (using 117.57 µg/ml transposon)

OE/AG

- pSK33 and frag 4 from EBL2 gel purified (accidentally didn't add ethanol to wash buffer)
 - image saved on computer
- Re-gibsoned the constructs and transformed into DH5a

AC

- Miniprep from overnight inoculation
 - From 2 different colonies
- Ran PCR on minipreps next to stock control
 - 25uL Apex
 - 1uL Forward Primer (Calin's pTKRED Ori primers)
 - 1uL Reverse Primer (Calin's pTKRED Ori primers)
 - 23uL dH2O
 - 1uL of diluted miniprep
 - ~6ng/uL
 - Looking for a band at ~1.2kb




Results: Stock was in the correct location, minipreps did not show up

August 5th

AC

- Prepared cPCR with grown colonies (6 colonies)
 - Created mastermix with primers and Apex

- For 1 colony
 - 0.5uL of 10nM Fwd primer (Calin's Ori primer)
 - 0.5uL of 10nM Rev primer (Calin's Ori primer)
 - 12.5uL of Apex
 - 10.5uL dH2O
- Cycling Conditions
 - Initial denaturation: 3min at 94C
 - 30x of the following
 - 30s----->94C
 - 30s----->49C
 - 1:10----->72C
 - Final extension time 5min at 72C
 - Infinite hold at 14C
- Also prepared to run circular plasmid
 - Use 200ng of stock and minipreps
 - cPCR failed
 - Most likely amplified segments of gDNA
 - Very faint bands show up next to the minipreps, not entirely sure why/ what it is



OE

- Performed a CPCR of colonies grown overnight
 - Tested for the presence of at least 3 fused fragments
 - H2O up to 25 μ L
 - 0.5 μ L forward primer
 - 0.5 μ L reverse primer
 - 1 μ L colony suspension (suspended in 100 μ L dH2O)
 - 12.5 μ L APEX
 - Cycle:
 - 95 C for 5 minutes
 - 27 X
 - 95 C for 30 seconds
 - Annealing Temp for 30 seconds
 - 72 C for extension time
 - 72 C for 5 minutes
 - 12 C hold for 1 hour
- Ran a gel with no hits....

August 7th

OE

- Theorized that the gibson had not been with enough attention to detail and stringent practice
- Redid Gibson with the same parameters but more care and attention

August 8th

AC

- Inoculated Enterobacter Aerogenes colony for use in making electro-competent cells tomorrow

AP

- Formed transposomes (stored at -20 C)
 - 2 µl Transposon DNA (117.57 µg/ml in TE Buffer)
 - 4 µl Transposase
 - 2 µl 100% glycerol
 - 8 µl Total reaction volume

OE

- Performed CPCR testing for the presence of four homology regions across 32 colonies...
 - All but two came back negative and even the two positives were not from the same colony
 - We theorized that perhaps uncut plasmid may be contaminating our samples and that may be why we are getting growth on our plates
 - Solution: DPN1 treat the backbone fragments to remove methylated plasmid (uncut plasmid):
 - Add 10 units of restriction enzyme per 1 µG of DNA insert (20 units for gDNA)
 - 1X CutSmartBuffer
 - Total reaction volume: 50 µL
 - 1hr incubation @ 37 C
 - Heat inactivate @ 80 C for 20 minutes

August 9th

AP

- Electroporated transposomes into Enterobacter
 - Added 1uL of transposomes to 50uL of homemade electro-competent Aerogenes (in 1.5 uL microcentrifuge tube)
 - Gently flicked tube to mix
 - Pipetted mixture into a chilled 10mm gap electroporation cuvette
 - Electroporated at 1.8kV
 - Arc time 3.8
 - Immediately added 950uL of SOC to cuvette
 - Pipette down, up, down, up
 - Transferred to microcentrifuge tube and shook at 800rpm at 30C for 30 min, then transferred to lightly capped culture tube at 200 rpm 37 C for 30 min
 - Centrifuged microcentrifuge tubes at 6000rpm for 4min on tabletop centrifuge
 - Discard supernatant and resuspend pellet in 50uL SOC
 - Plate 50uL of resuspended cells onto warm kan 25 plate
- Incubate overnight at 37C for 16h
- ***Need to:
 - Transform 1 µL transposomes into Dh5 alpha
 - Transform 1 µL pUC19 into Dh5 alpha (negative control)
 - Transform 1 µL pUC19 into Enterobacter (calculate transformation efficiency)

AC

- Made competent Enterobacter aerogenes
 - Take overnight culture and add 0.5mL of culture into 50mL of LB in a 250mL flask
 - Did twice, so have 2x 50mL of LB shaking in 37C
 - Let cells grow until reached an OD of 0.497
 - Immediately place flasks into ice-water slurry and swirl gently for ~5min
 - Moved everything to the cold room at 4C
 - All following steps were performed at 4C or lower
 - Aliquoted culture into 4x 50mL Falcon tubes
 - Spun tubes at 5000rpm on floortop centrifuge
 - Removed LB supernatant
 - Resuspended all pellets in 12.5mL Milli-Q Water

- Combined suspensions into 2 Falcon tubes (final volume in each was 25mL)
- Spin tubes at 5000rpm on floortop again
- Discard supernatant and resuspend in 25mL Milli-Q Water
- Repeat spin cycle and discard supernatant, but resuspend in 1mL of Milli-Q Water
- Transfer to microcentrifuge tubes
- Spin microcentrifuge tubes at 5000rcf on tabletop centrifuge
- Discarded supernatant and resuspended in 30uL of Milli-Q water
- Transformation
 - Added 2uL of 1.33ng/uL pTKRED solution to 30uL of homemade electro-competent Aerogenes (still in microcentrifuge tube)
 - Lightly mixed by flicking bottom of tube
 - Pipet mixture into a chilled electroporation cuvette (between the 2 metal plates)
 - Electroporate at 1.8kV
 - Arc time was 5.5
 - Immediately add 970uL of SOC to cuvette
 - Pipette down-up twice
 - Transfer to microcentrifuge tube and shake at 800rpm at 30C for 1hr
 - After 1hr, spin microcentrifuge tubes at 6000rpm for 4min on tabletop centrifuge
 - Discard supernatant and resuspend pellet in 50uL SOC
 - Plate 50uL of resuspension onto warmed up SpecT 25 plate
 - Incubate overnight at 30C


AG/OE:

- CPCR of region 1, 2, 3 and kan marker control of colonies (DH5a transformed with EBL1 TS and EBL2 andTS)
- Only one positive result for region 1 DH5a w/ EBL1 TS (neg for other two regions)
 - Inoculated for sequencing the next day

August 10th

AP


- No colonies after 16 h



18 h

- Transformed:
 - Dh5 alpha w/ transposomes
 - Dh5 alpha w/ pUC19 (negative control)
 - Dh5 alpha w/ pSK33 (positive control) (arctime 5.3)
 - 1.0-mm gap cuvette, 1.8 kV

Results from transformation w/ Dh5 alpha:16 h



No growth on all except Negative control was a lawn because flipped positive and negative control labels. Positive was actually w/ pSK33.

AG/OE


- Sequenced pSK33 w/ frag 1 attached
- Re-PCR'd pSK33, EBL1 frag4 and EBL1 frag 5 (Omar re-PCR'd all the EBL2 fragments)
 - dpnI treat the pSK33 after gel purification to cleave any leftover methylated template DNA

August 11th

AP

- Drop dialysis of transposase
- Formed transposomes with desalted transposases
- re-Transformed:
 - Dh5 alpha w/ transposomes (arctime 5.1)
 - Dh5 alpha w/ pUC19 (negative control) (arctime 5.00)
 - Dh5 alpha w/ pSK33 (positive control) (arctime 5.10)
 - 2.0-mm cuvette, 2.5 kV

No growth on transposome plates. Positive control fine.



AG

- Sent in (+) CPCR region 1 colony for sequencing using 12F_EBL1-p33-TS_SQ and 21F_EBL1-p33-TS_SQ → no priming
- Gel purified fragment 4 for EBL1
- DpnI treated pSK33 gel purified product and PCR purified
 - a. 9uL pSK33 sample
 - b. 5uL Cutsmart
 - c. 0.5uL enzyme
 - d. 35.5 uL H₂O
 - e. 37 C for ~1.5hr and 80 C for 20 min
- Gibsoned EBL1 and transformed DH5a (plus negative control without Gibson MM)
 - a. Frag 1: 25 ng → 0.449 uL
 - b. Frag 2: 10.34 ng → 0.179 uL
 - c. Frag 3: 17.90 ng → 0.210 uL
 - d. Frag 4: 17.17 ng → 0.191 uL
 - e. Frag 5: 13.83 ng → 0.110 uL
 - f. water up to 5uL
 - g. 5uL Gibson MM
- Inoculated 20mL more of (+) colony
- Streaked (+) colony
- Streaked enterobacter

August 12th

AP


Optimizing Electroporation of Transposomes:

- Form transposomes (JL)
- Desalt transposomes using drop dialysis
- Transform:
 - Dh5 alpha + 1 μL transposomes (arctime 5.2)
 - Dh5 alpha + 2 μL transposomes (arctime 5.1)
 - Dh5 alpha + 1 μL desalted transposomes (arctime 5.1)
 - Dh5 alpha + 2 μL desalted transposomes (arctime 5.1)

- 2.5 kV, 2.0-mm cuvette

AG

- Saw growth on neg control plate still (less than before)
- CPCR of region 3 for 32 colonies
 - No positives



ladder/ cpcr region 3 colony 1/ 2/ 3/ 4/ 5/ 6/ 7/ 8/ 9/ 10/ 11/ 12/ 13/ 14/ 15/ 16
ladder/ 17/ 18/ 19/ 20/ 21/ 22/ 23/ 24/ 25/ 26/ 27/ 28/ 29/ 30/ 31/ 32

Troubleshooting Ideas

- Try PCR of pSK33 with 0.1ng and 0.5ng of plasmid template in a 50uL PCR rxn and up cycle number to 30
- Digest PCR product with DpnI restriction enzyme to destroy plasmid template
 - NEB protocol:
 - 8uL PCR rxn
 - 1uL DpnI
 - 1uL Cutsmart
 - Incubate at 37C for 30 min and heat inactivate at 80C for 20 min
- Extended Gibson incubation shown to help efficiency
- PCR amplify Gibson product (RCA)
- Assemble two smaller plasmids and combine them


August 13th

AG/OE

- PCR of pSK33 with 0.1 ng and 0.5 ng of template instead of 1 ng
 - DpnI treated
 - gel purified

AP/ JL

16h at 37 C



- Desalted transposome 1 ul (arctime 5.10)
- Desalted transposome 2 ul (arctime 5.10)
- Transposome in regular buffer 1 ul (arctime 5.20)
- Transposome in regular buffer 2 ul (arctime 5.10)

No growth on all

Troubleshoot:

1. *Transposase concentration from Nate too low (Play around with 4ul in forming transposomes)*
2. *Transposomes/ transposon not forming properly→ assay for function by transforming transposed pUC19*
3. ~~Transposome not getting in (since 55kda per unit transposase)~~
4. ~~Transposome with multiple insertions in gDNA (kills it)~~

August 16th

AG

- inoculate positive colony to see if the toxin region is assembled

AP/ JL

- React transposase with pUC19, PCR purify
- Nanodrop concentration:
- Transform 1 ul with Dh5 alpha
 - 2.5 kV, arctime 4.8, 2.0-mm gap cuvette
 - Plate on kan 25

OE

- PCR amplify fragments together such that you have a new fragment consisting of the old fragments 1 and two. The same was done with fragments 4 and 5 and an additional fragment with 3, 4 and 5.
 - The idea is to try and PCR fuse as many fragments as possible to reduce the number of fragments in the gibson
 - Because fragment 3 is so much smaller than the other fragments, it may have to be added to the end of the fused fragment 4 and 5 if we cannot PCR amplify the three combined in the first place

August 17th

AG/OE

- Sent in pos colony for sequencing (primer on toxin region)
- Single cutter digestion of plasmid in pos colony (using NcoI-HF) to check for plasmid size
 - Expected: 11-12kb
 - Gel image saved on AGuo
- Fusion PCR step in protocol using already formed EBL2 frag 4 and DpnI treated pSK33 (made w/ 0.1 ng of template)
 - Ran 3 different reactions:


- Rxn 1:
 - 1ng F4, 0.7ng pSK33, 10.8uL H2O, 12.5uL Q5
- Rxn 2:
 - 2ng F4, 1.38ng pSK33, 9.12uL H2O, 12.5uL Q5
- Rxn 3:
 - 5ng F4, 3.45ng pSK33, 4.05uL H2O, 12.5uL Q5
- Ext time: 5min 28s
- Final ext 2 min
- 15 cycles
- Gel extract and PCR using Gibson primers (4F and MS 1R)

JL

- React formed transposomes + pUC19
 - 1.5 µl EZ-Tn5 10X Reaction Buffer
 - .33 µl target DNA (.2 µg)
 - 11.2 µl sterile water
 - 2 µl EZ-Tn5 Transposome
 - 15 µl Total reaction volume
- Transform that (assay if pre-formed transposomes are viable)



AP



Very few tiny colonies <kan25>

More colonies than 1:1 (???) <kan25>

More colonies/Satellites than kan plates <amp 50>


Fungal Contamination → threw away amp plates

- Make amp plates
- Tomorrow: Re transform
 - Dh5 alpha + transposed pUC19 (amp 50) positive control
 - Dh5 alpha + transposed pUC19 (kan 25)

AC

- Prepared electro-competent aerogenes
 - Took overnight culture and added 0.5mL to 50mL LB (Did 2x of this)
 - Shake in 37C until reached OD of 0.463
 - Take both flasks out and place in ice water slurry

- Swirl flasks in slurry for 5min
 - Let sit in slurry in the 4C cold room for 15min (All following work was performed at 4C or on ice)
 - Transfer culture to 4 falcon tubes (put 25mL of culture in each)
 - Spin tubes for 20min at 5000rcf
 - Remove supernatant
 - Resuspend each in 12.5mL Milli-Q Water at 4C
 - Mix suspensions so that 4x falcon tubes with 12.5mL each become 2x falcon tubes with 50mL each
 - Centrifuge at 5000rcf for 20min
 - Discard supernatant
 - Resuspend in Mili-Q water at 4C
 - Repeat this step once more
 - Resuspend pellet in 1mL water
 - Transfer into 1.5mL microcentrifuge tube
 - Centrifuge at 5000rcf for 20min
 - Discard supernatant
 - Resuspend using 30uL of Milli-Q water
- Electroporation
 - Add 2uL of pTKRED to 30uL of cells
 - Transform in 0.1cm cuvette
 - Arc time was 5.6
 - Add 2uL of pTKRED to 50uL of cells
 - Transform in 0.2cm cuvette
 - Arc time was 5.9
 - For both transformations
 - Immediately resuspend in 970uL of SOC
 - Incubate for 1hr at 30C
 - Centrifuge at 6000rpm on tabletop for 4min
 - Resuspend in 50uL of SOC
- Plating
 - Plate on Spect 25 plates
 - Use all of the 50uL SOC to plate
 - Incubate overnight at 30C



August 18th


AP

- Re transform Dh5 alpha + transposed pUC19
- Arc time: 5.4, .2 cm cuvette
- 1:10 and 1:100 on kan25 and amp100






Results: 16h 20min at 37C

AC

- Performed cPCR
 - Suspend 15 colonies in 30uL of dH₂O
 - For 1 reaction:
 - 0.5uL of Forward Primer
 - 0.5uL of Reverse Primer
 - 12.5 Apex
 - 1uL colony suspension (one PCR tube was reserved for pTKRED plasmid stock to ensure PCR was not faulty)
 - Fill rest up to 25uL with dH₂O
 - Cycling Conditions
 - 94C----->3min
 - 30x
 - 94C----->20s
 - 49C----->20s
 - 68C----->1min 10s
 - 68C----->5min
 - 14C----->Inifinite Hold
- Ran PCR products on 1% Agarose Gel with 1kb ladder and Sybrsafe



- Well 1 is ladder
- Well 2 is just pTKRED plasmid
- Wells 3-17 are colony PCRs
- Most likely PCR failed as not even the one with just the plasmid showed up
 - Will try different cycling conditions
 - Also do 50uL reaction as no bands showed up




August 19th

AC

- Redid cPCR of pTKRED
 - Screened 7 colonies
 - For 1 reaction:
 - 25uL Apex
 - 1uL Forward Primer
 - 1uL Reverse Primer
 - 23uL Water
 - ½ of a colony added directly to PCR tube (Add 2ng pTKRED plasmid-stock to check PCR works and is fine-in a separate PCR tube)

- Cycling Conditions
 - 94C----->5min
 - 30x of
 - 94C----->30s
 - 49C----->30s
 - 68C----->1min 10s
 - 68C----->5min
 - 14C Infinite Hold
- Run PCR with 1kb ladder in 1% Agarose Gel with Sybrsafe




Well 9 is just pTKRED; Wells 2-8 are cPCR of colonies; 2nd ladder from bottom is 1kb

- Expected to see 1kb fragment same size as well 9
 - Lots of different bands showed up, most likely gDNA or other DNA

JL

Transposome activity assay (unsuccessful)



From left: ladder, normal reaction, 2X transposase, 3X transposase, pUC19 control

August 20th

OE

- Ran a gel on the PCR conjoined fragments and nothing assembled correctly
-

August 22nd

AC

- Inoculate 3 colonies overnight in 5mL LB with 2.5uL of SpecT antibiotic

OE

- Ran a gel for the the gibson a
-

August 23rd

AG

- Sequencing results back → exonuclease chewed back the strand enough s.t. there was compatible region between 3' end of pSK33 and 5' end of toxin region → ~140 bp deletion on toxin region (3 fragment plasmid made)
 - Order 5 primers to gradually add back toxin sequence (use KAPA)
 - try doing 5 cycles matching annealing temp of actual primer and 25 cycles matching annealing temp of primer + overhang if the product doesn't form with the long primer
-

August 26th

AG

- PCR using psk33_fix_1R and EBL1_toxin_fix_1F
 - 12kb fragment successfully amplified using KAPA




- Gibsoned 11kb frag
 - CPCR region 3 amplify
 - 1uL of gibson
 - 0.5uL F
 - 0.5uL R
 - 12.5 Taq 2X MM
 - 10.5uL dH2O
 - Saw faint band at the right size



cpcr region 3/ 100bp ladder

- PCR fragment using gibson as template
 - First five cycles at 51C Tm and 20 cycles at 64C
 - No band at 11kb... probably bad primers
 - ordered new forward primer
- PCR using psk33_fix_1R and EBL1_toxin_fix_2F
 - Did one normal KAPA at Tm = 51 C
 - Did one with the first 5 cycles at 51 and 20 cycles at 65
 - worked better




ladder/ 11kb fragment/ normal KAPA protocol/ increased Tm protocol

August 29th

AC

- Ran gel with EBL2 gDNA primers to verify Aerogenes



Wells

#1: 1kb ladder


#2-6: AG wells

#7-9: gDNA extraction from 3 different colonies

Only 2 colonies were Aerogenes

AG

- PCR using psk33_fix_1R and EBL1_toxin_fix_3F
- PCR using new 11kb frag primers




1kb ladder/ toxin fix/ 11kb frag

August 30th

AG

- Inoculate 15 colonies from devices 1-3 and pos/neg control in 5mL LB and 5uL cat for interlab study
- EBL1 toxin fix 1R + 4F didn't work → saw large smear instead of band
 - troubleshoot tmr




August 31st

AG

- Troubleshoot PCR w/ EBL1_toxin_fix_1R and EBL1_toxin_fix_4F
 - Use new dH2O
 - Try reaction mix without template DNA (check for carryover contamination)
 - Try going back down to 25 cycles (did 28 cycles last time) to check for overcycling

- Try 5 cycles annealing temp of annealing region and 20 cycles annealing temp of entire primer including overhang (54C → 68C) to increase specificity
- Try running PCR using CPCPCR region 1 primers and OneTaq with template DNA to verify correct sequence is in there
 - Results:
 - Gel image under AGuo (1kb ladder; 25 cycles 54C; no template control; 25 cycles 54 → 68 C; CPCPCR region 1; 100bp ladder)
 - Two smears for 25 cycles 54 C and 25 cycles 54 → 68C, faint band at correct size observed for both
 - CPCPCR region 1 positive → template is fine
 - No template control has only primer dimers → no contamination
 - **4F primer is the problem**



1kb ladder/ 25 cycles 54C/ no template/ 5 cycles 54 → 20 cycles 68C/ CPCPCR primers with template/ 100bp ladder