


vector
 Φ_{RE}
 digest

RE fragment
 Φ_{RE}
 digest


$1ng$ of $3KB \approx 1x$
 vs
 $1ng$ of $1KB = 3x$



\Downarrow PCR
 cleanup



\Downarrow PCR
 cleanup



$3KB$



$1KB$

\sim Primers & GreI
IGEM website or GreI

\sim Kan or chl

[DNA] what ratio 1:3 or 1:5
 of vector to
 insert (molar ratio)



x bp $\times \frac{mol}{650-660nm \text{ per bp}}$

U

U



vector



fragment

1 μ g - 2 μ g of
DNA for
digestion


~~BBR1 MCS.2~~

BBR1 MCS.2

iGEM

- Kan

- chl



06 Jun 2018

Competent cells procedure:

made 1M stock /L NaOAc

created CCMBSO buffer

creates SOB

8 tubes

~~Observations~~
Observations

- pH 7.01

- pH 6.41

- pH 7.50

4 PCR products each in tetr + Hendingeyn
Artificial Sea water

U	2	Per Liter
		800 mL Deionized Water
110	55	27.5g NaCl
15.2	7.6	3.8g $MgCl_2 \cdot 6H_2O$
27.12	13.56	6.78g $MgSO_4 \cdot 7H_2O$
2.88	1.44	0.72g KCl
2.48	1.24	0.62g $NaHCO_3$
11.16	5.58	2.79g $CaCl_2 \cdot 2H_2O$
4	2	1g NH_4Cl
2	0.10	0.05g K_2HPO_4
4 mL	2 mL	1 mL Wolfe's Trace Mineral Soln.

PH \approx 6.1 - 6.5

bring to 1 Liter w/ Deionized H₂O

$$mmol/kg \rightarrow mM \quad m \frac{moles}{liter}$$

$$1.025 kg \text{ Sea } H_2O \rightarrow 1$$

$$1.025 kg = 1L \text{ Sea water}$$

$$1 \frac{mmol}{kg} \left| \frac{1.025 kg}{1L} \right| = 1.025 mmol/L$$

Concentration in (mM)

Na ⁺	K ⁺	Cl ⁻
481.75	10.455	563.75

$$Na^+ \quad \frac{470 mmol}{kg} \left| \frac{1.025 kg}{1L} \right| = 481.75 mM$$

$$K^+ \quad \frac{10.2 mmol}{kg} \left| \frac{1.025 kg}{1L} \right| = 10.455 mM$$

$$Cl^- \quad \frac{550 mmol}{kg} \left| \frac{1.025 kg}{1L} \right| = 563.75 mM$$

Making Chloramphenicol Plate

1) get 1 15mL Sterile tube.

2) weigh out 50mg of Chloramphenicol

0.05090g 50.90mg

$$50 \text{ mg} / 2 \text{ mL} = 25 \text{ mg/mL}$$

$$53.7 \text{ mg} / x \text{ mL} = 25 \text{ mg/mL}$$

$$x = 2.15 \text{ mL}$$

3) cool to 55°C

4) add chloramphenicol to LB Agar

5) Plate

Making LB Agar (Per liter)

- 10g of NaCl

- 10g tryptone

- 5g yeast extract

- 15g agar

- 1L of deionized H₂O

- Autoclave 1 hr. fluid

Test Plates

10 ng/mL

(50 ng/mL) \Rightarrow 5mL of nuclease, fresh water.

$$50 \text{ ng} / 5 \text{ mL} = 10 \text{ ng/mL}$$

$$51.7 \text{ ng} / x \text{ mL} = 10 \text{ ng/mL}$$

$$x = 5.17 \text{ mL}$$

1) cool LB Agar to 55°C

2) add the test into LB Agar

3) Plate

8/1/16

Spinning down primers.

- 1.) Max speed for 5 min (cent)
- 2.) SE by drate w/ nuclease-free H₂O to total 100μL.

8/2/16

PCR

For 1 rxn tube (50μL)

- 25μL Q5
- 0.25μL Forward primer
- 0.25μL Reverse primer
- 0.5 μL DNA
- 24 μL Nuclease-free H₂O

Setting up A and B primers.

A1

- 25μL Q5
- 0.25μL F+Pr
- 0.25μL Rev Pr
- 0.5μL DNA
- 24μL NF H₂O

A2

- 25μL Q5
- 0.25μL F+Pr
- 0.25μL Rev Pr
- 0.5μL DNA
- 24μL NF H₂O

* Refer to Fitz's Note book pg 154

8/3/16

Gel



Caddy's stuff.

I Gem


A = CP011927.1 / 729


Loading

- 5 μL of ~~Primer~~ primer
 - 1 μL of Purple dye
 - 7 μL of 10 kb ~~ladder~~ ladder
- } mix on Parafilm


well

- 6 - 10 kb ladder
- 7 - X BF Blank
- 8 - X BF1
- 9 - BF2



Gel #2

- A gel was run to see if ~~we~~ the sections of DNA we desired existed in Arc A and Arc B.
- once we figured out that the DNA we desired existed and showed a strong band, we did a PCR clean up. (see Fitz's lab notebook)




Ladder, blank a1, a1, blank a2, a2 (which evaporated), blank b, b1, b2, a2, again (still didn't show)

8/4/16

Gel 3

we rerun A2 and X2 (blank of A2) b/c it evaporated last time.



Ladder, blank a2, a2


ligation - refer to Fitz's notebook

8/8/16

Transformations

Retrieve competent cells from -80 freezer.

- 2) let it thaw for 30 min in ice.
- 3) After ligation, transformation occurred w/ the following protocol:
- 3) Remove competent cells from -80°C freezer
- 4) thaw on ice (20-30 min)
- 5) add 10 pg - 100 ng (1-5 μ l) of DNA to 50 μ l of competent cells (used 3 μ l of each)
- 6) incubate on ice for 30 min.
- 7) Heat shock at 42°C for 3 min
- 8) cool on ice for 5 min
- 9) Add 750 μ l of LB
- 10) Incubate and Shake for 1 hour
- 11) warm selection plates (37°C)
- 12) spread 100 μ l onto 1 plate (labeled C-A1 and C-B1) and spin down the rest @ 8000 xg for 5 min.
- 13) pour off supernatant, and resuspend pellet in 100 μ l of LB and plate (labeled C-A2 and C-B2)
- 14) Incubate plates at 37°C overnight.



6/8/16

Prep for Arc B

Set up for a PCR genes F and H using the high Fidelity Q5 polymerase and 50 μ L reactions (both will use heterostaggered primers)

F1

25 μ L Q5
0.25 μ L Fwd Pr. 1
0.25 μ L Rev Pr.
0.5 μ L DNA
24 μ L Nuclease Free H₂O

F2

25 μ L Q5
0.25 μ L Fwd Pr. 2
0.25 μ L Rev Pr.
0.5 μ L DNA
24 μ L N-F H₂O

H1

25 μ L Q5
0.25 μ L Fwd Pr.
0.25 μ L Rev. Pr. 1
0.5 μ L DNA
24 μ L N-F H₂O

H2

25 μ L Q5
0.25 μ L Fwd Pr.
0.25 μ L Rev. Pr. 2
0.5 μ L DNA
24 μ L N-F H₂O

Used the T_m calculator to find the temp for the Q5, which for all 4 is 72 °C (F is 1446 bp and H is 1489)

98 °C 30 sec
98 °C 10 sec
72 °C 30 sec
72 °C 1:30 sec
repeat x 34
72 °C 3 min
10 °C ∞

8/9/16 Gel run for Arc B

[1kb] [XH] [H1] [H2] [] [XF] [F1] [F2] [] []

- A gel for the F and H gene was run after the PCR. w/ 1 μ L of purple dye mixed w/ 5 μ L of each component above.

No results, so we reran the gel. this time w/ 1 kb.

1 kb

1 μ L 1 kb ladder

4 μ L H₂O (N-F)


1 μ L dye (purple)

- again no result.

[1kb] [XH] [H1] [H2] [XF] [F1] [F2]

- b/c No result we re did the PCR later to previous page for PCR procedure.

Gel #2



Still Zero bands.


↳ must redo PCR. (again)

8/10/16 Gel run for Arc B (attempt #3)

[1kb] [XF] [F1] [F2] [XH] [H1] [H2]
1 2 3 4 5 6 7

F bands appeared; however, H is still resisting amplification, but since F showed up this time maybe it will work next time?

redo PCR w/ gene H.



The previous gel run on genes F and H was showed results for only gene F, so we will be PCRing again for gene H and as well as genes A and B because our last transformation produced no results, so we used up all our DNA.

★ N-F H₂O =
Nuclease-free
H₂O.

For A (hetero staggered)

A1) 25 μ L Q5
0.25 μ L Fwd Primer 1
0.25 μ L rev primer
0.5 μ L DNA
24 μ L N-F H₂O

A2) 25 μ L Q5
0.25 μ L Fwd Pr. 2
0.25 μ L Rev Pr.
0.5 μ L DNA
24 μ L N-F H₂O

ten (71°)

tm (72°) = A1, B, H1, H2

B

25 μ L Q5
0.25 μ L Fwd pr
0.25 μ L Rev pr
0.5 μ L DNA
24 μ L N-F H₂O

H (hetero staggered)

H1) 25 μ L Q5
0.25 μ L Fwd pr
0.25 μ L Rev pr 1
~~24~~ 0.5 μ L DNA
24 μ L N-F H₂O

H2) 25 μ L Q5
0.25 μ L Fwd pr
0.25 μ L Rev pr 2
0.5 μ L DNA
24 μ L N-F H₂O

PCR Protocol

98°C 30 sec.

98°C 10 sec

gradient \rightarrow 71°/72°C 30 sec


72°C 1:30 sec

repeat x34

72°C 3min

10°C \checkmark

1 kb	XA	A1	A2	XB	B	XH	H1	H2	
	2	3	4	5	6	7	8	9	10




A and B - showed bands; however, H still gave no results.

- run gene H PCR again, but this time use different temperatures.

↳ set up a gradient to experiment w/ multiple different temperatures.


PCR thermocycler has 8 rows, so we can put 8 different temp.

gradient should be:



need to make total of 17 tubes.

- 1 - blank
- 8 - H1 label them
- 8 - H2



8/11/16 PCR again! (only for gene H [ancB])

TM gradient

H1 - 65°C	H5 - 69.2°C
H2 - 65.5°C	H6 - 70.7°C
H3 - 66.4°C	H7 - 71.5°C
H4 - 67.7°C	H8 - 72.0°C

H (hetero staggered)

HA	x8
25 mL - QS	200 mL
0.25 mL - Fwd Pr	2 mL
0.25 mL - Rev Pr 1	2 mL
0.5 mL - DNA	4 mL
24 mL - N-F H ₂ O	192 mL

HB	x8
25 mL - QS	200 mL
0.25 mL - Fwd Pr	2 mL
0.25 mL - Rev Pr 2	2 mL
0.5 mL - DNA	4 mL
24 mL - N-F H ₂ O	192 mL

Blank = 49.5 mL b/c NO DNA.

↓
total for each tube = 50 mL.

Gel #4

1 kb	KT	H _{A1}	H _{A2}	H _{A3}	H _{A4}	H _{A5}	H _{A6}	H _{A7}	H _{A8}
1 kb	KT	H_{A1}	H_{A2}	H_{A3}	H_{A4}	H_{A5}	H_{A6}	H_{A7}	H_{A8}

No results, however, the agar was very weak, so tomorrow I will re-run the gel w/ fresh Agars.

* ASK Brian to see if our genes show up on the Blast

6/12/16

Cell growth

A 10 ~~ML~~ μ L - No growth
 * A 100 μ L - 16 colonies.
 A 400 μ L - 88 colonies.

B 10 μ L - 5 colonies
 * B 100 μ L - 43 colonies
 B 400 μ L - 150 colonies

Colony PCR

12.5 μ L green mix
 0.5 μ L primer #1
 0.5 μ L primer #2
 11.5 μ L N-F H₂O
 1 colony.

Master Mix

100 μ L green mix
 1 μ L primer #1
 1 μ L primer #2
 92.1 ~~94.1~~ μ L N-F H₂O

25 μ L per rxn - (Put 25 μ L in each tube)

A

	1	2	3
D	\checkmark _S	\checkmark _S	\checkmark _M
E	\checkmark _M	\checkmark _M	\checkmark _L
F	\checkmark _L	\checkmark _L	

B

	1	2	3
D	\checkmark _S	\checkmark _S	\checkmark _M
E	\checkmark _M	\checkmark _M	\checkmark _L
F	\checkmark _L	\checkmark _L	

PCR Protocol

1.) 95°C for 3:00

2.) 95°C for 0:30

3.) X°C for 0:30. (1-2°C lower than lowest T_m of primers)

4.) 72°C for 1 kb/min (1:00 for 1000 bp for 1000 bp Plasmids)

5.) 2-4 34X

6.) 72°C for 2:00

7.) 10°C

X°C for A = ~~66.8~~ 66.8°C

X°C for B = 67.6°C

8/12/16


Marking Agarose300ml294ml of distilled H₂O

6ml of 50X TAE

1g / 100ml of agarose for 1%.

This is the

ladder

measurement used
for the gel
~~but~~ runs
from this point
forward

- boil until clear and agarose dissolves.

- cool to "uncomfortably warm"

- add EtBr (toxic) wear gloves. - 15ml

8/15/16



Cloning PCR gel run for B

- the bands are barely visible.

* = bands.

~~I read it~~

- did not add primers - will re-do.



- good clear bands.

- move on to next step.

D1, D2, E2 - best

Second best - (D3, F1A)



Qrel run for gene H (gradient)

Water | H¹ AB | H² B | H³ B | H⁴ B | H⁵ B | H⁶ B | H⁷ B | H⁸ B | HX B

still nothing


will have to do colony PCR using Nick's Labrenzia colonies.

- Ask Nick for Colony.

- Question taken, -

what should I do with
the grown colonies of gene
A and B in the Kan
vector?

Colony PCR redo for B



Colony PCR

12.5 mL green mix
0.5 mL B fwd pr
0.5 mL B rev pr.
11.5 mL N-F H₂O
1 colony.

Master mix

100 mL green mix
4 mL PR fwd
4 mL REV PR.
92 mL N-F H₂O.
25 mL in each
tube.

PCR Protocol

~~See Page 15.~~

~~Step 3~~ X °C for B = 67.6 °C
for step 3

$$\begin{array}{r} 67.6^\circ\text{C} \\ - 1 \\ \hline 66.6^\circ\text{C} \\ \approx 66^\circ\text{C} \end{array}$$

x34

- 95 °C — 3 min
- 95 °C — 30 sec
- X °C — 30 sec
- 72 °C — 30 sec
- 72 °C — 2 min
- 10 °C — ∞

6/16/16

PCR clean up

Run PCR clean up for gene F.

clean-up on F1, and F2 (heterostaggered)

F1 - 40 μ LF2 - 40 μ Lratio = $\frac{1}{5}$ $\begin{matrix} \leftarrow \text{PCR} \\ \leftarrow \text{Pb} \end{matrix}$

$$\frac{1}{5} = \frac{40}{x} \Rightarrow 200 = x$$

200 μ L of Pb -40 μ L of F -

PCR

clean up

→ F1

40 μ L F1200 μ L Pb

spin

discard

750 μ L PE

spin

discard

spin

discard

switch to collection tube

add 50 μ L EB

spin

toss filter

F2

40 μ L F2200 μ L Pb

spin

discard

750 μ L PE

spin

discard

spin

discard

add 50 μ L EB

spin

toss filter

Nano drop

	conc. (ng/ μ L)	260/280	260/230
F1	29.4	2.11	2.92
F2	28.4	2.09	3.65

Since heterostaggered - combine F1 and F2 into single tube label it F

Now a digest will be performed.

Get the previously extracted chloramphenicol plasmids from freezer (see F2's pg 42)

~~67.7 ng~~ \Rightarrow ~~1.00 ng~~ 14.8 μ L of plasmid

8/16/16

Recipes for Digest (500ng)

$$F_1) \frac{29.4 \text{ ng}}{1 \text{ mL}} = \frac{500 \text{ ng}}{X \text{ mL}} \quad X = 17.0068027 \text{ mL of } F_1$$

$$F_2) \frac{28.4 \text{ ng}}{1 \text{ mL}} = \frac{500 \text{ ng}}{X \text{ mL}} \quad X = 17.6056338 \text{ mL of } F_2$$

total = 1000ng

Since ratio = 1:3 for ligation, make 2 F_1 and F_2 mixtures.
1 for the chlor vector and one for the Kan vector.

Cl Plasmid

$$\frac{1000}{3} = 333.333 \text{ ng} \quad \frac{67.7 \text{ ng}}{1 \text{ mL}} = \frac{333.33 \text{ ng}}{X \text{ mL}} \Rightarrow X = 4.9236829 \text{ mL} \approx 4.92 \text{ mL}$$

Plasmid Cl

4.92 mL DNA

5 mL buffer (Cutsmart)

1 mL RE1 (EcoRI)

1 mL RE2 (PstI)

37.08 mL N-FH₂OF17 mL of F_1 17.61 mL of F_2

5 mL of (Cutsmart) buffer

1 mL of ~~RE1~~ (PstI)9.39 mL of N-FH₂O.

K Plasmid

$$\frac{1000}{3} = 333.33 \text{ ng} \quad \frac{12.7 \text{ ng}}{1 \text{ mL}} = \frac{333.33 \text{ ng}}{X \text{ mL}} \Rightarrow X = 26.246719 \text{ mL} \approx 26.3 \text{ mL}$$

Plasmid K

26.3 mL DNA

5 mL buffer (Cutsmart)

1 mL RE1 (EcoRI)

1 mL RE2 (PstI)

need 333.33 ng of K

$$12.7 \text{ ng/mL} \times ? = 333.33 \text{ ng}$$

X =

No need to digest K, already cut (check Igem box) (it says "cut")

for ligation


→ 26.3 mL of K

F

look above for Recipe

8/16/16 gel run colony PCR.

[ladder] [D₁B] [D₂B] [D₃B] [E₁B] [E₂B] [E₃B] [F₁B] [F₂B] []




Best band is E₃B colony PCR
Second - D₃B colony PCR.

This gel ran much better
than the initial IgM B colony
PCR run.

✶ — Digestion complete No do a PCR cleanup and
place in freezer for ligation and transformation
tomorrow.

ladder D1B D2B D3B E1B E2B E3B F1B F2B



Best band is E3B colony PCR
Second - D3B colony PCR.

This gel ran much better
than the initial E3B colony
PCR run.

* Digestion complete No do a PCR cleanup and
place in freezer for ligation and transformation
tomorrow.

Nanodrop

	conc (ng/ul)	260/280	260/230
C1	5.2	2.4	1.61
F1	11.4	2.2	1.58
F2	12.4	2.16	1.80

* tomorrow


ligate
transform

into

both

Kan

Chlor^r



7/16 ligation and transformation.

- heated and slowly cooled the F and C1 tube to get the correct sticky ends from the heterostaggered PCR. Put F in the thermal cycler for the following.

95°C 5min.
 85°C 30sec
 75°C 30sec
 65°C 30sec
 55°C 30sec
 45°C 30sec
 35°C 30sec
 25°C 30sec
 15°C 30sec
 5°C ∞

ligate K-F

2 mL T4 buffer
~~4.81 mL Plasmid C (10.8 ng/μL)~~ 2.3 μL V_K
 2.82 μL of F2 (12.4 ng/μL)
 11.88 μL ~~of H2O~~
 1 μL ligase

After this a ligation is performed for C-F and K-F

ligate E-F

2 mL T4 buffer
 4.81 μL Plasmid C (10.8 ng/μL)
 3.07 ~~3.82~~ μL F1 (11.4 ng/μL)
 9.12 ~~11.88~~ μL H2O
 1 μL ligase.

★ Add ligase
Last!!!

Amount of Plasmid needed for 20 μL ligation is 25 ng.

$$(10.8 \text{ ng}/\mu\text{L}) V_K = 25 \text{ ng}$$

$$V_K = 2.3 \mu\text{L}$$

Also assumed that since ~35 ng of insert B was needed in Kirstie's notebook, we will still need ~35 ng.

F

$$F_1 \rightarrow (11.4 \text{ ng}/\mu\text{L})$$

$$F_2 \rightarrow (12.4 \text{ ng}/\mu\text{L})$$

$$C_1 \rightarrow (5.2 \text{ ng}/\mu\text{L})$$

$$\frac{35 \text{ ng}}{11.4 \text{ ng}/\mu\text{L}} = 3.07 \mu\text{L of } F_1$$

$$\frac{35 \text{ ng}}{12.4 \text{ ng}/\mu\text{L}} = 2.82 \mu\text{L of } F_2$$

$$\frac{25 \text{ ng}}{5.2 \text{ ng}/\mu\text{L}} = 4.81 \mu\text{L of } C_1$$

for chlor vector ligation/transformation

for kan vector