

# Week 10

August 8, 2016

What We Did Today:

Digested pgRNA with BstX1 and xho1 for Gibson Vector

Sequencing

pick 4 colonies of WT and PTC colonies from each Gibson plate

- Split cell culture

Confluence: >90%

1:5 split

- Cell plating

Cell count: 339; 371 → Average = 355

- $1.42 \times 10^6$  cells/mL

Plating:


12-well plate:  $3 \times 10^5$  cells \* 13 wells → 2.75 mL cells

Purified pSB1C3

August 9, 2016

What We Did Today:

Gel purified digest from overnight  
Prepped art for A term



Gel – pgRNA Vector

Lane 1 –ladder  
Lane 2 – Empty  
Lane 3 – 30ul Sample  
Lane 4 – Empty  
Lane 5 – 30ul Sample

- Transfection  
APOBEC trial 2
    - Row 1: pRetro empty
    - Row 2: pRetro-APO2X-dCas9
    - Row 3: pRetro-APO3X-dCas9
  - Cell plating  
Cell count: 127; 112 → Average = 120
    - $0.48 \times 10^6$  cells/mLPlating:  
12-well plate:  $3 \times 10^5$  cells \* 13 wells → 8.125 mL cells
- Miniprep

Digest with Not1 to linearize  
neg control GFPn3  
16uL water, 1uL each enzyme, 2uL buffer  
gel:  
L, neg, WT1-4, PTC1-4

Test digested biobrick ATG/ACG eGFP, run on gel

# August 10, 2016

What We Did Today:

Western Blot  
Gibson gRNA

Western Blot Protocol

Used premade gels

15  $\mu$ L of each sample added to each well

Ran for an hour and a half on 120 V

Took gel out of plastic container and stored in water

In a Pyrex dish containing **Transfer Liquid**, made stacks using two sided clasp, sponge, filter paper, gel, transfer paper, filter paper, sponge and closed clasp, being careful to keep it always wet with transfer liquid and without air bubbles.


Put in blot box filled with transfer liquid, stir bar, and ice box. Ran at 100 V for 1 hour.

Removed from blot box, took out the membrane and put in water, added **Ponceau S Stain** liquid, drained liquid, rinsed with water and photographed.


Rinsed with **Tween Wash Buffer** until clean.

Moved membrane into bag, added 10 mL **Dry milk and Tween PBS wash buffer**.

Removed bubbles and sealed bag. Refrigerated overnight.





**Transfer Liquid**

50 ml 10/20x  
100 ml meOH  
850 water

**Ponceau S Stain liquid**

95 ml water added to 0.1g Ponceau S Power  
5 ml glacial acetic acid added last

**Tween Wash Buffer**

100 ml 10x PBS  
900 water  
5 ml 20% Tween 20

**Dry milk and Tween PBS wash buffer**

2g dry nonfat milk  
40 ml Tween Wash Buffer

**Gibson Protocol**

Scrambled pgRNA  
9ul pgRNA vector  
1ul Insert  
10ul Gibson mix

**ACG-GFP pgRNA**

9ul pgRNA vector  
1ul Insert

10 ul Gibson mix

Globin pgRNA  
9ul pgRNA vector  
1 ul Insert  
10ul Gibson mix

Negative Control  
9ul pgRNA vector  
10ul Gibson mix  
1ul Water

- Induction

Doxycycline → APOBEC plate

- Column 1: 0
- Column 2: 10 ng/mL
- Column 3: 25 ng/mL
- Column 4: 50 ng/mL

Transfected globin into cells

Sequenced WT/PTC

Rerun ACG eGFP biobrick samples on gel to confirm

August 11, 2016

### **Gibson Plate Counts/Liquid Cultures**

Neg. Control #1,2,3 – 0

Scrambled pgRNA

#1 – 1

#2 – 4

#3 – 0

Globin #1,2,3 – 0

ACG-GFP pgRNA

#1 – 0

#2 – 1

#3 – 1

### **Western Blot**

Took blots out of overnight milk

Added 15 mL milk and 1.5 uL mono-clonal anti-flag

Shook for 1 hour

Rinsed with wash buffer and shook for 10 mins (3x)


In a tube mixed 1 mL of each developer and add to blot

Shook for 5 minutes


Placed between two layers of plastic wrap

Measured on imaging machine using Chemi protocol at different exposures


### **Blot 1**




Exposure Time: 0.38 secs



Exposure Time: 2 secs




Exposure Time: 45 secs




Exposure Time: 60 secs


## **Blot 2**




Exposure Time: 2 secs



Exposure Time: 10 secs



Exposure Time: 60 secs



Exposure Time: 10 mins

### **Gibson Cloning Globin pgRNA**

Mixture

9ul Vector

2 ul Insert

10 ul Gibson mix

Took pictures

Miniprepmed more pgRNA

Set up+ sent off Biobrick ACG+ATG eGFP plasmids for sequencing

- Harvest APOBEC plate (for Western blot)
  - 1.5 mL 2xSDS sample buffer + 30  $\mu$ L DTT
  - Heat at 80 °C, 10 min
  - Store in freezer



August 12, 2016

**Gibson Globin Counts/Liquid Cultures**

Neg. Control #1,2,3 – 0

Globin

#1 – 0


#2 – 0

#2 – 1

**Miniprep**

**Digested Scrambled and ACG-GFP pgRNAs with BstX1 and xHo1**

**Ran on Gel**



Lane 1 – Ladder

Lane 2 – Neg. Control

Lane 3 – Scrambled #1

Lane 4 – Scrambled #2

Lane 5 – Scrambled #3

Lane 6 – Scrambled #4

Lane 7 – ACG-GFP #1

Lane 8 – ACG-GFP #2

**Analysis**

The Scrambled #1 and ACG-GFP #2 looked like they worked because the insert is 150bp and these are slightly bigger than the negative control