

Week 8

July 25, 2016

What We Did Today:

Ran Gibson Assembly with ADARs and Neg Control (Using 0.8 of the original concentration)

Made Liquid Cultures from Gibson Apobec Plates

Will Digest with EcoR1 and xHo tomorrow to run on gel

Gibson Assembly Mixtures

ADAR 1-1x

0.98 ul pRETRo-ON
1.52 ul dCas9 Insert
2.12 ul ADAR 1-1x Insert
8 ul Mastermix
3.38 ul Water

ADAR 1-2x

0.98 ul pRETRo-ON
1.52 ul dCas9 Insert
3.56 ul ADAR 1-2x Insert
8 ul Mastermix
1.94 ul Water

ADAR 1-3x

0.98 ul pRETRo-ON
1.52 ul dCas9 Insert
1.62 ul ADAR 1-3x Insert
8 ul Mastermix
3.88 ul Water

ADAR 2-1x

0.98 ul pRETRo-ON
1.52 ul dCas9 Insert
3.65 ul ADAR 2-1x Insert
8 ul Mastermix
1.85 ul Water

ADAR 2-2x

0.98 ul pRETRo-ON
1.52 ul dCas9 Insert
3.48 ul ADAR 2-2x Insert
8 ul Mastermix
2.02 ul Water

ADAR 2-3x

0.98 ul pRETRo-ON
1.52 ul dCas9 Insert
3.35 ul ADAR 2-3x Insert
8 ul Mastermix
2.15 ul Water

Negative Control

0.98 ul pRETRO-ON
8 ul Mastermix
7.02 ul Water


APOBEC 1x

0.98 ul pRETRO-ON
1.52 ul dCas9 Insert
3.34 ul APOBEC 1x Insert
8 ul Mastermix
2.16 ul Water

Colony PCR again

One successful wild type clone with a clear insert

Picked it and grew in liquid culture



Pick colonies from 10uL ligation for liquid cultures

Transformed eGFP samples

Set up 12 well plate of HEK293T cells for transfection

- Digest
dCas9-ADAR ligations (7 samples from 7/22)
- Gel
Run for >1.5 hours

July 26, 2016

What We Did Today:

Counted Gibson Plates

Made liquid cultures of Gibson plates

Mini-prepped Apobec 2x and 3x Gibson

Digested Apobec 2x and 3x with EcoR1 and xHo1

Ran Digests on gel

Interviews

Gibson Plate Counts

Adar 2-1x

1.) 2

2.) 0

3.) 3

Adar 2-2x

1.) 1

2.) 4

3.) 2

Adar 1-2x

1.) 10

2.) 7

3.) 18

Adar 2-3x

1.) 11

2.) 9

3.) 2

Neg Control

1.) 0

2.) 0

3.) 0

Apobec 1x

1.) 5 (none single)

2.) 8

3.) 1

Adar 1-1x

1.) 11

2.) 7

3.) 3


Adar 1-3x

1.) 6

2.) 8

3.) 12


Gel 1



Gel 1: Apobec 2x
Lane 1 – Empty
Lane 2 – Ladder
Lane 3 – Neg. Control
Lane 4 – Apobec 2x #1
Lane 5 – Apobec 2x #2
Lane 6 – Apobec 2x #3
Lane 7 – Apobec 2x #4
Lane 8 – Apobec 2x #5

Analysis: The Apobecs ran correctly!! They are around where we would expect the 12kb marker to be. We realized that the EcoR1 site was cut out during the Gibson Assembly.

Gel 2




Gel 2: Apobec 3x
Lane 1 – Ladder
Lane 2 – Neg. Control
Lane 3 – Apobec 3x #1
Lane 4 – Apobec 3x #2
Lane 5 – Apobec 3x #3
Lane 6 – Apobec 3x #4
Lane 7 – Apobec 3x #5
Lane 8 – Empty
Lane 9 – Frederick Neg
Lane 10 – Frederick test


Analysis: The Apobecs ran correctly!! They are around where we would expect the 12kb marker to be. We realized that the EcoR1 site was cut out during the Gibson Assembly.


Wild type colony digested with XhoI and BamHI

AND THE INSERT DISAPPEARED
Lanes: Negative, WT clone



screened 21 more PTC colonies to try to find a clone
didn't find one





picked WT colonies 1-7 and PTC 7 just in case something got mixed up

Test digest:

14uL water

2uL DNA

2uL buffer

1uL each enzyme

Made liquid cultures of eGFP colonies in 2x AMP

- Glycerol stocks
- Miniprep
 - dCas9-ADAR/APOBEC ligations (8 samples)
 - Concentration: ~130-330 ng/ μ L
- Digest
 - Enzymes: ECoRI, XbaI
 - 0.5 μ L each enzyme
 - 4 μ L DNA (dCas9-ADAR/APOBEC)
 - 1 μ L CutSmart buffer
 - (10 μ L total)
- Gel

July 27, 2016


What We Did Today:

Mini-prepped Apobec 1x and Adar 1-1x, 1-2x, 1-3x, 2-1x, 2-2x, 2-3x Gibson

Digested Apobec 1x and Adar 1-1x, 1-2x, 1-3x, 2-1x, 2-2x, 2-3x xH01

Ran Digests on gel

Made 2x Amp Plates



Gel 1

Lane 1 – Ladder

Lane 2 – Neg. Control

Lane 3 – Adar 1-1x #1

Lane 4 – Adar 1-1x #2

Lane 5 – Adar 1-1x #3

Lane 6 – Adar 1-2x #1


Lane 7 – Adar 1-2x #2

Lane 8 – Adar 1-2x #3

Lane 9 – Apobec 1x #1

Analysis

Did not run well



Gel 2

Lane 1 – Ladder

Lane 2 – Neg. Control

Lane 3 – Adar 1-3x #1

Lane 4 – Adar 1-3x #2

Lane 5 – Adar 1-3x #3

Lane 6 – Adar 2-1x #1

Lane 7 – Adar 2-1x #2

Lane 8 – Adar 2-1x #3

Lane 9 – Apobec 1x #2

Analysis

No clones



Gel 3

Lane 1 – Ladder
 Lane 2 – Neg. Control
 Lane 3 – Adar 2-2x #1
 Lane 4 – Adar 2-2x #2
 Lane 5 – Adar 2-2x #3
 Lane 6 – Adar 2-3x #1
 Lane 7 – Adar 2-3x #2
 Lane 8 – Adar 2-3x #3
 Lane 9 – Apobec 1x #3

Analysis

No clones


Colony PCR the 8 liquid cultures, 10uL of liquid culture brought up to 25uL with 15uL water, boiled then typical procedure

digested 2uL of each (after miniprep) with XhoI and BamHI


Gel 1:

Digest:

Ladder, WT 1-7, PTC 7



Gel 2:
Colony PCR
Ladder, Positive control, WT 1-7, PTC 7




The globin came back, then was gone again

July 28th, 2016

What We Did Today:
PCR dCas

PCR Mixture

- 1.) Primers were diluted 1:10 with water (9ul water, 1ul Primer)
- 2.) Components were added to PCR tubes
 - 1.) 22ul Water
 - 2.) 1ul Template (ADAR 1, ADAR 2, and APOBEC with all extensions)
 - 3.) 1ul Primer 1 (3' XTen)
 - 4.) 1ul Primer 2 (5' dCas)
 - 5.) 25 ul MasterMix
- 3.) Tubes were run in the PCR Machine under the iGEM Colony Protocol for 50ul tubes
- 4.) 5ul of each sample was run on a gel



Gel – dCas9 Insert Testing

Lane 1 – Empty

Lane 2 – Ladder

Lane 3 – Empty

Lane 4 – Old dCas9 Insert

Lane 5 – Empty

Lane 6 – New dCas9 Insert

Analysis

The old dCas9 insert ran correctly however the new one did not. We can only see the primer dimers for this sample.

- Fixing cells
 - eGFP trial 1
 - 6 slides
 - Positive, negative, ATG, ACG
- Imaging
 - 2 fields/cover slip

- Green/red fluorescence channel each
- 5X objective

48 images total

Sequenced globin trials, editing enzymes

Examined transfected cells

Fixed and mounted transfected cells

Created standard curve for Interlab using FITC and made standardization samples for Interlab

Measured samples using plate reader at gateway

July 29, 2016

Gibson Mixtures

Adar 1-1x

1.23 ul pRETRO-ON
1.90 ul dCas Insert
2.65 ul Adar Insert
10 ul Mastermix
4.22 ul Water

Adar 1-2x

1.23 pRETRO-ON
1.90 dCas9 Insert
4.44 ul Adar Insert
10 ul Mastermix
2.43 ul Water

Adar 1-3x

1.23 ul pRETRO-ON
1.90 dCas9 Insert
2.03 ul Adar Insert
10 ul Mastermix
4.84 u. Water

Adar 2-1x

1.23 pRETRO-ON
1.90 dCas9 Insert
4.57 ul Adar Insert
10 ul Mastermix
2.3 ul Water

Adar 2-2x

1.23 ul pRETRO-ON
1.90 dCas9 Insert
4.35 ul Adar Insert
10 ul Mastermix
2.52 ul Water

Adar 2-3x

1.23 ul pRETRO-ON
1.90 dCas9 Insert
4.18 ul Adar Insert
10 ul Mastermix
2.69 ul Water


Apobec 1x

1.23 ul pRETRO-ON
1.90 dCas9 Insert
4.18 ul Apobec Insert
10 ul Mastermix
2.69 ul Water

Neg. Control

1.23 pRETRO-ON
10 ul Mastermix
8.77 ul Water

Test digest of globin vector with BglII, XhoI, HindIII, BamHI, EcoRI, PstI to test for efficiency of the cutting of the enzyme
Digest 1 hr, SAP heat kill ran 10uL on a gel



ligate 5 mins 1uL

Gel:

L, neg, Bam, Bgl, Eco, Hind, Pst, Xho

Mini prepped additional pgRNA

Maintenance of HEK293T cells