

# June 20, 2016

## Notes on NEGEM Presentation

- Next slide should explain what rice crispr is
- Background slides before defending
- In text citations

sgRNA is produced by RNA Pol III  
mCherry on separate promotor  
Transient transfection in HEK 293T  
rtTa tetracyclin repressor

## Outreach activities

- ~3 slides
- Say that we have notes to share with W&M and will put on wiki
- Add in touch tomorrow slides, pictures, graphs

Add acknowledgements

Wednesday 3pm run through slides

Draft e-mail to UMASS Zhang, CRISPR Technologies, Team at Harvard/MIT

Start stable integration

Transform and make mini prep

Transfect + add antibiotic

What We Did Today

- Split Cells

- TouchTomorrow Data

- NEGEM Slides

- Transformation for interlab study

Transformation Protocol for interlab study

- Used Plate 2- cell 6F- Bba-R0040 (Cam<sup>r</sup>)

  - 10 microliter water and pipetted up and down, put in microcentrifuge tube

  - 1 microliter pRetro on (amp<sup>r</sup>)

- Used 50 microliters of dh5a cells and added all of the above to 200 microliters of LB

- Capped tubes and let shake for an hour and a half

- Plated cells and left overnight

GM, JW

Worked on NEGEM PowerPoint

# June 21, 2016

What we did today:

- NEGEM Slides

- Interlab study:

  - Checked transformation

  - Made liquid cultures

Liquid Cultures:

- Chose 6 cultures from the pRetro Plates, avoided major satellite cells

- 5mL LB

- 2.5 microliters Ampicillin

- Added cultures to 15 mL conical tube

- Put them in the shaking incubator

GM, JW

Split cell culture (large flask), 1:5

# June 22, 2016

## What We Did Today:

- Finished NEGEM presentation

- Looked up people to contact for Integrated Practices

- Practiced presentation for NEGEM

Mini Prepped P-retro


June 23, 2016

What We Did Today:

Ran Gel with Restriction Enzymes


Gel 1: Singles

- Lane 1 – Empty
- Lane 2 – Ladder
- Lane 3 – Neg. Control
- Lane 4 – BamH1
- Lane 5 – Not1
- Lane 6 – Pac 1
- Lane 7 – Spe 1 + Cal1



Gel 2: Multiple

- Lane 1 – Empty
- Lane 2 – Ladder
- Lane 3 – Neg. Control
- Lane 4 – EcoR1 + BamH1
- Lane 5 – xBa1 + HindIII
- Lane 6 – EcoRI + Pst1
- Lane 7 – xHo1 + Spe1



### Restriction Digest of Plasmids

- BamHI+EcoRI
- PacI
- BamHI
- NotI

15-16uL water  
2uL SmartCut  
1uq plasmid  
1uL enzyme

Heat incubate 1-2 hours  
Run gel

Big flask - bacterial contamination  
Small - fungal  
Threw out all media

Streaked humanized dcas9 on amp and 2x amp plates

**June 24, 2016**

**What We Did Today:**

Attended NEGEM Meeting #1